

ChIP-seq 実験にベリタス取扱製品を使用 研究者の声【16】

東京大学 山中 総一郎 先生

プロトコル紹介

【実際に行っている Dynabeads Protein G を利用した ChIP-seq のプロトコル】

1. 細胞の固定

- ↓ FACS で生殖細胞を 10000 細胞分取する (in low-binding tube)
- ↓ remove supernatant
- ↓ suspend with 1ml DMEM (10% FBS)
- ↓ add 67 μ l, 16% Formaldehyde (Pierce) (final 1%)
- ↓ RT, 10min, rotation
- ↓ add 62.5 μ l 2M glycine (final 125mM)
- ↓ RT, 5min, rotation
- ↓ centrifuge, 6000rpm, 1min x2
- ↓ take supernatant (~900 μ l)
- ↓ centrifuge, 6000rpm, 1min
- ↓ take supernatant completely > stock in -80°C freezer

2. Dynabeads のブロッキング

- ↓ take 20 μ l Dynabeads Protein G into two tubes per sample (one for IP and one for pre-clear)
- ↓ wash with 800 μ l, 0.5% BSA in PBS, twice
- ↓ suspend with 1ml, 0.5% BSA in PBS
- ↓ 4°C, rotation, 1h

3. IP 用 Dynabeads の調製

- ↓ add 2 μ g anti-body (H3K27me3)
- ↓ 4°C, rotation, o/n

4. ソニケーション

Frozen cells

- ↓ add 1ml Swelling Buffer, 4°C, rotation, 20min
- ↓ centrifuge, 8000rpm, 1min x2
- ↓ take supernatant (~900 μ l)

↓ centrifuge, 8000rpm, 1min
↓ take supernatant completely
↓ suspend with 100μl 1x Sonication buffer (Covaris)
↓ sonication [Branson: 30%, total ON time=2min, ON/OFF=2sec/3sec]
↓ centrifuge, 13500rpm, 5min
↓ save supernatant to low-binding tube
↓ centrifuge, 13500rpm, 10min
↓ save supernatant to low-binding tube
↓ dilute with 1ml of Abcam RIPA buffer

5. pre-clear

2 で調製した“pre-clear beads”

↓ wash with 1ml, Abcam RIPA buffer x2 回
↓ suspend with 20μl Abcam RIPA buffer
↓ add to cell extract
↓ 4°C, rotation, 1h
↓ save supernatant
↓ aliquot 100μl for [WCE] > the rest of the cell extract(~1ml) is for IP

6. IP

3 で調製した“beads conjugated with anti-body”

↓ wash with 1ml, Abcam RIPA buffer x4 回
↓ suspend with 20μl Abcam RIPA buffer
↓ add beads to 1ml cell extract
↓ 4°C, rotation, o/n

7. Beads wash > elution

IP を行った beads

↓ wash with 1ml, LOW buffer x3 回
↓ wash with 1ml, HIGH buffer x3 回
↓ add 200μl direct elution buffer
↓ incubate at 65°C, 15min, with constant vortexing
↓ save supernatant

WCE

100µl of WCE

↓ add 100µl direct elution buffer

↓ incubate at 65°C, 15min, with constant vortexing

8. ProK treatment > reverse crosslinking

IPed sample and WCE sample

↓ transfer to screw tube

↓ add 5µl proteinase K (Roche)

↓ incubate 37 °C, >6h

↓ 65°C, >6h

9. DNA 回收

↓ add 20µl 3M NaOAc, 200µl Phenol:Chloroform (Thermo)

↓ centrifuge, 13500rpm, 10min

↓ save supernatant to low-binding tube

↓ +1µl pellet paint (Merck)

↓ -80°C, 15min

↓ centrifuge, 14000rpm, 30min

↓ remove supernatant

↓ add 150µl, 70% ethanol

↓ centrifuge, 14000rpm, 5min

↓ remove supernatant

↓ dry up RT, 10min

suspend with 35µl 0.1M Tris-HCl.


10. ChIP-seq library preparation

Check the quality of IPed sample with qPCR

↓ prepare the DNA library for sequencing with QIAGEN QIAseq Ultralow Input library Kit

According to the manufacturer's instruction

↓ check the DNA size with Bioanalyzer



使用した Buffer の組成

- Swelling buffer

20 mM Hepes(pH7.9)

1.5 mM MgCl₂

10 mM KCl

0.1% NP-40

1mM DTT (add before use)

- Abcam RIPA buffer

(RIPA buffer in Abcam protocol)

50 mM Tris-HCl, pH8.0

150 mM NaCl

2 mM EDTA pH8

1% NP-40

0.5% Sodium Deoxycholate (DOC)

0.1% SDS

- Low buffer

(Wash buffer in Abcam protocol)

0.1% SDS

1% Triton X-100

2 mM EDTA pH8

150 mM NaCl

20 mM Tris-HCl, pH8.0

- High buffer

(Final wash buffer in Abcam protocol)

0.1% SDS

1% Triton X-100

2 mM EDTA pH8

500 mM NaCl

20 mM Tris-HCl, pH8.0

- Direct elution buffer

10 mM Tris-HCl, pH8.0

5 mM EDTA

300 mM NaCl

0.5% SDS

株式会社ベリタス 〒105-0013 東京都港区浜松町1-10-14 住友東新橋ビル3号館5F

TEL 03-5776-0078 FAX 03-5776-0076

技術的なお問い合わせは：TEL 03-5776-0040 E-mail techservice@veritastk.co.jp