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
HLA Fusion™

Create and Update Demographic/Allele Frequency

Create Demographic Group with Alleles and Frequency Delimited

Demographic Frequency File : FUSION\data/catalog/ver1_2020_Jan_allelefilter.csv

Select Name	Allele/Sero : 2020_Ja
<input checked="" type="checkbox"/> ver1_Japanese_2020_Jan_GF2020JSHI	A*01:01 0.439
	A*01:0... 0.439
	A*01:0... 0.439
	A*01:0...
	A*01:0...
	A*01:0...

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
New Duplicate Delete Export Save Close

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
HLA Fusion™

Product Type: LABType

Code: NMDP Local Code P Group G Group No Code

Cross Code: Enable Cross Code (A,B,C,DRB1,DRB345,DQ) Enable Cross Code (DP)

Demographic: ver_1_2020_Jan_allelefilter [Edit] Number of False Rns: 1

Possible Allele Code: G1 G2 G3

Bw4/Bw6 Information

Computer Assigned Serology

Flag Sero Ambiguity

Skip Batch Analysis of Samples with

Min Positive Control *: 1000

Min Bead Count *: 100

Min Bead Failure Threshold (%): 10

Close Bead Rxn Threshold *: 3

Display Popup message for Low Bead Count and Low Positive Control

Allow Auto Accept All

Save Force 1 Pair (If not checked, the Force 1

Display Converted HD Bead Count in Raw Dat

Exclude Exon 3 Probes for DQ/DP Analysis

Global Bead Adjustment (Impact All Samples)

Include Generic Test Comments (Exon 4-7)

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
HLA Fusion™

Product Type: MicroSSP

Code: NMDP Local Code P Group G Group No Code

Cross Code: Enable Cross Code (A,B,C,DRB1,DRB345,DQ) Enable Cross Code (DP)

Demographic: ver_1_Japanese_2017_Jan_G [Edit] Number of False Rns: 1

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
①「Upgrade prior versions of HLA Fusion database to 4.4」をクリックします

... Create/Select Audit Log Database
... Merge Audit Log database
... Current Database Details
... Select/Connect to HistoTrac Database

Data Migration

... Upgrade prior versions of HLA Fusion database to 4.4

Database

Database Name : FUSION44_0820

Version : 4.4.0.13925, Created on: 2020/08/20

Size : Used: approximately 1% - 13 MB of 10240 MB DB size.

Settings : us_english#SQL_Latin1_General_CP1_CI_AS
SQL Express Edition

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Database Name : FUSION44_0820

Version : 4.3.1.13925, Created on: 2020/05/25

Size : Used: approximately 1% - 62 MB of 10240 MB DB size.

Settings : us_english#SQL_Latin1_General_CP1_CI_AS
SQL Express Edition

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Database

Database Name : FUSION44_0820

Version : 4.4.0.13925, Created on: 2020/08/20

Size : Used: approximately 1% - 100 MB of 10240 MB DB size.

Settings : us_english#SQL_Latin1_General_CP1_CI_AS
SQL Express Edition

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