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

[LABScreen Single Antigen Class I](#)
[LABScreen Single Antigen ExPlex Class I](#)
[LABScreen Single Antigen Supplement Class I](#)
[LABScreen Single Antigen Class II](#)

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


④

⑤

Check Control Import Delete Patient Close

Well	Sample	Sample Date	Sample Treatment	Sample Source	Dilution Factor	Secondary Ab	Luminex Min Bead Cnt	NS	Exist In DB	Patient ID	First Name	Last Name	Ethnicity	Patient/Donor
78(1.A1)_NC							100	<input checked="" type="checkbox"/>	N					
74(1.B1)_SH8101							100	<input type="checkbox"/>	N					
76(1.C1)_SH8102							100	<input type="checkbox"/>	N					
76(1.D1)_SH8103							100	<input type="checkbox"/>	N					
77(1.E1)_SH8104							100	<input type="checkbox"/>	N					
78(1.F1)_SH8105							100	<input type="checkbox"/>	N					
79(1.G1)_SH8106							100	<input type="checkbox"/>	N					

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

PC: (002) 9094.41
 NC: (001) 609.91
 PC/NC: 14.91

NC Inrshold (CL1): 50(50)
 NC Threshold (CL2): 50(50)

Class I Computer Assignment Positive
 Final Assignment Positive Negative Undetermined

Class II Computer Assignment Positive
 Final Assignment Positive Negative Undetermined

MIC Computer Assignment Negative

Comments/ (System): NC Raw >=500.NC Raw >=500Adjusted Cutoff.

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


Analyze Data **Reports** Data Sample Patient Info Profile Utilities Help Exit

Home, Print, Search, etc. icons

Patient Generic Typing LABType MicroSSP Generic Antibody LABScreen Specialty Statis

Sample Summary LSM Details LSM Summary LSM Overview Product Compare

Patient or Donor ID: * Session: * Batch: * Sample ID: *

SessionID WellPosition SampleID PatientID

LSM Detail

Approved By: _____ Date: _____

Session ID: 200907_LSM022_NC023_20200907_124154
Catalog: LSM12NC23_022_00

NS Sample: NC

Background Value	005	007	132.31	008	129.79	013	114.74	018	128.56	030	129.09	041	123.49	043	133.09	053	133.52	059	126.96
061	128.84	062	143.13	078	92	081	109.28	083	107.59	084	66.02	096	95.5	099	90.6	100	90.09	NC	114.76
PC	1212247																		

Class I Positive Cutoff: 1.5 Negative Cutoff: 1.2 NC Threshold: 50 Class I Sensitivity/MFI: 0
Class II Positive Cutoff: 1.5 Negative Cutoff: 1.2 NC Threshold: 50 Class II Sensitivity/MFI: 0
MIC Positive Cutoff: 1.5 Negative Cutoff: 1.2 NC Threshold: 50 HNA Sensitivity/MFI: 0

Patient ID: _____ Patient Name: _____ Local ID: _____

Patient ID:	Status	Transplant Type						Class II						Sample Date									
		Class I			Class II			Class II			Class II												
Overall Positive	Bead	Result	Raw Data	Ratio	Count	Overall Positive	Bead	Result	Raw Data	Ratio	Count	Overall Positive	Bead	Result	Raw Data	Ratio	Count	Overall Positive	Bead	Result	Raw Data	Ratio	Count
	005	Positive	831.5	16.21	112		078	Undetermined	95.78	1.49	157		081	Positive	232.02	4.22	163		083	Negative	68.14	0.94	155
	007	Positive	977.6	19.13	192		084	Positive	137.78	2.33	152		096	Negative	75.5	1.09	151						
	008	Positive	622.27	12.02	147																		
	013	Positive	849.77	16.57	101																		
	018	Positive	951.22	18.60	100																		
	030	Positive	775.82	15.09	148																		
	041	Positive	703.83	13.65	169																		
	043	Positive	885.08	17.28	162																		
	053	Positive	706.84	13.71	197																		
	059	Positive	698.43	13.55	142																		
	061	Positive	650.73	12.59	159																		
	062	Positive	942.17	18.42	184																		

MIC

Overall Negative	Bead	Result	Raw Data	Ratio	Count
	099	Negative	24.83	0.07	147
	100	Negative	22.72	0.03	170

Sec. Ab: _____ Treatment: _____
Saved By: 1,1 Saved Date: 11/22/2020
Confirmed By: _____ Confirmed Date: _____
Comment: _____
User Comment: _____

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Analyze Data **Reports** Data Sample Patient Info Profile Utilities Help Exit

Patient Generic Typing LABType MicroSSP Generic Antibody LABScreen Specialty Statistics

Sample Summary

- LSM Details
- LSM Summary**
- LSM Overview
- Product Compare

SessionID
WellPosition
SampleID
CatalogID

LSM Summary

Approved By: _____ Date: _____

Session ID: 200907_LSM022_NC023_20200907_124154
Catalog ID: LSM12NC23_022_00

NS Sample: NC

Background Values

005	121.21	007	132.31	008	129.79	013	114.74	018	128.56	030	129.09	041	123.49	043	133.09	053	133.52	059	126.96
061	128.84	062	143.13	078	92	081	109.28	083	107.59	084	66.02	096	95.5	099	90.6	100	90.09	NC	114.76

PC 1212247

Class I Positive Cutoff: 1.5 Negative Cutoff: 1.2 NC Threshold: 50 Class I Sensitivity/MFI: 0
Class II Positive Cutoff: 1.5 Negative Cutoff: 1.2 NC Threshold: 50 Class II Sensitivity/MFI: 0
MIC Positive Cutoff: 1.5 Negative Cutoff: 1.2 NC Threshold: 50 HNA Sensitivity/MFI: 0

[WellPosition]	Sample ID	Class I	Class II	MIC	Raw	Count	Raw	Count	PC/NC
[2(1,B1)]	FL1-64	Positive	Positive	Negative	21.16	168	9278.56	148	438.5

Patient Name: _____
Patient ID: _____
Sec. Ab: _____ Treatment: _____
Saved By: 1,1 Saved Date: 11/22/2020
Confirmed By: _____ Confirmed Date: _____
Comment: _____
User Comment: _____

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
The screenshot shows the HLA Fusion software interface. The 'Reports' menu is highlighted in the top navigation bar. Below it, the 'Generic Antibody' sub-menu is open, showing options like 'Antibody Custom', 'Antibody Screening/Identification', and 'Antibody Screening/Single Antigen'. The 'Custom Antibody Screening Report Setup' dialog box is open, displaying various configuration options for the report. A red box highlights the 'Setup' button in the bottom right corner of the dialog box. A text box with the Japanese text 'レポートに出力する項目を選択' (Select items to output to the report) is overlaid on the dialog box.

Reports Data Sample Patient Info Profile Utilities Help Exit

Generic Typing LABType MicroSSP Generic Antibody LABScreen Specialty Statistical Miscellaneous My Favorite Tools

Antibody Custom
Antibody Screening/Identification
Antibody Screening/Single Antigen

Custom Antibody Screening Report Setup

HLA Fusion™

Type or enter the report name:
LABScreen

Lab Information

Patient Information
Check All **Uncheck All**
 Patient ID, Name (REQUIRED for this group)
 Project Type
 SSN
 Ethnicity, DOB, Gender
 Disease, Status, Blood, Transplant Type
 Address, Phone, Email
 Spouse, Emg. Contact, Employer
 DCN, Hospital, Division
 UNOS CPRA
 Donor Info
 Molecular Typing Sero Typing
 Antibody Assignment
 Match / Mismatch
 Add MFI Baseline
 Acceptable Antigens
 Unacceptable Antigens
 AM DFC ETKAS DFC
 Unacceptable Epitopes

Sample Information
Check All **Uncheck All**
 Sample ID (REQUIRED) Barcode
 System Comments User Comments
 Sample Source Dilution Factor

Session Information
Check All **Uncheck All**
 Session ID/Date (REQUIRED for this group)
 Lumindex Info
 Session Details: Catalog ID, Locus Type, Test Position, Date, NS Sample
 More Tests and Tests Remarks
 Quantiplex/SFI

Test Configuration
Check All **Uncheck All**
 User Settings:
Excluded Antigens, Positive Region Threshold, Normalization Formula, Allele Level
 NU and PG Values, % PHA-or- %SA and Cutoffs by Region

Overall Results/Assignments
Check All **Uncheck All**
 Overall Pos/Neq Other
 % PRA % Donor PRA
 Antibody New Antibody
 Possible Epitope Amino Acid Position

Epitope Analysis Results
Check All **Uncheck All**
 Antibody Specificity (REQUIRED for this group)
 CREG Analysis Details:
TP, FN, Mean (Raw) of Positives
 CREG Chart, Circle DSA

Tail Analysis Results
Check All **Uncheck All**
 Antibody Specificity (REQUIRED for this group)
 Tail Analysis Details:
TP/FP/FN/TN, R Value, Average Score, % Inclusion, Strength Index
 Manual Tail Assignment

Test and Catalog Details
Check All **Uncheck All**
 Bead ID or Test Well (REQUIRED for this group)
 Test Values: Rxn, Raw, Normalized, Count
 SFI Raw, SFI Normal
 Sero Specificity **Sort by:** Bead ID
 ASC DESC
 Allele Specificity

Cutoff Summary
 Bead MFI Chart
 Raw Data Baseline
 Sort Antigen Using HLA Order
 Rxn Score Sort Alphanumerically

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ONE LAMBDA

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Custom Antibody Screening Report Setup

HLA Fusion™

Type or enter the report name*: LABScreen

Lab Information

Patient Information

- Patient ID, Name (REQUIRED for this group)
- Project Type
- SSN
- Ethnicity, DOB, Gender
- Disease, Status, Blood, Transplant Type
- Address, Phone, Email
- Spouse, Emg. Contact, Employer
- DCN, Hospital, Division
- UNOS CPRA
- Donor Info
- Molecular Typing Sero Typing
- Antibody Assignment
- Match / Mismatch
- Add MFI Baseline
- Acceptable Antigens
- Unacceptable Antigens
- AM DFC ETKAS DFC
- Unacceptable Epitopes

Test Configuration

- User Settings
- Excluded Antigen Threshold, N Level
- NO and PC V and Cutoffs

Overall Results/Assessment

- Overall Pos/ % PRA
- Antibody
- Possible Epitope

Epitope Analysis Report

- Antibody Specificity
- CREG Analysis: TP, FN, Mear
- CREG Chart

Sample Information

- Sample ID (REQUIRED)
- System Comments
- Sample Source
- Barcode
- User Comments
- Dilution Factor

Session Information

- Session ID/Date (REQUIRED for this group)
- Luminex Info
- Session Details: Catalog ID, Locus Type, Test Position, Date, NS Sample
- More Tests and Tests Remarks
- Quantiplex/SFI

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Approved By:

Lab Name 1
Institute:
Street: 1
Lab City:
Region:
Country:
Notes:

Lab Code
Contact: 1
Email: chiefujwara@veritastk.co.jp
Phone:
Fax:

Province:
Postal Code:

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Patient Name:
Project Type
Ethnicity:
DOB: Gender: UNK SSN:
Disease
Transplant Type
City:
Address
Phone:
Email:
Blood Type
Spouse:
Emergency Contact:
Employer: Division:
Hospital
UNOS CPRA: as of
Donor PRA:
Donor Type

Donor Center ID:
From Other Facility: NO

Patient/Donor:
AMDFC:
ETKAS DFC:
Molecular Typing
Sero Typing
Antibody Assignment Class I:
Class II:
MIC

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Patient Donor Match/ Mismatch

Acceptable Antigens
Unacceptable Antigens
Unacceptable Epitope- Class I:
Unacceptable Epitope- Class II:
Unacceptable Epitope- MICA:

Session ID: 200908_LS1A04011_NC023-2New_20200908_142249 Luminex: xPONENT/FlexMAP 5D-4.5.2290 / SN BMBDD14251002
Catalog: LS1A04NC23_011_00 Test Pos: 74(1,B10)

Quantiplex/SFI Used

Sample ID: FL1-PC Local ID:

Sec. Ab: Sample Date:  Test Date: 9/8, 2020

Treatment

Sample Source: Dilution Factor:

Saved By: 1,1 Date: 9/8, 2020 NS Sample: NC

Confirmed By: Date: Test Type: Single Class I

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Custom Antibody Screening Report Setup

HLA Fusion™

Type or enter the report name*: LABScreen

Lab Information

Patient Information

Patient ID, Name (REQUIRED for this group)
 Project Type
 SSN
 Ethnicity, DOB, Gender
 Disease, Status, Blood, Transplant Type
 Address, Phone, Email
 Spouse, Emg. Contact, Employer
 DCN, Hospital, Division
 UNOS CPRA
 Donor Info
 Molecular Typing Sero Typing
 Antibody Assignment
 Match / Mismatch
 Add MFI Baseline
 Acceptable Antigens
 Unacceptable Antigens
 AM DFC ETKAS DFC
 Unacceptable Epitopes

Test Configuration

User Settings:
 Excluded Antigens, Positive Threshold, Normalization Factor
 NC and PC Values, % PRA and Cutoffs by Region

Overall Results/Assignments

Overall Pos/Neg Overall %
 % PRA % Negative
 Antibody Negative
 Possible Epitope Antigen

Epitope Analysis Results

Antibody Specificity (REQUIRED)
 CREG Analysis Details: TP, FN, Mean (Raw) of Positives
 CREG Chart, Circle DSA

Sample Information

Sample ID (REQUIRED) Barcode
 System Comments User Comments
 Sample Source Dilution Factor

Session Information

Session ID/Date (REQUIRED for this group)
 Luminex Info
 Session Details: Catalog ID, Locus Type, Test Position, Date, NS Sample
 More Tests and Tests Remarks
 Quantiplex/SFI

Donor ID: 222
 Serology Typing: A2,A26,B8,B60,Bw6,Cw7,Cw8,DR4,DR12,DR52,DR53,DQ7,DQ8,DP2,DP6
 Molecular Typing: A*02:06,A*26:01,B*08:01,B*40:01,C*07:02,C*08:01,DRB1*04:07,DRB1*12:01,DQA1*03:01,DQA1*05:05,DQB1*03:01,DQB1*03:02,DPA1*01:03,DPA1*01:03,DPB1*02:01,DPB1*06:01

Serology Typing	Molecular Typing	MFI	(Baseline)	Match/ Mis Match
A2				Mismatch not in Final Assignment
	A*02:01	14.89	0	
	A*02:03	12.47	0	
	A*02:06	20.62	0	Mismatch not in Final Assignment
A26				Mismatch not in Final Assignment
	A*26:01	32.77	0	Mismatch not in Final Assignment
B8				Mismatch in Final Assignment
	B*08:01	2315.36	2253	Mismatch not in Final Assignment
B60				Mismatch in Final Assignment
	B*40:01	9004.45	8930.45	Mismatch not in Final Assignment
Cw7				Match
	C*07:02	135.77	2.55	Match
Cw8				Mismatch in Final Assignment
	C*08:01	4698.46	4609.17	Mismatch not in Final Assignment

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ONE LAMBDA

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Custom Antibody Screening Report Setup

HLA Fusion™

Type or enter the report name*: LABScreen

Lab Information

Patient Information

Patient ID, Name (REQUIRED for this group)

Project Type

SSN

Ethnicity, DOB, Gender

Disease, Status, Blood, Transplant Type

Address, Phone, Email

Spouse, Emg. Contact, Employer

DCN, Hospital, Division

UNOS CPRA

Donor Info

Molecular Typing Sero Typing

Antibody Assignment

Match / Mismatch

Add MFI Baseline

Acceptable Antigens

Unacceptable Antigens

AM DFC ETKAS DFC

Unacceptable Epitopes

Sample Information

Sample ID (REQUIRED)

System Comments

Sample Source

Barcode

User Comments

Dilution Factor

Session Information

Session ID/Date (REQUIRED for this group)

Luminex Info

Session Details: Catalog ID, Locus Type, Test Position, Date, NS Sample

More Tests and Tests Remarks

Quantiplex/SFI

Test Configuration

User Settings:

Excluded Antigens, Positive Region Threshold, Normalization Formula, Allele Level

NC and PC Values, % PRA -or- %SA and Cutoffs by Region

Overall Results/Assignments

Overall Pos/Neg

% PRA

Antibody

Possible Epitope

Other

% Donor PRA

New Antibody

Amino Acid Position

Epitope Analysis Results

Antibody Specificity (REQUIRED for this group)

CREG Analysis Details: TP, FN, Mean (Raw) of Positives

CREG Chart, Circle DSA

Test Results

	Class I	Class II
Overall (Pos/Neg):	Positive	
%PRA - or - %SA:	15	
Antibody Assignment:	B48, B41, B81, B45, B50, B60, B72, B61, B78, B76, B62, B71, B18, B39, B35, B75, B*48:01, B*40:01	

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CREG Table

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Custom Antibody Screening Report Setup

HLA Fusion™

Type or enter the report name*: LABScreen

Lab Information

Patient Information

Patient ID, Name (REQUIRED for this group)

Project Type

SSN

Ethnicity, DOB, Gender

Disease, Status, Blood, Transplant Type

Address, Phone, Email

Spouse, Emg. Contact, Employer

DGN, Hospital, Division

UNOS CPRA

Donor Info

Molecular Typing Sero Typing

Antibody Assignment

Match / Mismatch

Add MFI Baseline

Acceptable Antigens

Unacceptable Antigens

AM DFC ETKAS DFC

Unacceptable Epitopes

Sample Information

Sample ID (REQUIRED) Barcode

System Comments User Comments

Sample Source Dilution Factor

Session Information

Session ID/Date (REQUIRED for this group)

Luminex Info

Session Details: Catalog ID, Locus Type, Test Position, Date, NS Sample

More Tests and Tests Remarks

Quantiplex/SFI

Test Configuration

User Settings:

Excluded Antigens, Positive Region Threshold, Normalization Formula, Allele Level

NC and PG Values, % PRA --or-- %SA and Cutoffs by Region

Overall Results/Assignments

Overall Pos/Neg Other

% PRA % Donor PRA

Antibody New Antibody

Possible Epitope Amino Acid Position

Epitope Analysis Results

Antibody Specificity (REQUIRED for this group)

CREG Analysis Details: TP, FN, Mean (Raw) of Positives

CREG Chart, Circle DSA

Tail Analysis Results

Antibody Specificity (REQUIRED for this)

Tail Analysis Details: TP/FP/FN/TN, R Value, Average Score, % Inclusion, Strength Index

Manual Tail Assignment

Test and Catalog Details

Bead ID or Test Well (REQUIRED for this)

Test Values: Rxn, Raw, Normalized, Count

SFI Raw, SFI Normal

Sero Specificity

Allele Specificity

Sort by: Bead ID

ASC DESC

Cutoff Summary



Bead MFI Chart

Raw Data Baseline

Sort Antigen Using HLA Order

Rxn Score Sort Alphabetically

Buttons: Check All, Uncheck All, Save, Delete, Cancel

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Custom Antibody Screening Report Setup

HLA Fusion™

Type or enter the report name*:
test

Lab Information

Patient Information
 Lab Information
Check **Uncheck**

- Patient ID, Name (REQUIRED for this group)
- Project Type
- SSN
- Ethnicity, DOB, Gender
- Disease, Status, Blood, Transplant Type
- Address, Phone, Email
- Spouse, Emg. Contact, Employer
- DCN, Hospital, Division
- UNOS CPRA
- Donor Info
- Molecular Typing Sero Typing
- Antibody Assignment
- Match / Mismatch
- Add MFI Baseline
- Acceptable Antigens
- Unacceptable Antigens
- AM DFC ETKAS DFC
- Unacceptable Epitopes

Sample Information
Check **Uncheck**

- Sample ID (REQUIRED) Barcode
- System Comments User Comments
- Sample Source Dilution Factor

Session Information
Check **Uncheck**

- Session ID/Date (REQUIRED for this group)
- Luminex Info
- Session Details: Catalog ID, Locus Type, Test Position, Date, NS Sample
- More Tests and Tests Remarks
- Quantiplex/SFI

Test Configuration
Check **Uncheck**

- User Settings:
Excluded Antigens, Positive Region Threshold, Normalization Formula, Allele Level
- NC and PC Values, % PRA -or- %SA and Cutoffs by Region

Overall Results/Assignments
Check **Uncheck**

- Overall Pos/Neg Other
- % PRA % Donor PRA
- Antibody New Antibody
- Possible Epitope Amino Acid Position

Epitope Analysis Results
Check **Uncheck**

- Antibody Specificity (REQUIRED for this group)
- OREG Analysis Details:
TP, FN, Mean (Raw) of Positives
- OREG Chart, Circle DSA

Tail Analysis Results
Check **Uncheck**

- Antibody Specificity (REQUIRED for this group)
- Tail Analysis Details:
TP/FP/FN/TN, R Value, Average Score, % Inclusion, Strength Index
- Manual Tail Assignment

Test and Catalog Details
Check **Uncheck**

- Bead ID or Test Well (REQUIRED for this group)
- Test Values: Rxn, Raw, Normalized, Count
- SFI Raw, SFI Normal
- Sero Specificity **Sort by** Sero Specificity ASC DESC
- Allele Specificity


Cutoff Summary
 Cutoff Summary

- Bead MFI Chart
- Raw Data Baseline
- Sort Antigen Using HLA Order
- Rxn Score Sort Alphanumerically

Check All **Uncheck** **Save** **Delete** **Cancel**

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
LABScreen Combined Analysis
Current Sample ID: 111
Patient ID: [redacted]
Analysis Date: 2019/03/04 - 2021/08/16


Select	Sample ID	Sample Date	Well Position	Session Name	Analysis Date	Final/Assign	W6-32	Quantplex	Sec Ab	Treatment
<input checked="" type="checkbox"/>	111		22(1,F3)	[redacted]			<input type="checkbox"/>	<input type="checkbox"/>		

HLA Fusion™
The Rxn of LABScreen Supplemental and SA beads will be recomputed.
Please note that the combine Rxn pattern will not be saved.
OK

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Side by Side Compare for Sample: SH3102

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Analyze Data Reports Data Sample Patient Info Profile Utilities Help Exit

- Update Reference
- Catalog Template Association
- Molecular Product Configuration
- Antibody Product Configuration
 - Create Combined Products
 - Set Default Negative Value
 - Set Mixed Product Configuration
 - Set Analysis Configuration
 - NS File Import
- General Settings
- Products Selection
- Validation
- Switch User
- Switch Database

Analysis Configuration Settings

HLA Fusion™

Antibody Product: MM

Product Type: LABScreen Single Antigen

Threshold: X Baseline

Low Bead Count *: 50 GREG *: OLI

Low NC Raw *: Low NC % *: Low PC *: Low PC/NC Ratio

Negative Sample Criteria

Normal Value of First Bead < 500

Normal Value of Third Bead < 300

Do Not Display Warning Messages

Swap Bead if NC is higher than all Beads

Swap Bead if NC is high

Swap Bead if Low NC

Set Max

Baseline: 10000 Ratio: 5 Raw: 10000

User Defined CutOff:

Exclude Patient Typing

Show DNA Mode

Use All Beads to Calculate Mean in Epitope Analysis

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
Epitopeの表示

The screenshot shows the Epitope software interface. At the top, there are checkboxes for 'DQA/DPA' and 'DNA', and a dropdown menu for 'Epitope - HMM Class2 0817'. Below this, there are radio buttons for 'All', 'Ab Verified', and 'Other', and a dropdown for 'iCr3D Latest Version'. There are also checkboxes for 'DNA Mode' and 'Assign epitope when it is double clicked'. The main part of the interface is a table with columns for 'Group Name', 'AA Position', 'Polymorphic Residues', and several columns of HLA alleles. The alleles are color-coded: blue (e.g., DRB1*01:01), pink (e.g., DRB1*01:02), orange (e.g., DRB1*14:02), yellow (e.g., DRB1*04:05), and green (e.g., DRB1*04:03). The table also includes a 'D' column on the right side.

Group Name	AA Position	Polymorphic Residues											
70QRA	70QRA	70Q71R73A	DRB1*01:01	DRB1*01:02	DRB1*14:02	DRB1*04:04	DRB1*04:05	DRB1*04:03	DRB1*01:31	DRB1*01:32	DRB1*01:34	D	
96Y2	96Y2	96Y98E120N(180L181T)	DRB1*04:01	DRB1*04:04	DRB1*04:05	DRB1*04:03	DRB1*04:02	DRB1*04:10	DRB1*04:06	DRB1*04:92	DRB1*04:34	D	
13FEL	13FEL	13F14E26L	DRB1*01:01	DRB1*01:02	DRB1*10:01	DRB1*01:03	DRB1*01:28	DRB1*01:29	DRB1*01:30	DRB1*01:31	DRB1*01:32	D	
96ES2	96ES2	96E98K120S(28E30C31)	DRB1*01:01	DRB1*01:02	DRB1*01:03	DRB1*01:04	DRB1*12:17	DRB1*01:13	DRB1*01:20				
57A	57A	57A58A60H	DRB1*14:01	DRB1*14:54	DRB1*14:120	DRB1*14:122	DRB1*14:124	DRB1*14:125	DRB1*14:16	DRB1*14:26	DRB1*14:28	D	
111R	111R	111R	DPA1*02:01	DPA1*02:02									
74SR3	74SR3	71A75S77R(116I)(125S)	DQB1*05:01	DQB1*05:02	DQB1*05:03	DQB1*05:04	DQB1*05:05	DQB1*05:06	DQB1*05:07	DQB1*05:08	DQB1*05:09	D	
55EA	55EA	55E56A57E	DPB1*05:01	DPB1*19:01	DPB1*47:01	DPB1*38:01	DPB1*36:01	DPB1*30:01	DPB1*24:01	DPB1*22:01	DPB1*21:01	D	
31Q	31Q	31Q	DPA1*02:01	DPA1*02:02	DPA1*01:06	DPA1*02:04							
125SQ	125SQ	125S126Q	DQB1*05:01	DQB1*05:02	DQB1*05:03	DQB1*05:04	DQB1*05:05	DQB1*05:06	DQB1*05:07	DQB1*05:08	DQB1*05:09	D	
178M	178M												
56PS	55 57												
112V	112V												

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