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
HLA Fusion™

Database Tasks

- Create Fusion Database
- Select/Connect to Fusion Database
- Delete Database
- Attach Database

Create Fusion Database

SQL Server

SQL Server : (local)\#FUSION_SQL14EXP

Settings : us_english#SQL_Latin1_General_CP1_CI_AS
SQL Express Edition

Database

Database Name : FUSION4.4_230119

File Size (GB) : 10.0

Fusion will create 4 files each with the size specified above

Create Help Exit

Message

Please note that some database functions such as Create, Delete, Backup, Restore, Attach, Detach etc. requires administrator level privilege. Please see your systems administrator if you need to perform these actions. The maximum allowed database name is 35 characters.

Database Connection

Database connected successfully.

OK

Version: 4.4.0.13925

User Name: FUSION_USER Server Name: (local)\#FUSION_SQL14EXP Database Name: FUSION44_220613 Database Mode: Server/Standalone

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Password*:

Confirm Password*:

What is the name of your best childhood friend? ②

What is the name of the city in which you were born?

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Contact Email:

Next >>

Note: The Lab information will be printed on HLA Fusion™ reports.

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The screenshot shows the 'HLA Fusion - Security Login' window. It has a title bar with the application name and standard window controls. The main area has a blue header with 'HLA Fusion™'. Below this are fields for 'User Name*' and 'Password*', both with red rectangular callouts. There are links for 'Forgot User Name' and 'Forgot Password'. A 'Log In' button is highlighted with a red rectangular callout and an arrow. Below the login fields, it shows 'SQL Server: 2014 Express Edition' and 'Database: (local)\FUSION_SQL14EXP\FUSIONNaa' with version and size information. At the bottom, it shows 'Regional Settings: Client: 日本語 (日本)'.

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