

# ViewRNA™ ISH Tissue Assay

## USER GUIDE



**Manufacturer:** Affymetrix Inc. | 3450 Central Expressway | Santa Clara, CA 95051 | USA

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision	Date	Description
A.0	January 2019	New document

**Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

**TRADEMARKS:** All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2019 Thermo Fisher Scientific Inc. All rights reserved.

# Contents

General Information .....	1
About this manual.....	1
How the ViewRNA ISH Tissue Assay works.....	1
Performance Highlights .....	2
Product description.....	2
Safety Warnings and Precautions.....	3
Required equipment and materials not provided.....	3
Microscopy and Imaging Equipment Guidelines .....	4
Contents and storage.....	6
ViewRNA™ ISH Tissue Assay, Core Kit.....	6
ViewRNA ISH Tissue Assay, Blue Module .....	6
Before you begin.....	8
Tissue Preparation Guidelines.....	8
Guidelines for Working with Tissue Microarrays .....	9
Experiment Design Guidelines .....	10
Sample Pretreatment Optimization .....	10
Assessment of Endogenous Alkaline Phosphatase .....	11
Probe Set Considerations.....	12
Assigning Colors to Target mRNAs in 1 vs. 2-plex Assays .....	12
Assay Procedure.....	15
About the ViewRNA™ ISH Tissue Assay Procedure .....	15
Part 1: Sample Preparation and Target Probe Hybridization.....	16
Part 2: Signal Amplification and Detection .....	21
Troubleshooting .....	27
Weak or No Signals.....	27
High Background.....	29
Diffused Signals .....	30
Endogenous Alkaline Phosphatase Activity .....	31
Tissue Detachment from Slide.....	31

Poor Cell Morphology .....	32
High Non-Specific Binding on Glass Slide .....	32
Pink Non-Specific Background Where Paraffin Was.....	33
Hydrophobic Barrier Falls Off .....	33
Fast Red Signal for TYPE 1 Target is Weak or Different in 2-Plex vs. 1-Pex .....	34
TYPE 1 Target Signals Observed in the Channel for TYPE 6 Target.....	34
Co-localized Fast Blue and Fast Red Signals When Using Only TYPE 6 Probe Set in a 2-plex Assay .....	34
Appendix A: Sample Pretreatment Optimization Procedures.....	35
About Pretreatment Optimization.....	35
Sample Pretreatment Optimization Setup.....	35
Sample Preparation and Target Probe Hybridization .....	36
Appendix B: Sample Pretreatment Lookup Table.....	39
Appendix C: Evaluating Results .....	42
Assessing Pretreatment Conditions.....	42
Analyzing Target Expression .....	43
Appendix D: Using Frozen Tissues with ViewRNA ISH Tissue Assay .....	45
About this Appendix.....	45
Important Procedural Notes.....	45
Modifications to Part 1: Sample Preparation and Target Probe Hybridization for Frozen Tissues.....	45

# General Information


## About this manual

This manual provides complete instructions for performing the ViewRNA™ ISH Tissue Assay for visualization of one or two target RNAs in formalin-fixed paraffin-embedded (FFPE) samples prepared in accordance with the guidelines provided. Appendix D provides a modified protocol for OCT-embedded frozen tissue sections.

## How the ViewRNA ISH Tissue Assay works

*In situ* hybridization (ISH) techniques are used to visualize DNA or localize RNAs within cells. Whereas *in situ* analysis of RNA has traditionally been limited by low sensitivity and complicated probe synthesis, the ViewRNA ISH Tissue Assay, based on highly specific branched DNA signal amplification technology, provides robust *in situ* detection of one or two target mRNAs within FFPE tissue sections with single-copy sensitivity.

**Figure 1: ViewRNA Assay workflow**



**Sample Preparation.** FFPE tissue sections are deparaffinized and pretreated to allow unmasking of RNA and probe accessibility.

**Target Hybridization.** Target specific probe pairs hybridize to the target RNA. A typical mRNA probe set contains 20 oligonucleotide pairs. For simplicity, only one pair per mRNA target is shown in Figure 1. TYPE 1 and TYPE 6 probe sets are designed to generate red and blue signals, respectively. These separate yet compatible signal amplification systems provide the assay with multiplex capability.

**Signal Amplification and Detection.** Signal amplification is achieved via a series of sequential hybridization steps. Pre-amplifiers hybridize to their respective pair of bound probe set oligonucleotides, then multiple amplifiers hybridize to their respective preamplifier. Next, TYPE-specific label probe oligonucleotides, conjugated to alkaline phosphatase, are sequentially hybridized to their corresponding amplifier molecules to provide up to 3,000-fold amplification per target RNA.

**Visualization.** Sequential hybridization of TYPE 6 label probe followed by addition of Fast Blue substrate and TYPE 1 label probe followed by addition of Fast Red substrate, produces blue and red precipitates (dots) respectively. The target mRNAs are visualized using a standard brightfield or fluorescent microscope.

## Performance Highlights

Specification	Description
Sample types	Formalin-fixed paraffin-embedded (FFPE) tissue section or microarray; OCT-embedded frozen tissue sections
Sensitivity	Single RNA molecule (one dot = one RNA molecule)
Targets	Detection of two target RNAs
Detection	Chromogenic and fluorescence
Nuclear stain	Hematoxylin and/or DAPI
Instrumentation	Brightfield and/or fluorescence microscope or scanner

## Product description

The ViewRNA™ ISH Tissue Assay consists of 3 modules, each sold separately:

- **ViewRNA ISH Tissue Assay, Core Kit:** contains the reagents required for the detection of a single mRNA target using the Fast Red substrate.
- **ViewRNA ISH Tissue Assay, Blue Module (for 2-plex assays only):** contains reagents for the detection of a second RNA target using the Fast Blue Substrate. Both the Blue Module and the Core Kit are required to perform a ViewRNA ISH Tissue 2-plex assay.
- **ViewRNA TYPE 1 and TYPE 6 Probe Sets:** contains oligonucleotides specific to your RNA target of interest. Two different amplification

system “TYPES” are used (consisting of Pre-amplifier, Amplifier, Label Probe, and Substrate). TYPE 1 probe sets and amplification systems are compatible with Fast Red, and TYPE 6 systems with Fast Blue.

## Safety Warnings and Precautions

- Formaldehyde is a poison and an irritant. Avoid contact with skin and mucous membranes. Use in a fume hood.
- Ammonium hydroxide is highly volatile. Use in a fume hood.
- Xylene is both flammable and an irritant. Avoid inhalation and contact with skin. Use in a fume hood.
- Probe Set Diluent QT, Pre-amplifier Mix QT, and Amplifier Mix QT contain formamide, a teratogen, irritant and possible carcinogen. Avoid contact with mucous membranes.
- DAPI is a possible mutagen. Avoid contact with skin and mucous membranes.
- Perform all procedural steps in a well-ventilated area at room temperature (RT) unless otherwise noted.
- Discard all reagents in accordance with local, state, and federal laws.

## Required equipment and materials not provided

Unless otherwise indicated, all materials are available through thermofisher.com. MLS: Fisher Scientific (fisherscientific.com) or other major laboratory supplier.

Materials	Source
Tissue Tek Staining Dish (clear color), 3	Thermo Fisher Scientific, QVC0502 American Master Tech Scientific, LWS20WH
Tissue Tek Clearing Agent Dish (green color), 2	American Master Tech Scientific, LWS20GR
Tissue Tek Vertical 24 Slide Rack, 1	Thermo Fisher Scientific, QVC0503
1000 mL glass beaker	MLS
Forceps	Thermo Fisher Scientific, QVC0501
Pipettes – P20, P200, P1000	MLS
Hydrophobic Barrier Pen	Thermo Fisher Scientific, QVC0500 Vector Laboratories, H4000
Mounting media <ul style="list-style-type: none"> <li>• Dako UltraMount Permanent Mounting Medium</li> <li>• HistoMount Mounting Solution (used only in conjunction with UltraMount)</li> <li>• ADVANTAGE Mounting Medium</li> </ul>	Agilent, S196430-2 Thermo Fisher Scientific, 008030 Innovex, NB300
Rectangular cover glass, 24 mm x 55 mm	Thermo Fisher Scientific, QVC0506
Aluminum foil	MLS
Double-distilled water (ddH <sub>2</sub> O)	MLS
100% ethanol (200 proof)	MLS

10X PBS, pH 7.2-7.4	Thermo Fisher Scientific, QVC0508
Gill's Hematoxylin I	American Master Tech Scientific, HXGHE1LT
xylene or Histo-Clear	National Diagnostics, HS-200 or equivalent
37% formaldehyde	EMD Millipore, FX0410-1 or equivalent
27-30% ammonium hydroxide	VWR, JT-9726-5 or equivalent
DAPI (optional, for fluorescence detection)	Thermo Fisher Scientific, QVC0515
<b>Equipment</b>	<b>Source</b>
Either of the following hybridization systems: <ul style="list-style-type: none"> <li>ThermoBrite System and ThermoBrite Humidity Strips</li> <li>Tissue culture incubator with &gt;85% humidity and 0% CO<sub>2</sub> and 3 aluminum slide racks for transferring slides to incubator during hybridization</li> </ul>	Abbott, 07J91-010 (110V) & 07J68-001  MLS
ViewRNA Temperature Validation Kit	Thermo Fisher Scientific, QV0523
Water-proof remote probe thermometers, validated for 90-100°C	VWR, 46610-024
Fume hood	MLS
Fisherbrand™ Isotemp™ Advanced Stirring Hotplate or equivalent	Fisher Scientific SP88857290
Table-top microtube centrifuge	MLS
Water bath capable of maintaining 40±1°C	MLS
Vortexer	MLS
Dry incubator or oven capable of maintaining 60°C for baking slides	Thermo Fisher Scientific, QS0720 (120V) or QS0721 (220V)
Microplate shaker (optional, for washing steps)	Thermo Fisher Scientific, 88880023 or 88880024
Microscope and imaging equipment	(See next section)

## Microscopy and Imaging Equipment Guidelines

The stains used to label RNA in the ViewRNA ISH Tissue Assay can be visualized using brightfield or fluorescence microscopy.

### Stains for ViewRNA ISH Tissue Assay:

Detect	Staining Reagent	Stain Color	
		Brightfield View	Fluorescent View
RNA 1 using TYPE 1 probe	Fast Red	Red	Red
RNA 2 using TYPE 6 probe	Fast Blue	Aqua blue	Far red
Nucleus	Hematoxylin/DAPI	Light purple-blue	Blue



### ViewRNA ISH Tissue Assay Viewing and Digital Capture Options:

	Brightfield viewing	Fluorescence viewing and image capture	Automated image capture (brightfield or fluorescence)
<b>Microscope Type</b>	Standard brightfield microscope	Microscope with camera and fluorescence options. Verify that the camera does not have infrared blocking filter.	Digital pathology scanner system
<b>Recommended Microscope</b>	<ul style="list-style-type: none"> <li>Leica DM series</li> <li>Nikon E series</li> <li>Olympus BX series</li> <li>Zeiss Axio Lab/Scope/Imager</li> <li>(or equivalent)</li> </ul>		<ul style="list-style-type: none"> <li>Aperio ScanScope AT/XT/CS, use FL version for fluorescence</li> <li>Leica SCN400-F</li> <li>Olympus Nanozoomer RS</li> </ul>
<b>Required Optics</b>	Requires 20X and 40X objectives	<ul style="list-style-type: none"> <li>Requires 20 and 40X objectives</li> <li>Numerical aperture (NA) &gt; 0.5</li> </ul>	Recommend scanning at 40X when expression is low
<b>Recommended Filters</b>	Requires neutral density filters and/or color filters for white balancing	<p>For Fast Red Substrate, use Cy3/TRITC filter set: Excitation: 530 ±20 nm Emission: 590 ±20 nm Dichroic: 562 nm</p> <p>For Fast Blue Substrate, use custom filter set<sup>1</sup>: Excitation: 630 ± 20 nm Emission: 775 ± 25 nm Dichroic: 750 nm</p> <p>For DAPI filter set: Excitation: 387/11 nm Emission: 447/60 nm</p>	(see left)

[1] Recommended vendor: Semrock Cy7-B/Alexa 750 filter modified with excitation filter FF02-28/40-25

# Contents and storage

## ViewRNA™ ISH Tissue Assay, Core Kit

Each ViewRNA Tissue Assay Core Kit contains the following components, supplied in 2 boxes based on storage temperature. Reagents have a minimum shelf life of 6 months from date of receipt when stored as recommended. Please refer to kit labels for expiration dates and to the Package Insert for component quantities. Kits are configured for processing a minimum of 6 assays per experiment.

Component	Description	Storage
100x Pretreatment Solution	Aqueous buffered solution	2-8 °C
Probe Set Diluent QT	Aqueous solution containing formamide, detergent, and blocker	
Pre-Amplifier Mix QT	DNA in aqueous solution containing formamide and detergent	
Amplifier Mix QT	DNA in aqueous solution containing formamide and detergent	
Label Probe Diluent QF	Aqueous solution containing detergent	
AP Enhancer Solution	Aqueous buffered solution	
Naphthol Buffer	Buffer required for preparation of Red Substrate	
Protease QF <sup>1</sup>	Enzyme in aqueous buffered solution	
Fast Red Tablets	Red precipitating substrate for the detection of alkaline phosphatase activity	
Label Probe 1-AP <sup>1</sup>	Alkaline phosphatase-conjugated oligonucleotide in aqueous buffered saline	
Wash Buffer Component 1	Aqueous buffered solution containing detergent	15–30 °C
Wash Buffer Component 2	Aqueous buffered solution	

[1] DO NOT FREEZE!

## ViewRNA ISH Tissue Assay, Blue Module

Each ViewRNA Tissue Assay Blue Module contains the following components. Reagents have a minimum shelf life of 6 months from date of receipt when stored as recommended. Please refer to kit labels for expiration dates and to the Package Insert for component quantities.

Component	Description	Storage
Label Probe 6-AP <sup>1</sup>	Alkaline phosphatase-conjugated oligonucleotide in aqueous buffered saline	2-8 °C
Blue Reagent 1	Blue precipitating substrate component 1 for the detection of alkaline phosphatase activity	
Blue Reagent 2	Blue precipitating substrate component 2 for the detection of alkaline phosphatase activity	
Blue Reagent 3	Blue precipitating substrate component 3 for the detection of alkaline phosphatase activity	
Blue Buffer	Buffer required for preparation of Blue Substrate	
AP Stop QT	Aqueous buffered solution intended for the inactivation of residual Label Probe 6-AP activity after the Fast Blue Substrate development	

[1] DO NOT FREEZE!

# Before you begin

## Tissue Preparation Guidelines

This section provides guidelines for preparation of FFPE tissue blocks, FFPE tissue slides, and tissue microarray (TMA) slides. Samples prepared outside of these guidelines may produce inadequate results.


### FFPE/TMA Tissue Block Preparation

- Immediately place freshly dissected tissues in  $\geq 20$  volumes of fresh 10% neutral buffered formalin (NBF) or 4% paraformaldehyde (PFA) at room temperature (RT) for 16-24 hr. Trim larger specimens to  $\leq 3$  mm thickness to ensure faster diffusion of the fixative into the tissue.
- Rinse, dehydrate, and embed in a paraffin block.
- Store FFPE tissue blocks at RT.

### FFPE/TMA Tissue Slide Preparation

- Section FFPE tissue to a thickness of  $5 \pm 1$   $\mu\text{m}$ .
- If working with TMAs, core size should be  $\geq 1.0$  mm diameter.
- Maximum tissue area is 20 mm x 30 mm and should fit within the hydrophobic barrier.
- Mount sections as shown in Figure 2.1 on page 9 onto one of the following types of positively-charged glass slides:
  - Leica Non-Clipped X-tra® Slides, 1 mm white (P/N 3800200 in U.S., Canada, and Asia Pacific regions or P/N 3800210 in Europe).
  - Fisherbrand™ Superfrost™ Plus Slides, white label (Fisher Scientific, P/N 12-550-15). Avoid other colored labels as they tend to give high background.
- Air dry freshly-mounted sections at RT overnight or at 37 °C for 5 hr.
- Bake slides at 60 °C for 1 hr to immobilize tissue sections.
- Storage:
  - Short-term: Store sections in a slide box at RT for up to 2 weeks.

- Long-term: Store sections in a slide box at  $-20^{\circ}\text{C}$  for up to 1 year (avoid freeze/thaw).
- Slides can be shipped at the temperature at which they were originally stored.



## Guidelines for Working with Tissue Microarrays

Process TMA slides using the same assay procedures but with the following two modifications:

- Increase the initial baking step time from 60 to 90 min. This additional baking time will increase the tissue attachment to the slide, reducing the risk of small ( $> 1\text{ mm}$ ) core sections falling off during assay procedure.
- Increase the volume/slide of the protease working solution to prevent tissues at the edge of the TMA from drying out.

When designing TMAs to be used in the ViewRNA™ ISH Tissue Assay, it is important to understand that only one optimized condition can be used when running the assay. Therefore, if you want multiple tissue types within the same TMA block, we recommend running an optimization procedure on each individual FFPE tissue type to identify the most favorable pretreatment boiling and protease condition. Based on the optimal condition of the tissue morphology, signal strength, and residual cores, you can judge whether there may be one optimization condition that will be suitable for all of the sample types.

# Experiment Design Guidelines

## Assay Controls

We recommend running one positive and one negative control slide in each assay, based on your sample type. This will allow you to qualify and interpret your results.

### Negative Control

This slide undergoes the entire assay procedure and assesses the assay background from different levels.

The negative control can be one of the following:

- Omit the target probe set. A no probe negative control.
- Use a probe set designed to the sense strand of the target – A more target-specific negative control used to subtract assay background when assessing results.
- Use a probe set for a target not present in your tissue sample – A more general negative control used to subtract assay background when assessing results, for example, the bacterial gene *dapB*.

### Positive Control

This slide undergoes the entire assay procedure using a probe set against an ubiquitous or tissue-specific target that has consistent, medium-high to high, but not saturating, expression level. A positive control ensures that the assay procedure has been successfully run. Examples of positive control targets include:

- Housekeeping Genes: *ACTB*, *GAPD*, or *UBC*
- Housekeeping Gene Panel: A panel of several housekeeping genes can be pooled and used as a positive control whenever the expression level of any one given housekeeping gene is unknown in the tissue of interest. For example, pool *ACTB*, *GAPD* and *PPIB* probe sets at equal volumes to form a panel, and then dilute the panel of probe sets 1:40 to create a working probe set solution for use as a positive control.

### Replicates

We recommend running all assays in duplicate.

## Sample Pretreatment Optimization

The pretreatment of tissue sections is critical for the success of all in situ assays, and consists of heat treatment and protease digestion. These pretreatment steps

help to unmask the RNA targets, allowing for better probe accessibility and thereby increasing assay signal. However, excessive pretreatment can have a negative effect on tissue morphology. Thus, we recommend optimizing the pretreatment conditions for each new tissue type (see Appendix A, Sample

Pretreatment Optimization Procedures). Once the optimal pretreatment conditions are determined, they can generally be used for most targets within the particular tissue. In instances when the transcript is particularly rare or expressed at an extremely low level, the optimal pretreatment condition may need to be one that favors signal over morphology.

Refer to the Sample Pretreatment Lookup Table in Appendix B for heat treatment and protease conditions that we have found to be optimal for a number of tissues prepared according to the recommended guidelines in this manual using 10% NBF. This table serves as a reference or starting point only and may not be applicable to tissues prepared using 4% PFA. If you do not obtain the desired results, we recommend performing either the full or limited Sample Pretreatment Optimization Procedures on page 31, depending on availability of your samples.

When optimizing pretreatment conditions for TMAs, it is important to understand that it is impossible to identify one condition that is ideal for every tissue type on the array. The optimal pretreatment conditions in such a case would be one that maximizes the number of cores with assay signal and minimizes the number of cores lost due to excessive heat treatment and protease digestion. Due to their high cost and limited quantity, TMAs would greatly benefit from the limited pretreatment optimization procedure, since only as few as three slides might be necessary (see Table B.2 in Appendix B).

## Assessment of Endogenous Alkaline Phosphatase

The ViewRNA™ ISH Tissue Assay uses alkaline phosphatase to convert a chromogenic substrate into a colored signal. For this reason it is important to assess the level of endogenous alkaline phosphatase (AP) activity in your tissue of interest prior to performing the assay.

Certain types of tissue (such as stomach, intestine, placenta and mouse embryo) are known to possess high levels of endogenous AP activity that can interfere with the assay. While the problem is more prevalent in fresh frozen tissues, it has also been observed in some FFPE samples.

To empirically determine the level of endogenous AP activity in your tissue type, perform the pretreatment protocol as instructed for fresh frozen or FFPE tissue. After the protease treatment and fixation in 10% NBF, wash the samples in 1X TBS (Sigma, T5912-1L) and incubate the sections with either Fast Blue Substrate or Fast Red Substrate.

If present, endogenous AP can be inactivated with 0.2 M HCl/300 mM NaCl at RT for 15 min just before the probe hybridization but after the sample has undergone protease treatment, 10% NBF fixation and 2 washes in 1X PBS.

## Probe Set Considerations

Probe sets of the same TYPE can be combined to create a target panel or cocktail. For example, identifying epithelial cells could be easily accomplished by pooling different cytokeratin probe sets of the same type, such as TYPE 1, KRT5, KRT7, KRT8, KRT10, KRT19, KRT19 and KRT20, into a single assay. However, we do not recommend combining more than 10 targets for any one signal amplification system, be it TYPE 1 or TYPE 6.

How the probe sets are diluted to generate a panel depends on the application. For example, if the goal is to identify all of the epithelial cells or to assess RNA integrity, then each probe set can be diluted 1:40. However, when using a panel of housekeeping gene probe sets for optimizing pretreatment conditions, the probe sets (e.g., ACTB, GAPD and PPIB) should be pooled at equal volumes to form the panel, and then diluted 1:40 to create the working probe set solution. This ensures that the panel expression is sufficiently high but not saturated so that the differences in signal between pretreatment conditions can be distinguished.

The typical design for a ViewRNA Probe Set consists of 40 unlabeled oligos, or 20 pairs of oligos per RNA target, and spans approximately 1000 bases of the target transcript to achieve maximal sensitivity. The binding of these oligo pairs side-by-side to the target sequence serves as a base upon which the signal amplification is built, and is the core of the assay's sensitivity and specificity. Using multiple pairs of oligos in a single probe set ensures that there are many opportunities for the probe to bind to the target's unmasked/accessible regions so as to achieve the maximal signal amplification possible for that particular RNA target molecule. When working with smaller targets or applications such as splice variants or RNA fusions, the available number of oligo pairs in the probe set is naturally reduced, and this will directly impact the sensitivity of the assay. That is, the probes will have fewer opportunities to find the unmasked areas of the target in order to generate signal at that location. In these cases, increasing the probe set concentration used in the assay from 1:40 to 1:30 or 1:20 might increase the sensitivity. However, note that there is always a general trade-off between sensitivity and specificity.

## Assigning Colors to Target mRNAs in 1 vs. 2-plex Assays

The ViewRNA™ ISH Tissue Assay allows in situ detection of up to two mRNA targets simultaneously, using the ViewRNA TYPE 1 and/or TYPE 6 probe sets. The standard workflow of the assay is designed to automatically assign Fast Red



signal to TYPE 1 and Fast Blue signal to TYPE 6 probe sets. While both the Fast Red and Fast Blue signals that form are easily visible under brightfield, the red dots generally have a much higher contrast than the blue dots, especially in the presence of hematoxylin. Thus, when the detection of only one target (1-plex assay) is desired, we recommend using either TYPE 1 or TYPE 6 probe set and developing the signal as Fast Red. The Core Kit (SKU 19931) is configured for use with TYPE 1 probe sets.

When performing a 2-plex assay, we recommend assigning the TYPE 1 probe set (Fast Red) to the more important target of the two. Reserve the TYPE 6 probe set (Fast Blue) for the less critical target, such as a housekeeping gene. Due to the nature of the chromogenic assay and the sequential development of Fast Blue before Fast Red signals, large quantities of blue precipitate that are deposited, particularly when a TYPE 6 target is expressed homogeneously at high level, have the potential to partially block subsequent hybridization of the TYPE 1 Label Probe and consequently the development of the Fast Red signal. For this reason, the target assigned to Fast Blue should preferably have lower expression than the one assigned to Fast Red to ensure against potential interference with Fast Red signal development downstream.

If only medium and high expressing housekeeping targets are available in a particular tissue type and the critical target of interest has low to medium expression, a 2-plex assay can still be performed by assigning Fast Red to the housekeeping target and Fast Blue to the second target. Brightfield detection of the Fast Blue signal for a medium expressing transcript could still be easily done, while fluorescent detection would provide a more sensitive alternative for detecting a low expressing target tagged with Fast Blue.

### **Fluorescent Mode Guidelines**

The advantage of using alkaline phosphatase-conjugated label probe for the enzymatic signal amplification is the availability of substrates with dual property, such as Fast Red and Fast Blue, which allows for both chromogenic and fluorescent detection of the targets. However, for a 2-plex assay in which both Label Probe 1 and Label Probe 6 are conjugated to the same alkaline phosphatase, the enzymes conjugates are unable to differentiate between Fast Red and Fast Blue if both substrates are added simultaneously. As a result, the enzymatic signal amplification has to be performed sequentially in order to direct substrate/color specificity to each target. Additionally, complete inactivation of the first alkaline phosphatase-conjugated label probe (LP6-AP) is necessary, especially when employing fluorescence mode for the detection of the targets. Otherwise, the residual LP6-AP activity can also convert Fast Red substrate in subsequent step into a red signal even at locations where TYPE 1 target is not present, giving a false impression that the Fast Blue and Fast Red signals are co-localized. For this reason, it is absolutely necessary to quench any residual LP6-AP activity with the ViewRNA AP Stop QT prior to proceeding

with the second label probe hybridization and development of the Fast Red color as this will ensure specific signals in fluorescent mode and brighter aqua blue dots in chromogenic mode.

Fast Red has a very broad emission spectrum and its bright signal that can bleed into adjacent Cy5 channel if one uses the standard Cy3/Cy5 filter sets for imaging. For this reason, it is critical that the recommended filter set for Fast Blue detection be used to avoid spectral bleed through of the Fast Red signal into the Fast Blue channel and interfering with Fast Blue detection. Please refer to the Viewing and Digital Capture Options table in the Microscopy and Imaging Equipment Guidelines for exact filter set specifications.

### **Limitations of Chromogenic In Situ Assay in Co-localization Studies**

When employing the ViewRNA™ ISH Tissue 2-Plex Assay for co-localization studies, it is crucial to understand the assay's strengths and limitations. By definition, a requisite for in situ detection is target accessibility. While the assay, with its branched DNA technology, has the capability to detect RNA molecules down to single-copy sensitivity and the probe sets are designed to maximize the binding opportunities to all accessible regions of the targets, the overall detection for any given target is only as good as the unmasking of the target site is able to provide. This essentially means that in situ assays in general are only capable of relative and not absolute detection. That is, not every single molecule of a given target can be detected. So in practice, even if two RNA targets are theoretically expected to be colocalized, only a subset these two transcripts will be detected as being so due to lack of complete target accessibility.

Another factor that can limit the use of this assay for co-localization studies is the nature of chromogenic assay and the sequential development of Fast Blue then Fast Red signals. In chromogenic assay, the enzyme converts the substrate into color precipitates and deposits them at the site where the RNA molecule is localized. Because the Fast Blue and Fast Red substrates are sequentially developed in the ViewRNA ISH Tissue 2-Plex Assay, the Fast Blue precipitates that are formed first and deposited have the potential to partially block subsequent hybridization of the TYPE 1 Label Probe, by masking its binding sites on a nearby/co-localized target and consequently affecting the development of the Fast Red signal. This is yet another form of accessibility issue that needs to be considered when performing colocalization studies and analyzing the data obtained from such studies. Consequently, even when two targets are co-localized, only a subpopulation of the two is actually observed as such because of target accessibility, be it at the probe hybridization step due to incomplete unmasking or at the label probe hybridization step due to masking of the binding site by the deposition of the Fast Blue precipitates.

# Assay Procedure

## About the ViewRNA™ ISH Tissue Assay Procedure

The ViewRNA ISH Tissue Assay can be run in a single long day or broken up over two days for added flexibility. The procedure includes two parts:

- Part 1: Sample Preparation and Target Probe Set Hybridization (optional stopping point)
- Part 2: Signal Amplification and Detection

### Important Procedural Notes and Guidelines

- The procedure assumes running a maximum of 12 slides at a time and that the size of the section does not exceed the maximum coverage area recommended.
- Do not mix and match kit components from different lots.
- Before beginning the procedure, know the optimized conditions (heat treatment time and protease digestion time) for your sample type. If you do not know these optimized conditions, refer to Appendix A: Sample Pretreatment Optimization Procedures.
- Throughout the procedure, dedicate the Tissue Tek staining dishes as follows:
  - Clear staining dish for Formaldehyde.
  - Green staining dish for Gill's hematoxylin.
  - Green staining dish for xylene/Histo-Clear.
  - The remaining two clear staining dishes can be used interchangeably for 1X PBS, 100% ethanol, Wash Buffer, ddH<sub>2</sub>O, Storage Buffer, and DAPI. Rinse staining dishes between steps with ddH<sub>2</sub>O.
- If using a humidified tissue culture incubator (without CO<sub>2</sub>) as the hybridization system:
  - Verify that the water jacket or bottom tray is filled with water.
  - Use an aluminum slide rack to transfer slides to the incubator.
  - Do not leave the incubator door open longer than necessary when transferring slides, particularly during the protease optimization procedure. This will help maintain the required temperature.
- Typical processing times included in the assay procedure assume that the preparations for the following step are being done during the incubation periods.

### Essential Keys for a Successful Assay

- Prepare samples following the Tissue Preparation Guidelines under “Before You Begin”.
- Organize the preparation of the assay before you start:
  - Verify that all materials and equipment are available.
  - Be mindful of the incubation times/temperatures, as variations can negatively affect assay signal and background.
  - Double-check all reagent calculations, as correct reagent volumes and concentrations are critical.
- Employ good washing techniques. Frequently, washing is performed too gently. Adequate washing is important for consistent low backgrounds.
- Calibrate temperatures for hybridization system (to 40°C) and dry oven (to 60°C) using the ViewRNA™ Temperature Validation Kit.
- Ensure that hybridization system is appropriately humidified.
- DO NOT let tissues dry out where indicated in the procedure.
- Incorporate controls, both positive and negative, so that results are unambiguous and can be interpreted. See Experiment Design Guidelines under “Before You Begin”.

## Part 1: Sample Preparation and Target Probe Hybridization

1. Bake slides (65 minutes)
  - a. Set the dry oven or hybridization system to 60±1°C.
  - b. Label the slides with a pencil
  - c. Bake the slides. With a dry oven, insert slides into the slide rack and bake for 60 minutes. With a ThermoBrite System, keep the lid open and bake for 60 minutes. Make sure that the temperature of the ThermoBrite is validated with the lid open.
2. Prepare buffers, reagents, and equipment while slides bake
  - a. Verify that the hybridization system is set to 40± 1°C and that it is appropriately humidified.
  - b. Prepare 3 L 1X PBS: add 300 mL 10X PBS and 2.7 L ddH<sub>2</sub>O to a 3 L capacity container.
  - c. Prepare 200 mL 10% NBF, working in fume hood: add 178 mL 1X PBS + 22 mL 37% formaldehyde to a 200 mL capacity container and mix well.
  - d. Prepare 4 L Wash Buffer: add the components below in the order listed to a 4 L capacity container and mix well:
    - 3 L ddH<sub>2</sub>O
    - 36 mL Wash Comp 1
    - 10 mL Wash Comp 2Adjust the total volume to 4 L with ddH<sub>2</sub>O.
  - e. Prepare 500 mL 1X Pretreatment Solution: add 5 mL 100X Pretreatment Solution and 495 mL ddH<sub>2</sub>O to a 1 L glass beaker.

- f. If using the optional stopping point, prepare 200 mL Storage Buffer: add 60 mL Wash Comp 2 and 140 mL ddH<sub>2</sub>O to a 200 mL capacity container.
  - g. Prepare 1 L of 0.01% ammonium hydroxide, working in a fume hood: add 0.33 mL 30% ammonium hydroxide and 999.67 mL ddH<sub>2</sub>O in a 1 L capacity container.
  - h. Ensure the availability of:
    - 600 mL 100% ethanol
    - 1.4 L ddH<sub>2</sub>O
    - 600 mL xylene or 400 mL Histo-Clear
    - 200 mL Gill's Hematoxylin I
    - 200 mL of 3 µg/mL DAPI in 1X PBS (optional, for fluorescence detection). Store in the dark at 4 °C until use.
  - i. Thaw probe set(s). Mix, briefly centrifuge to collect contents, and place on ice until use.
  - j. Pre-warm 40 mL 1X PBS and Probe Set Diluent QT to 40±1°C.
  - k. If performing both parts of the assay in 1 day:
    - Pre-warm Pre-amplifier Mix QT, Amplifier Mix QT, and Label Probe Diluent QF to 40°C .
    - Briefly spin down the Label Probe 1-AP, Label probe 6-AP, and Blue reagents (for 2-plex), then place on ice until use.
    - Bring Fast Red Tablets, Napthol Buffer, Blue Buffer (for 2-plex), and AP Enhancer Solution to RT.
  - l. Optional: if using a microplate shaker for the washes, set the speed to 550 rpm. Place a slide rack in a clear staining dish containing the appropriate reagent and insert the slides into the rack. Manually lift the rack up and down 10 times. Put the lid on the staining dish and place it on a microplate shaker platform that is equipped with a non-skid pad. Shake for the recommended amount of time.
3. Deparaffinization (30 minutes)
- If using xylene (work in a fume hood):**
- a. Pour 200 mL of xylene into a green clearing agent dish.
  - b. Transfer the rack of baked slides to the green clearing dish containing the xylene.
  - c. Incubate the slides at RT for 5 min. Agitate frequently by moving the rack up and down.
  - d. Discard the used xylene and refill with another 200 mL of fresh xylene. Incubate slides at RT for 5 min with frequent agitation.
  - e. Repeat Step D above.
  - f. Remove the slide rack from the xylene and wash the slides twice, each time with 200 mL of 100% ethanol for 5 min with frequent agitation.
  - g. Remove the slides from the rack and place them face up on a paper towel to air dry for 5 min at RT.

**If using Histo-Clear:**

- a. Pour 200 mL of Histo-Clear into a green clearing dish and insert an empty slide rack.
  - b. Set the dry oven or hybridization system to  $80\pm 1^{\circ}\text{C}$ .
  - c. Bake the slide for 3 min to melt the paraffin.
  - d. Immediately insert the warm slides into the Histo-Clear and agitate frequently by moving the rack up and down for 5 min at RT.
  - e. Discard the used Histo-Clear and refill the dish with another 200 mL of fresh Histo-Clear. Agitate frequently by moving the rack up and down for another 5 min at RT.
  - f. Remove the slide rack from the Histo-Clear and wash the slides twice, each time with 200 mL of 100% ethanol for 5 min with frequent agitation.
  - g. Remove the slides from the rack and place them face up on a paper towel to air dry at RT for 5 min.
4. Draw hydrophobic barrier (40 minutes)
- a. Dab the hydrophobic barrier pen on a paper towel several times before use to ensure proper flow of the hydrophobic solution.
  - b. To create a hydrophobic barrier:
    - Place the slide over the template image, making sure that the tissue sections fall inside the blue rectangle.
    - Lightly trace the thick blue rectangle 2-4 times with the hydrophobic barrier pen to ensure a solid seal.
    - Allow for barrier to dry at RT for 20-30 min. Begin the next step while the barrier is drying.



5. Heat pretreatment (10-25 minutes, depending on optimized time)
- a. Tightly cover the beaker containing the 500 mL of 1X Pretreatment Solution with aluminum foil, place it on a hot plate, and heat the solution to a temperature of  $90-95^{\circ}\text{C}$ . Use a waterproof probe thermometer to measure and maintain the temperature of the solution at  $90-95^{\circ}\text{C}$  during the pretreatment period.
  - b. Load the slides into the vertical slide rack.
  - c. Using a pair of forceps, submerge the slide rack into the heated 1X Pretreatment Solution. Cover the glass beaker with aluminum foil and incubate at  $90-95^{\circ}\text{C}$  for the optimal time as determined in Appendix A: Sample Pretreatment Optimization Procedure.
  - d. After pretreatment, remove the slide rack with forceps, submerge it into a clear staining dish containing 200 mL of ddH<sub>2</sub>O, and wash for 1 min with frequent agitation.

- e. Repeat the wash one more time with 200 mL of fresh ddH<sub>2</sub>O.
- f. Transfer the slide rack to a clear staining dish containing 1X PBS.

**IMPORTANT:** Do not let the tissue sections dry out from this point forward. After heat pretreatment, sections can be stored covered in 1X PBS at RT overnight.

6. Protease digestion and fixation (30-50 minutes, depending on optimized time)
  - a. Prepare the working protease solution by diluting the Protease QF 1:100 in prewarmed 1X PBS (e.g. 4µL Protease QF added to 396 µL 1X PBS prewarmed to 40°C) and briefly vortex to mix. Scale reagents according to the number of assays to be run. Include one slide volume overage.
  - b. Remove each slide and flick it to remove excess 1X PBS. Without completely drying out the sections, tap the slides on the edge and then wipe the backside on a laboratory wipe.
  - c. Place the slides face up on a flat, elevated platform (e.g., Eppendorf tube rack for easier handling) and immediately add 400 µL of the working protease solution onto the tissue section. Make sure that the tissue section is covered with working protease solution. It may be necessary to spread the solution with a pipette tip.
  - d. Transfer the slides to the hybridization system and incubate at 40 °C for the optimal time as determined in Appendix A: Sample Pretreatment Optimization Procedure.
  - e. Pour 200 mL of 1X PBS into a clear staining dish and insert an empty slide rack into the dish.
  - f. After the incubation, decant the working protease solution from the slides, insert the slides into the rack and wash gently but thoroughly by moving the rack up and down for 1 min.
  - g. Repeat the wash one more time with another 200 mL of fresh 1X PBS.
  - h. Transfer the slide rack to a clear staining dish containing 200 mL of 10% NBF and fix for 5 min at RT under a fume hood.
  - i. Wash the slides twice, each time with 200 mL of fresh 1X PBS for 1 min with frequent agitation.
7. Target probe set hybridization (2 hours and 10 minutes)
  - a. Prepare the working probe set solution using the table below as a guide. Dilute the ViewRNA Probe Set 1:40 in prewarmed Probe Set Diluent QT and briefly vortex to mix. Scale reagents according to the number of assays to be run and include one slide volume overage.  
**NOTE:** Add only 400 µL of Probe Set Diluent QT to the "negative control" or probe negative control slide.

1-plex (400 µL total volume)	Volume
Probe Set Diluent QT (prewarmed to 40°C)	390 µL
ViewRNA TYPE 1	10 µL

2-plex (400 $\mu$ L total volume)	Volume
Probe Set Diluent QT (prewarmed to 40°C)	380 $\mu$ L
ViewRNA TYPE 1 Probe Set	10 $\mu$ L
ViewRNA TYPE 6 Probe Set	10 $\mu$ L

- b. Remove each slide and flick it to remove excess 1X PBS. Without completely drying out the sections, tap the slides on the edge and then wipe the backside on a laboratory wipe
  - c. Place the slides face up on a flat, elevated platform and immediately add 400  $\mu$ L of pre-warmed Probe Set Diluent QT to the negative probe control and 400  $\mu$ L of working probe set solution to each test sample.
  - d. Transfer the slides to the hybridization system and incubate at 40°C for 2 hr.
8. Wash slides (8 minutes)
  - a. Insert an empty slide rack into a clear staining dish containing 200 mL of Wash Buffer.
  - b. After incubation, decant the working probe set solution from the slides and insert them into the slide rack.
  - c. Wash the slides 3 times, each time with 200 mL of fresh Wash Buffer at RT for 2 min with constant and vigorous agitation.
  - d. If you plan to perform the assay over the course of two days, proceed to Step 9, "Optional Stopping Point". Otherwise, proceed to Part 2 Step 12, "Pre-amplifier Hybridization" to complete the entire assay in one day.
9. Optional stopping point (1 minute)
  - a. Store slides in a clear staining dish containing 200 mL of Storage Buffer at RT for up to 24 hr. Cover the dish with a lid or sealing film to prevent evaporation.
  - b. Discard 1X Pretreatment Solution, 10% NBF, remaining protease and probe set working solutions.
  - c. Store the remaining 1X PBS and Wash Buffer at RT for use in Part 2: Signal Amplification and Detection.
  - d. If using a ThermoBrite System, rewet the ThermoBrite Humidity Strips in ddH<sub>2</sub>O.
  - e. Proceed to Step 10, "Prepare Additional Buffers and Reagents" when you are ready to continue the assay.



## Part 2: Signal Amplification and Detection

10. Prepare additional buffers and reagents (5 minutes)
  - a. Pour Gill's Hematoxylin into a clear staining dish and store at RT protected from light until use.
  - b. Pre-warm Pre-amplifier Mix QT, Amplifier Mix QT, and Label Probe Diluent QF buffers to 40°C.
  - c. Briefly spin down the Label Probe 1-AP, Label probe 6-AP (for 2-plex), and Blue reagents (for 2-plex). Place them on ice.
  - d. Bring Fast Red Tablets, Naphthol Buffer, AP Enhancer Solution, and Blue Buffer (for 2-plex) to RT.
11. Wash slides (5 minutes)
  - a. Remove the slides from Storage Buffer.
  - b. Wash the slides 2 times, each time with 200 mL of fresh Wash Buffer at RT for 2 min with constant and vigorous agitation.
12. Pre-amplifier hybridization (35 minutes)
  - a. Swirl the Pre-amplifier Mix QT bottle briefly to mix the solution.
  - b. Remove each slide and flick it to remove the Wash Buffer. Without completely drying out the sections, tap the slide on its edge and then wipe the backside on a laboratory wipe. Place the slides face up on a flat, elevated platform and immediately add 400 µL of Pre-amplifier Mix QT to each tissue section.
  - c. Transfer slides to the hybridization system and incubate at 40°C for 25 min.
13. Wash slides (8 minutes)
  - a. Insert an empty slide rack into a clear staining dish containing 200 mL of Wash Buffer.
  - b. After incubation, decant the Pre-amplifier Mix QT from the slides and insert them into the slide rack.
  - c. Wash the slides 3 times, each time with 200 mL of fresh Wash Buffer at RT for 2 min with constant and vigorous agitation.
14. Amplifier hybridization (20 minutes)
  - a. Swirl the Amplifier Mix QT bottle briefly to mix the solution.
  - b. Remove each slide and flick it to remove the Wash Buffer. Without completely drying out the sections, tap the slide on its edge and then wipe the backside on a laboratory wipe. Place the slides face up on a flat, elevated platform and immediately add 400 µL of Amplifier Mix QT to each tissue section.
  - c. Transfer slides to the hybridization system and incubate at 40°C for 15 min.

15. Wash slides (8 minutes)
- Insert an empty slide rack into a clear staining dish containing 200 mL of Wash Buffer.
  - After incubation, decant the Amplifier Mix QT from the slides and insert them into the slide rack.
  - Wash the slides 3 times, each time with 200 mL of fresh Wash Buffer at RT for 2 min with constant and vigorous agitation.

**IMPORTANT: For 1-plex assays, proceed directly to step 21: Label Probe 1-AP Hybridization, and omit steps 16-20. For 2-plex assays, continue with step 16.**

16. Label Probe 6-AP hybridization (20 minutes)
- Briefly vortex and spin down Label Probe 6-AP before using.
  - Prepare the Working Label Probe 6-AP Solution using the table below as a guide. Dilute Label Probe 6-AP 1:1000 in pre-warmed Label Probe Diluent QF and briefly vortex to mix. Scale reagents according to the number of assays to be run and include overage of one slide's volume.

Working Label Probe 6-AP Solution per slide (400 $\mu$ L total volume)	Volume
Label Probe Diluent QF (pre-warmed to 40°C)	399.6 $\mu$ L
Label Probe 6-AP	0.4 $\mu$ L

- Remove each slide and flick to remove the Wash Buffer. Without completely drying out the sections, tap the slide on its edge and then wipe the backside on a laboratory wipe. Place the slides face up on a flat, elevated platform and immediately add 400  $\mu$ L of Working Label Probe 6-AP Solution to each tissue section.
  - Transfer the slides to the hybridization system and incubate at 40°C for 15 min.
17. Wash slides (12 minutes)
- Insert an empty slide rack into a clear staining dish containing 200 mL of Wash Buffer.
  - After incubation, decant the Working Label Probe 6-AP Solution from the slides and insert them into the slide rack.
  - Wash the slides 3 times, each time with 200 mL of fresh Wash Buffer at RT for 3 min with constant and vigorous agitation.
18. Apply Fast Blue Substrate (35 minutes)
- Prepare the Fast Blue Substrate: Add 5 mL of Blue Buffer and 105  $\mu$ L of Blue Reagent 1 to a 15 mL conical tube and vortex. Add 105  $\mu$ L of Blue Reagent 2 and vortex. Add 105  $\mu$ L Blue Reagent 3 and briefly vortex. Protect from light by wrapping in aluminum foil until use.
  - Remove each slide and flick it to remove the Wash Buffer. Without completely drying out the sections, tap the slide on its edge and then wipe the backside on a laboratory wipe. Place the slides face up on a flat, elevated platform and immediately add 400  $\mu$ L of Fast Blue Substrate.

- c. Transfer the slides to the hybridization system and incubate in the dark at RT for 30 min.
19. Wash slides (12 minutes)
    - a. Insert an empty slide rack into a clear staining dish containing 200 mL of Wash Buffer.
    - b. After incubation, decant the Fast Blue Substrate from the slides and insert them into the slide rack.
    - c. Wash the slides 3 times, each time with 200 mL of fresh Wash Buffer at RT for 3 min with constant and vigorous agitation.
  20. Quench Label Probe 6-AP (35 minutes)
    - a. Remove each slide and flick it to remove the Wash Buffer. Without completely drying out the sections, tap the slide on its edge and then wipe the backside on a laboratory wipe. Place the slides face up on a flat, elevated platform and immediately add 400  $\mu$ L of AP Stop QT. Incubate in the dark at RT for 30 min.
    - b. Insert an empty slide rack into a clear staining dish containing 200 mL of 1X PBS.
    - c. After incubation, decant the AP Stop Buffer from the slides and insert them into the slide rack.
    - d. Wash the slides twice, each time in 200 mL of fresh 1X PBS at RT for 1 min with frequent agitation.
    - e. Replace the 1X PBS with 200 mL of fresh Wash Buffer and rinse any residual PBS from the slides by moving the slide rack up and down for 1 min.
  21. Label Probe 1-AP hybridization (20 minutes)
    - a. Briefly vortex and spin down Label Probe 1-AP before using.
    - b. Prepare the Working Label Probe 1-AP Solution using the table below as a guide. Dilute Label Probe 1-AP 1:1000 in pre-warmed Label Probe Diluent QF and briefly vortex to mix. Scale reagents according to the number of assays to be run and include overage of one slide's volume.

<b>Working Label Probe 6-AP Solution per slide (400 <math>\mu</math>L total volume)</b>	<b>Volume</b>
Label Probe Diluent QF (pre-warmed to 40°C)	399.6 $\mu$ L
Label Probe 1-AP	0.4 $\mu$ L

- c. Remove each slide and flick to remove the Wash Buffer. Without completely drying out the sections, tap the slide on its edge and then wipe the backside on a laboratory wipe. Place the slides face up on a flat, elevated platform and immediately add 400  $\mu$ L of Working Label Probe 1-AP Solution to each tissue section.
  - d. Transfer the slides to the hybridization system and incubate at 40°C for 15 min.
22. Wash slides (12 minutes)
    - a. Insert an empty slide rack into a clear staining dish containing 200 mL of Wash Buffer.

- b. After incubation, decant the Working Label Probe 1-AP Solution from the slides and insert them into the slide rack.
  - c. Wash the slides 3 times, each time with 200 mL of fresh Wash Buffer at RT for 3 min with constant and vigorous agitation.
23. Apply Fast Red Substrate (45 minutes)
- a. Remove each slide and flick it to remove the Wash Buffer. Without completely drying out the sections, tap the slide on its edge and then wipe the backside on a laboratory wipe. Place slides face up on a flat, elevated platform.
  - b. Immediately add 400  $\mu$ L of the AP Enhancer Solution to each tissue section and incubate at RT for 5 min while preparing the Fast Red Substrate.
  - c. Prepare the Fast Red Substrate: Add 5 ml of Naphthol Buffer and one Fast Red Tablet to a 15 ml conical tube. Vortex at high speed to completely dissolve the tablet. Protect from light until use by wrapping the tube in aluminum foil.
  - d. Decant the AP Enhancer Solution and flick the slide twice to completely remove any excess AP Enhancer Solution. Tap the slide on its edge then wipe the backside on a laboratory wipe. Immediately add 400  $\mu$ L of Fast Red Substrate onto each tissue section.
  - e. Transfer the slides to the hybridization system and incubate at 40°C for 30 min.
  - f. Insert an empty slide rack into a clear staining dish containing 200 mL of 1X PBS.
  - g. After incubation, decant the Fast Red Substrate solution from the slides and insert them into the slide rack.
  - h. Rinse off the excess Fast Red Substrate from the slides by moving the slide rack up and down for 1 min.
24. Counterstain (25 minutes)
- a. Transfer the slide rack to the clear staining dish containing the 200 mL of Gill's hematoxylin and stain for 5-10 sec at RT.
  - b. Wash the slides 3 times, each time with 200 mL of fresh ddH<sub>2</sub>O for 1 min by moving the slide rack up and down.
  - c. Pour off the ddH<sub>2</sub>O, refill with 200 mL of 0.01% ammonium hydroxide and incubate the slides for 10 seconds. Unused 0.01% ammonium hydroxide can be stored at RT for up to 1 month.
  - d. Wash the slides once more in 200 mL of fresh ddH<sub>2</sub>O by moving the rack up and down for 1 minute.
  - e. Optional, if you plan to view slides using a fluorescent microscope: Move the slide rack into a clear staining dish containing 200 mL DAPI (3  $\mu$ g/mL). Stain the slides for 1 min, then rinse them in 200 mL of fresh ddH<sub>2</sub>O by moving the slide rack up and down for 1 min.

- f. Remove the slides from the slide rack and flick to remove the excess ddH<sub>2</sub>O. Tap the slide on its edge then wipe the backside on a laboratory wipe. Place them face up onto a paper towel to air dry in the dark.
  - g. Ensure that slide sections are completely dry before mounting (~20 min).
25. Mount and image (40 minutes)

**If using DAKO Ultramount mounting medium, with no coverslipping (20X viewing or imaging):**

- a. Place the slide flat on a counter top with specimen facing up.
- b. Dab the first 2-3 drops of Ultramount onto a paper towel to remove bubbles.
- c. Apply a sufficient amount of Ultramount to completely cover the specimen with a thin layer (3-4 drops) of mounting medium.
- d. Place slides horizontally in a 70°C oven/incubator to dry the mounting medium. Allow 10-30 minutes for the mounting medium to harden completely. The drying time depends on the amount of mounting medium applied.
- e. Image or store slides at RT.

**If using DAKO Ultramount mounting medium, and post mounting with coverslip (crisper 20X or 40X viewing or imaging):**

- a. Work in a fume hood and follow the procedure above for no coverslipping. Make sure that the Ultramount is completely hardened, and allow the slides to come to RT.
- b. Apply HistoMount directly on top of the dried Ultramount.
- c. Place coverslip on top and allow to air dry at RT for 15 min.
- d. Image or store slides at RT.

**If using ADVANTAGE mounting medium:**

- a. Place a 24 mm x 55 mm cover glass horizontally onto a clean, flat surface.
- b. Dab the first 2-3 drops of mounting media onto a paper towel to remove bubbles.
- c. Add 2 drops of the ADVANTAGE medium directly onto the middle of the cover glass.
- d. Use a pipette tip to draw out any air bubbles in the droplets.
- e. Invert the specimen slide and slowly place it onto the mounting medium at an angle. Make sure that the tissue comes into contact with the mounting medium first before completely letting go of the glass slide to overlap with the cover glass.
- f. After mounting, flip the slide over and place it on its edge on a laboratory wipe to soak up and remove excess mounting medium. Allow the slide to dry at RT in the dark for 15 min. Do not bake the slides to speed up the drying process.

- g. To prevent bubble formation, seal all 4 edges of the cover glass with a flat black-colored nail polish (iridescent or colored nail polish can autofluoresce and interfere with fluorescent imaging).
- h. Image the results using a brightfield and/or fluorescence microscope. Store slides at RT.

# Troubleshooting

## Weak or No Signals

Probable Cause	Recommended Action
Incorrect pretreatment conditions	<p>Repeat pretreatment assay optimization procedure to determine optimal heat treatment time and protease digestion time that will strike a balance between morphology and signal.</p> <ul style="list-style-type: none"> <li>• Under-pretreatment yields good morphology but poor signal due to insufficient unmasking of target.</li> <li>• Over-pretreatment yields poor morphology and loss of signal due to over digestion.</li> </ul>
Sample preparation	<p>Immediately place freshly dissected tissues in <math>\geq 20</math> volumes of fresh 10% neutral buffered formalin (NBF) or 4% paraformaldehyde (PFA) at RT for 16-24 hours.</p>
Tissue over-fixed after protease digestion	<p>Make sure the tissue sections are not fixed for more than 5 min in 10% NBF after protease digestion.</p>
RNA in tissue is degraded	<p>Verify tissue fixation:</p> <ul style="list-style-type: none"> <li>• Immediately place freshly dissected tissues in <math>\geq 20</math> volumes of fresh 10% neutral buffered formalin (NBF) or 4% paraformaldehyde (PFA) for 16-24 hours at RT.</li> <li>• If fixation cannot be performed immediately, be sure that the tissue is placed on dry ice or in liquid nitrogen to prevent RNA degradation.</li> <li>• Use positive control probe set(s) such as one for a housekeeping gene or a housekeeping gene panel (ACTB, GAPD and UBC) to assess RNA integrity.</li> </ul>
Reagents applied in wrong sequence	<p>Apply target probe sets, Pre-amplifier Mix QT, Amplifier Mix QT, Label Probe-AP and substrates in the correct order.</p>
Gene of interest not expressed	<ul style="list-style-type: none"> <li>• Verify expression using other tissue lysate methods such as QuantiGene™ Singleplex or QuantiGene Plex Assays.</li> <li>• Run the same probe set on known samples that have been validated to express the target of interest.</li> </ul>

Incorrect storage condition	Store the components at the storage condition as written on the component label or kit boxes.
Hybridization temperature not optimal	Calibrate the hybridization system at 40°C using a ViewRNA™ Temperature Validation Kit (QV0523).
Mounting solution contained alcohol	Use the recommended mounting media to mount your tissue (see Step 25: Mount and Image). Avoid any mounting solution containing alcohol.
Tissue dries up during hybridization steps	<p>Recommendations for hybridization systems:</p> <ul style="list-style-type: none"> <li>• Ensure the hybridization system is appropriately humidified and that door/lid is closed during hybridization steps.</li> <li>• Make sure the hybridization system is placed on a level bench.</li> <li>• Calibrate the hybridization system to 40 °C using the ViewRNA Temperature Validation Kit (QV0523).</li> </ul> <p>Prevent sections from drying out:</p> <ul style="list-style-type: none"> <li>• Prepare enough reagents and use the recommended volumes for each step of the assay.</li> <li>• Ensure that you have a solid seal when drawing your hydrophobic barriers.</li> <li>• Add all working reagents onto the slides before moving them to the 40°C hybridization system.</li> </ul>
Tissue dries up during processing	<p>Keep tissue sections moist starting from the heat pretreatment step:</p> <ul style="list-style-type: none"> <li>• Add respective reagents immediately after decanting solution from the slides.</li> <li>• Keep tissue exposure to air as short as possible before adding hybridization reagents.</li> <li>• Add all working reagents onto the slides before moving them to the 40°C hybridization system.</li> </ul>
Fast Red and Fast Blue Substrate solutions not freshly prepared	Prepare Fast Red and Fast Blue Substrate solutions immediately before use.
Small targets, splice variants or RNA fusions	<p>Doing one or both of the following may increase sensitivity, but it should be noted that there is always a general trade-off between specificity and sensitivity:</p> <ul style="list-style-type: none"> <li>• Increase probe set concentration by diluting target probe set 1:30 instead of 1:40 and hybridize for 2 hr.</li> <li>• Decrease hybridization temperature from 40 to 38°C.</li> <li>• Increase Fast Red incubation time to 45 min.</li> </ul>



Probe set hybridization temperature, time, and/or concentration not optimal	Decrease hybridization temperature from 40°C to 38°C and increase the probe set concentration by diluting the target probe set 1:30 instead of 1:40. Hybridize for 2 hr.
Label Probe-AP concentration too low	<ul style="list-style-type: none"> <li>• Verify that the correct concentrations were used.</li> <li>• Increase the recommended concentration for Label Probe-AP. If this is necessary, it may result in higher background.</li> </ul>
Dark hematoxylin stain reduces visibility of the blue dots	<ul style="list-style-type: none"> <li>• Tissues with lower cell density require longer hematoxylin incubation than tissues with higher cell density. It may be helpful to titrate incubation times.</li> <li>• Increase the lamp brightness during viewing.</li> <li>• View under a 40X objective.</li> <li>• Image using fluorescent mode.</li> </ul>

## High Background

Probable Cause	Recommended Action
Tissue dries up during processing	<p>Prevent tissue sections from drying out after the pretreatment step:</p> <ul style="list-style-type: none"> <li>• Ensure that you have a solid seal when drawing your hydrophobic barrier.</li> <li>• Prepare enough reagents and use the recommended volume for each step of the assay.</li> <li>• Add respective reagents immediately after decanting solution from the slides.</li> <li>• Keep tissue exposure to air as short as possible before adding hybridization reagents.</li> <li>• Make sure that the hybridization system is appropriately humidified.</li> <li>• Make sure the hybridization system is set at 40 °C and that the lid/door is closed during hybridization steps.</li> <li>• Process as few or as many slides at a time as you are comfortable doing.</li> </ul>
Incomplete removal of paraffin	<ul style="list-style-type: none"> <li>• Use fresh xylene or Histo-Clear solution.</li> <li>• Immediately submerge the warm slides into the Histo-Clear solution after baking.</li> </ul>
Insufficient washing	<ul style="list-style-type: none"> <li>• Move the slide rack up and down with constant and vigorous agitation. Click for a helpful video.</li> <li>• Increase wash incubation time by 1 min per wash.</li> </ul>
Hybridization temperature not optimal	Calibrate the hybridization system at 40°C using the ViewRNA™ Temperature Validation Kit (QV0523).

Concentration of hybridization reagents too high	Double check the dilution calculation for all working solutions.
Suboptimal pretreatment conditions	Perform the pretreatment optimization procedure to determine the optimal heat treatment and protease digestion time.
Label Probe-AP concentration too high	<ul style="list-style-type: none"> <li>• Verify that the correct concentrations were used.</li> <li>• Decrease the recommended concentration for Label Probe-AP.</li> </ul>

## Diffused Signals

Probable Cause	Recommended Action
Tissue dries up during processing	<p>Prevent tissue sections from drying out after the pretreatment step:</p> <ul style="list-style-type: none"> <li>• Ensure that you have a solid seal when drawing your hydrophobic barrier.</li> <li>• Prepare enough reagents and use the recommended volume for each step of the assay.</li> <li>• Add respective reagents immediately after decanting solution from the slides.</li> <li>• Limit tissue exposure to air before adding hybridization reagents.</li> <li>• Make sure that the hybridization system is appropriately humidified.</li> <li>• Make sure the hybridization system is set at 40 °C and that the lid/door is closed during hybridization steps.</li> <li>• Process as few or as many slides at a time as you are comfortable doing.</li> </ul>
Incomplete removal of AP Enhancer	Ensure that excess AP Enhancer is removed by decanting the AP Enhancer and flicking the slides twice prior to adding Fast Red Substrate.
Insufficient washing	Make sure tissues are washed twice in 1X PBS after protease digestion and twice again after subsequent fixing in 10% NBF.
Fast Red Substrate and Fast Blue Substrate solutions not freshly prepared	Prepare Fast Red and Fast Blue Substrate solutions immediately before use.
Slides are not dried before mounting	Ensure that the sections are completely dry (~20 min) before mounting.
Mounting solution contained alcohol	Use the recommended mounting media to mount your tissue (see Step 25: Mount and Image). Avoid any mounting medium containing alcohol or any cover slipping method requiring alcohol dehydration.

## Endogenous Alkaline Phosphatase Activity

Probable Cause	Recommended Action
Endogenous alkaline phosphatase activity	Verify alkaline phosphatase activity by incubating protease-treated sample with Fast Red Substrate or Fast Blue Substrate. If endogenous AP activity is present, diffused signals (which can be weak or strong) will appear. Inactivate endogenous AP with 0.2 M HCl at RT for 10 min before the protease step. Wash samples twice with 1X PBS before proceeding to protease digestion.

## Tissue Detachment from Slide

Probable Cause	Recommended Action
Improper tissue preparation	Make sure that the tissue is prepared as recommended in the Tissue Preparation Guidelines, including fixation time and reagent, thickness of sections, brand of positively charged glass slide, and baking of the sections at 60 °C for 1 hr before storing at –20 °C.
Insufficient baking of slides	Verify that the 60 min at 60 °C baking step was performed prior to storage of slides at –20 °C and again just before the deparaffinization step to ensure adhesion of tissue to slide.
Incorrect pretreatment conditions	Perform full pretreatment optimization procedure to determine optimal heat treatment and protease digestion time.
Temperature of heat pretreatment condition too high	Make sure the temperature is within the tolerance range of 90-95°C. For fatty soft tissue such as breast, adjust to 90°C.
Proteinase treatment is too long or at too high of a concentration.	Reduce proteinase concentration and/or incubation time.

## Poor Cell Morphology

Probable Cause	Recommended Action
Incorrect pretreatment conditions	Perform full pretreatment optimization procedure to determine optimal heat treatment and protease digestion time. See Appendix A.
Tissue sample not fixed properly	Make sure that freshly dissected tissues are fixed in 10% NBF or 4% PFA for 16-24 hr.
Section thickness is variable or not optimal	Make sure microtome is calibrated and tissue is sectioned at $5 \pm 1 \mu\text{m}$ .

## High Non-Specific Binding on Glass Slide

Probable Cause	Recommended Action
Incompatible glass slide	<ul style="list-style-type: none"> <li>• Use the recommended glass slides:</li> <li>• Leica Non-Clipped X-tra® Slide, 1 mm White P/N 3800200 or 3800210</li> <li>• Fisherbrand™ Superfrost™ Plus Slides, white label (Fisher Scientific, P/N12-550-15); avoid other colored labels as they tend to give high background.</li> <li>• Prevalidate each new batch of slides by running the entire assay, including probe set on empty slides with hydrophobic barriers (without fixed tissues) to determine if the slides are suitable for the assay.</li> </ul>
Insufficient washing	<ul style="list-style-type: none"> <li>• Move the slide rack up and down with constant and vigorous agitation.</li> <li>• Increase wash incubation time by 1 min per wash.</li> </ul>
Concentration of hybridization reagents was too high	Confirm that the dilution calculations are correct for all working solutions.

## Pink Non-Specific Background Where Paraffin Was

Probable Cause	Recommended Action
Incomplete removal of paraffin	<ul style="list-style-type: none"> <li>• Be sure to use fresh Histo-Clear or xylene for the indicated amount of time during the dewaxing step.</li> <li>• Use 3 changes of Histo-Clear instead of 2 changes.</li> </ul>
Polymerization of poor quality paraffin	<ul style="list-style-type: none"> <li>• Melt the paraffin at 80°C for 3 min and remove paraffin using 3 changes of fresh Histo-Clear.</li> <li>• Do not bake the slides at a temperature higher than 60°C.</li> </ul>

## Hydrophobic Barrier Falls Off

Probable Cause	Recommended Action
Incompatible glass slide	<ul style="list-style-type: none"> <li>• Use the recommended glass slides: <ul style="list-style-type: none"> <li>○ Leica Non-Clipped X-tra® Slide, 1 mm White P/N 3800200 or 3800210</li> <li>○ Fisherbrand™ Superfrost™ Plus Slides, white label (Fisher Scientific, P/N12-550-15); avoid other colored labels as they tend to give high background.</li> </ul> </li> <li>• Prevalidate each new batch of slides by drawing a hydrophobic barrier onto an empty slide (without fixed tissue), allow it to dry for 20-30 min, boil in pretreatment solution for 40 min to determine if the hydrophobic barrier is intact and the slides are suitable for the assay.</li> </ul>
Incorrect hydrophobic pen	Use the recommended Hydrophobic Barrier Pen (QVC0500 or Vector Laboratories H4000).
Hydrophobic barrier was not completely dried	Be sure that the hydrophobic barrier is completely dry before proceeding to the next step. This can be 20-30 min or longer depending on how heavily the barrier is created.

## Fast Red Signal for TYPE 1 Target is Weak or Different in 2-Plex vs. 1-Plex

Probable Cause	Recommended Action
Cross-inhibition of LP1-AP by Fast Blue precipitate	Assign lower expressing target to TYPE 6 (Fast Blue) and higher expressing target to TYPE 1 (Fast Red).
Co-localization of TYPE 1 and TYPE 6 targets	<ul style="list-style-type: none"> <li>Perform a 1-plex assay for each target.</li> <li>Assign lower expressing target to TYPE 6 (Fast Blue) and higher expressing target to TYPE 1 (Fast Red).</li> <li>If co-localization study is desired, try reducing development time for Fast Blue from 30 min to 10-15 min.</li> </ul>

## TYPE 1 Target Signals Observed in the Channel for TYPE 6 Target

Probable Cause	Recommended Action
Spectral bleed through of Fast Red signal	Check to make sure that the filter set for Fast Blue is as recommended.
Incorrect filter set for Fast Blue signal	Use the correct filter set. See the Microscopy and Imaging Equipment Guidelines under “General Information” for recommended filter set specifications for Fast Blue.

## Co-localized Fast Blue and Fast Red Signals When Using Only TYPE 6 Probe Set in a 2-plex Assay

Probable Cause	Recommended Action
Residual LP6-AP activity	<ul style="list-style-type: none"> <li>Do not omit Step 20: Quench Label Probe 6-AP.</li> <li>Be sure to quench LP6-AP activity with AP Stop QT for the entire 30 min.</li> </ul>

# Appendix A: Sample Pretreatment Optimization Procedures

## About Pretreatment Optimization

Critical to any *in situ* assay is the balance between the adhesion of the tissue to the glass surface, crosslinking of the target molecules to the cellular structures by chemical fixatives and the subsequent unmasking of the RNA targets by heat treatment and protease digestion for the probes to hybridize. For the ViewRNA ISH Tissue Assay, this balance between signal strength and tissue morphology is largely sample dependent (tissue types as well as the modes of fixation and sample preparation) and can be achieved by optimizing the pretreatment conditions to empirically determine the optimal time for heat treatment and protease digestion.

When optimizing the pretreatment conditions for your tissue type, choose a target that is known to be expressed in the tissue of interest with medium to medium-high levels of expression. This will avoid possible signal saturation that may be associated with extremely high expressing targets and allow for detectable changes in the signals to be assessed as a function of the different pretreatment conditions. In general, a housekeeping gene with medium-high expression, such as GAPD or ACTB, can be used for this purpose. Once the optimal pretreatment conditions are determined, they can generally be used for most targets within the particular tissue. If the transcript is expressed at an extremely low level, the optimal pretreatment condition may need to be one that favors signal over morphology.

## Sample Pretreatment Optimization Setup

Ten FFPE tissue sections from the same block are treated with different set of pretreatment conditions prior to target probe hybridization step. Slide 7 serves as a "no probe control", while the remaining 9 slides are processed with the control target probe set.

Table B.1 in Appendix B provides sample pretreatment conditions for some common tissues. If samples are limited, see Table B.2 in Appendix B.

**Table A.1:** Pretreatment Optimization Setup

Protease Incubation Time (min)	Heat Pretreatment Time (min)			
	0	5	10	20
0	Slide 1 Morphology reference			
10		Slide 2	Slide 5	Slide 9
20		Slide 3	Slide 6 Slide 7 No Probe Control	Slide 10
40		Slide 4	Slide 8	

Before starting the pretreatment optimization protocol, please read the sections on Important Procedural Notes and Guidelines and Essential Keys for a Successful Assay, in the “About the ViewRNA ISH Tissue Assay Procedure” section under “Assay Procedure”.

The pretreatment optimization procedure for the ViewRNA ISH Tissue Assay is divided into two parts that can be performed in a single day or over two days:

- Part 1: Sample Preparation and Target Probe Set Hybridization (optional stopping point).
- Part 2: Signal Amplification and Detection.

We do not recommend stopping the procedure at any point in the assay unless specifically indicated.

## Sample Preparation and Target Probe Hybridization

1. Bake slides
  - a. See Step 1 to Step 4, under “Assay Procedure”.
2. Heat Pretreatment (10-25 minutes)
  - a. Tightly cover the beaker containing the 500 mL of 1X Pretreatment Solution with aluminum foil, place it on a hot plate and heat the solution to a temperature of 90-95°C. Use a waterproof probe thermometer to measure and maintain the temperature of the solution at 90-95°C during the pretreatment period.
  - b. Set slide 1 aside on the lab bench.
  - c. Load slides 9 and 10 into the vertical slide rack.
  - d. Using a pair of forceps, submerge the slide rack into the heated 1X Pretreatment Solution. Cover the glass beaker with aluminum foil and incubate at 90-95°C for 10 min.
  - e. At the end of the 10 min, add slides 5, 6, 7 and 8 to the rack in the 90-95°C 1X Pretreatment Solution. Cover the glass beaker with aluminum foil and incubate for 5 min.



- f. At the end of the 5 min, add slides 2, 3, 4 into the rack in the 90-95°C 1X Pretreatment Solution. Cover the glass beaker with aluminum foil and incubate for 5 min.
  - g. After the pretreatment, remove the slide rack with forceps, submerge it into a clear staining dish containing 200 mL of ddH<sub>2</sub>O and wash for 1 min with frequent agitation.
  - h. Repeat the wash one more time with another 200 mL of fresh ddH<sub>2</sub>O.
  - i. Transfer the slide rack to a clear staining dish containing 1X PBS.
  - j. **IMPORTANT:** From this point forward, do not let the tissue sections dry out. Tissue sections that have been heat treated can be stored covered in 1X PBS at RT for up to one week. Continue with Step 3 below when ready.
3. Protease Digestion and Fixation (30-50 minutes)
- a. Prepare the Working Protease Solution using the table below as a guide. Dilute the Protease QF 1:100 in pre-warmed 1X PBS and briefly vortex to mix. Scale reagents according to the number of assays to be run. Include one slide volume overage.

Reagent	Volume
Protease QF	4 µL
1X PBS (pre-warmed to 40°C)	396 µL

- b. Leave slide 1 on the lab bench as it is excluded from this step.
- c. Begin by removing slides 4 and 8 and flicking each to remove excess 1X PBS. Tap the slide on its edge then wipe the backside on a laboratory wipe. Leave remaining slides in 1X PBS.
- d. Place slides 4 and 8 face up on a flat, elevated platform (e.g., an Eppendorf tube rack for ease of handling) and immediately add 400 µL of the working protease solution onto the tissue section. It may be necessary to spread the solution with a pipette tip.
- e. Transfer the slides to the hybridization system and incubate at 40°C for 20 min.
- f. After 19 min, remove slides 3, 6, 7 and 10 from the clear staining dish and flick off excess 1X PBS. Without completely drying out the sections, tap the slides on their edges and then wipe the backsides on a laboratory wipe.
- g. Place slides 3, 6, 7, and 10 face up on a flat, elevated platform and immediately add 400 µL of the working protease solution onto the tissue section.
- h. Transfer the slides to the hybridization system and incubate at 40°C for 10 min.
- i. After 9 min, remove slides 2, 5 and 9 from the clear staining dish and flick off excess 1X PBS. Without completely drying out the sections,

tap the slides on their edges and then wipe the backsides on a laboratory wipe.

- j. Place slides 2, 5, and 9 face up on a flat, elevated platform and immediately add 400  $\mu$ L of the working protease solution onto the tissue section.
- k. Transfer the slides to the hybridization system and incubate at 40°C for 10 min.
- l. Pour 200 mL of 1X PBS into a clear staining dish and insert an empty slide rack into it.
- m. At the end of 10 min (40 min total of incubation time), decant the working protease solution from the slides, insert the slides into the rack and wash gently but thoroughly by moving the rack up and down for 1 min.
- n. Repeat the wash one more time with another 200 mL of fresh 1X PBS before adding slide 1 to the rack.
- o. Transfer the slide rack, containing all 10 slides, to a clear staining dish containing 200 mL of 10% NBF and fix at RT for 5 min under a fume hood.
- p. Wash the slides twice, each time with 200 mL of fresh 1X PBS for 1 min with frequent agitation.
- q. Proceed to Step 7: Target Probe Set Hybridization under “Assay Procedure” Part 1 to continue the assay.

# Appendix B: Sample Pretreatment Lookup Table

Table B.1 shows a list tissues that were prepared according to the guidelines outlined in this manual (Tissue Preparation Guidelines under “Before You Begin”) and optimized using the recommended pretreatment assay optimization procedure. This table provides a reference or a starting point to minimize the number of test conditions if you do not have sufficient slides to perform the full recommended pretreatment optimization procedure.

Please note that the conditions listed here are specific to tissues prepared using 10% NBF and may not be applicable to tissue prepared using 4% PFA. If you chose to use any of the pretreatment conditions listed in the lookup table, it is important to include a "negative control" slide to assess whether the assay background is clean and cellular morphology is well-defined.

**Table B.1: Sample Pretreatment Optimization Lookup Table**

Species	Tissue	Optimal Conditions (minutes)		Range of Tolerance (Heat pretreatment, Protease)
		Heat pretreatment at 90-95°C	Protease at 40°C	
Human	Brain	20	10	(10, 10) (10, 20)
	Breast	20	15	(25, 15) (30, 20) (25, 20)
	Colon	5	20	(5, 10)
	Kidney	20	10	
	Liver	20	20	(10, 20)
	Lung	10	20	
	Lymph Node	10	20	
	Nasal polyp	5	5	
	Osteoarthritic tissue	20	20	
	Pancreas	10	10	(10, 20) (5, 10)
	Prostate	10	20	(5, 10) (20, 10) (10, 10)
	Salivary gland	10	10	(5, 10)

	Skin	5	10	
	Tonsil	10	20	
	Thyroid	10	20	
<b>Rat</b>	Kidney	10	20	(10, 10) (20, 20)
	Liver	10	20	
	Spleen	20	10	
	Thyroid	10	20	
<b>Mouse</b>	Bone	20	20	
	Brain	10	10	
	Heart	10	40	(20, 20)
	Kidney	20	20	(10, 20)
	Liver	20	20	(5, 40) (10, 20)
	Lung	10	20	
	Retina	10	10	
<b>Salmon</b>	Heart	10	10	
	Muscle	10	20	
<b>Monkey</b>	Mucosal rectum	10	20	

If your tissue type is not listed in Table B.1, and you have only limited slides available for the pretreatment optimization, Table B.2 provides the recommended heat treatment and protease incubation times that will likely give the best chance of achieving an acceptable pretreatment conditions for your ViewRNA™ ISH Tissue Assay.

**Table B.2: Heat Treatment and Protease Incubation Times for Limited Optimization**



Number of Available Slides	Heat Pretreatment Time (min)	Protease Time (min)
3	5	10
	10	10
	10	20
5	5	10
	5	20
	10	10
	10	20
	20	10
7	5	10
	5	20
	10	10
	10	2
	20	10
	20	20
	0	0




# Appendix C: Evaluating Results

## Assessing Pretreatment Conditions

This section provides sample images obtained from the ViewRNA ISH Tissue Assay, performed on rat kidney tissue, to illustrate the effects of optimal and suboptimal pretreatment conditions on Arbp signal strength versus morphology and to demonstrate how data gathered from the in situ assay can be analyzed to determine target expression.

**Table C.1: Assessing Pretreatment Conditions: Synpo and SPP1 Expression in Rat Kidney FFPE Tissue**

Heat Pretreatment Time (min)	Protease Digestion Time (min)	Brightfield Image	Results Interpretation
0	0		<b>Untreated Morphology Reference Slide; +Probes (Synpo and SPP1)</b> <ul style="list-style-type: none"> <li>• Good morphology</li> <li>• Intact cellular structure</li> <li>• Good hematoxylin counterstaining of nuclei</li> <li>• Little or no signal (dots) observed</li> </ul>
5	10		<b>Insufficient Pretreatment or Over Fixation of Tissue; +Probes (Synpo and SPP1)</b> <ul style="list-style-type: none"> <li>• Good morphology</li> <li>• Intact cellular structure</li> <li>• Strong hematoxylin counterstaining of nuclei</li> <li>• Weak, diffused and non-ubiquitous signal</li> <li>• Few number of dots</li> </ul>

10	20		<b>Optimal Pretreatment and Sample Preparations; +Probes (Synpo and SPP1)</b> <ul style="list-style-type: none"> <li>• Good morphology</li> <li>• Cellular structures and boundaries are retained and still identifiable</li> <li>• Good hematoxylin counterstaining of nuclei</li> <li>• Strong, punctated and ubiquitous signals in (+) probe sample and clean background in (-) probe sample</li> </ul>
10	20		<b>Optimal Pretreatment and Sample Preparations; No Probes</b> <ul style="list-style-type: none"> <li>• Clean background</li> <li>• Acceptable morphology and cellular architecture</li> <li>• Good hematoxylin counterstaining of nuclei</li> </ul>
10	40		<b>Over Pretreatment or Under Fixation; + Probes (Synpo and SPP1)</b> <ul style="list-style-type: none"> <li>• Poor morphology</li> <li>• Loss of cellular structure and boundaries due to excessive heat treatment and protease digestion</li> <li>• Poor hematoxylin counterstaining of nuclei</li> <li>• Weak signal and fewer number of dots</li> </ul>

## Analyzing Target Expression

Each observable punctated dot represents a single RNA molecule within the cell that the ViewRNA ISH Tissue Assay is able to detect, assuming the RNA target is intact and properly unmasked for the probe to access. These dots are typically uniform in size. However, smaller than average size dots can also be present, and this usually

indicates that the transcript is not properly unmasked, or that the RNA target is not intact, resulting in the binding of only one or a few pairs of oligonucleotides from the probe set.

Conversely, a larger than average size dot can occur when multiple targets are found clustered in the same physical area. Naturally, with everything being equal, an RNA target with a low expression will yield fewer numbers of dots than one with a high expression.

In quantifying the results to assess the RNA target expression, it is important to consider the pattern and number of dots observed in the “Negative Control”, such as bacterial *dapB* or sense strand of the target, in order to confidently differentiate between low expressing targets and non-specific background dots. The ViewRNA ISH Tissue Assay typically has an average background of < 1 dot/10 cells. Consequently, as long as your target is consistently showing an expression level above the “Negative Control” threshold, even if the RNA target expression is extremely low (e.g., 1 dot/every 2 cells), you can trust that the detection is real.



# Appendix D: Using Frozen Tissues with ViewRNA ISH Tissue Assay

## About this Appendix

This appendix provides procedural modification for running the ViewRNA ISH Tissue Assay on fresh frozen or OCT-embedded frozen tissue sections. This modified assay protocol has been tested on the following OCT-embedded frozen tissue samples:

- Bovine – Ovary
- Human – Colon, skin, testis
- Mouse – Brain, duodenum, eye, liver, lung, pancreas, skin, spinal cord
- Rat – Brain, spinal cord

## Important Procedural Notes

- This protocol requires overnight fixation in chilled 4% NBF prior to starting Day 1 of the assay. See Step 2 under the “Modifications to Part 1” below.
- Samples should be freshly sectioned at  $12 \pm 1 \mu\text{m}$  and mounted onto one of the following positively charged glass slides.
  - Leica Non-Clipped X-tra® Slides, 1 mm white (P/N 3800200 in US, Canada, and Asia Pacific regions or P/N 3800210 in Europe).
  - Fisherbrand™ Superfrost™ Plus Slides, white label (Fisher Scientific Cat# 12-550-15); avoid other colored labels as they tend to give high background.
- Prepared frozen tissue slides should be used immediately in the assay or can be stored at  $-80^{\circ}\text{C}$  for up to 6 months.
- Perform optimization for Protease only (see Table A.1 in Appendix A for recommended times). No heat treatment step required.

## Modifications to Part 1: Sample Preparation and Target Probe Hybridization for Frozen Tissues

The following procedural steps replace Step 1 to Step 5 in Part 1: Sample Preparation and Target Probe Hybridization under “Assay Procedure”:

1. Prepare and chill 10% NBF: Add 178 mL 1X PBS and 22 mL 37% formaldehyde to a 200 mL capacity container. Mix well and chill on ice for 1 hr.
2. Fix tissue overnight: Pour chilled 10% NBF into a clear staining dish and insert an empty slide rack into the solution. Insert frozen tissue slides into the slide rack and incubate at 4°C for 16-18 hr.
3. Prepare buffers, reagents, and equipment:
  - a. Verify that the hybridization system is set to  $40 \pm 1^\circ\text{C}$  and appropriately humidified.
  - b. Prepare 2 L 1X PBS: add 200 mL 10X PBS and 1.8 L ddH<sub>2</sub>O to a 2 L capacity container.
  - c. Prepare 200 mL of 50% ethanol: add 100 mL 100% ethanol and 100 mL ddH<sub>2</sub>O to a 200 mL capacity container.
  - d. Prepare 200 mL of 70% ethanol: add 60 mL 100% ethanol and 140 mL ddH<sub>2</sub>O to a 200 mL capacity container.
  - e. Prepare 4 L Wash Buffer: add the components below in the order listed to a 4 L capacity container and mix well:
    - 3 L ddH<sub>2</sub>O
    - 36 mL Wash Comp 1
    - 10 mL Wash Comp 2
 Adjust the total volume to 4 L with ddH<sub>2</sub>O.
  - f. Prepare 200 mL Storage Buffer (for optional stopping point): add 60 mL Wash Comp 2 and 140 mL ddH<sub>2</sub>O to a 200 mL capacity container.
  - g. Ensure availability of:
    - 1000 mL ddH<sub>2</sub>O
    - 200 mL Gill's Hematoxylin I
    - 200 mL of 3 µg/mL DAPI in 1X PBS (optional for fluorescent detection), store in the dark at 4°C until use.
  - h. Thaw probe set(s). Mix, briefly centrifuge to collect content and place on ice until use.
  - i. Pre-warm 10 mL of 1X PBS and Probe Set Diluent QF to 40°C.
  - j. Optional for 1-day assay:
    - Pre-warm Pre-amplifier Mix QT, Amplifier Mix QT, and Label Probe Diluent QF to 40°C.
    - Briefly spin down the Label Probe 1-AP, Label Probe 6-AP, and Blue Reagents. Place on ice.
    - Bring Fast Red Tablets, Naphthol Buffer, Blue Buffer, and AP Enhancer Solution to RT.
    - Working in a fume hood, prepare 1 L of 0.01% ammonium hydroxide. Add 0.33 mL 30% ammonium hydroxide and 999.67 mL ddH<sub>2</sub>O in a 1 L capacity container.
  - k. Optional: If using a microplate shaker for the washes, set the speed to 550 rpm.

4. Wash slides: Remove slide rack from the 10% NBF and wash the slides twice, each time with 200 mL of 1X PBS for 1 min with frequent agitation.
5. Tissue dehydration:
  - a. Dehydrate the tissue by sequentially soaking the rack of slides in 50%, 70% and then 100% ethanol in a clear staining dish, each time at RT for 10 min without agitation.
  - b. Remove the slide rack from the 100% ethanol and drain the excess on a paper towel.
  - c. Transfer the entire rack of slides to a 60°C dry incubator/oven and bake the slides for 60 min.

**IMPORTANT:** Following the baking step, continue with Part 1 of the assay procedure within 1 hour, beginning with Step 4: Draw Hydrophobic Barrier and skipping to Step 6 Protease Digestion and Fixation. The heat treatment step is NOT REQUIRED for frozen tissues.