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| 製品コード | 製品名 | 梱包単位 | 保存温度 |
|--------------|---|--------|------|
| IVT-X008061 | <i>InVitroCYP</i> H-class 10-Donor Pooled Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-X008062 | <i>InVitroCYP</i> H-class 10-Donor Pooled Human Liver Microsomes | 20 mg | 液体窒素 |
| IVT-X008063 | <i>InVitroCYP</i> H-class 10-Donor Pooled Human Liver Microsomes | 100 mg | 液体窒素 |
| IVT-X008064 | <i>InVitroCYP</i> H-class 25-Donor Pooled Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-X008065 | <i>InVitroCYP</i> H-class 25-Donor Pooled Human Liver Microsomes | 20 mg | 液体窒素 |
| IVT-X008066 | <i>InVitroCYP</i> H-class 25-Donor Pooled Human Liver Microsomes | 100 mg | 液体窒素 |
| IVT-FX008061 | <i>InVitroCYP</i> H-class 10-Donor Pooled Female Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-FX008062 | <i>InVitroCYP</i> H-class 10-Donor Pooled Female Human Liver Microsomes | 20 mg | 液体窒素 |
| IVT-MX008061 | <i>InVitroCYP</i> H-class 10-Donor Pooled Male Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-MX008062 | <i>InVitroCYP</i> H-class 10-Donor Pooled Male Human Liver Microsomes | 20 mg | 液体窒素 |
| IVT-FX008064 | <i>InVitroCYP</i> H-class 25-Donor Pooled Female Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-FX008065 | <i>InVitroCYP</i> H-class 25-Donor Pooled Female Human Liver Microsomes | 20 mg | 液体窒素 |
| IVT-MX008064 | <i>InVitroCYP</i> H-class 25-Donor Pooled Male Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-MX008065 | <i>InVitroCYP</i> H-class 25-Donor Pooled Male Human Liver Microsomes | 20 mg | 液体窒素 |
| IVT-F008084 | <i>InVitroCYP</i> H-class Individual Female Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-F008085 | <i>InVitroCYP</i> H-class Individual Female Human Liver Microsomes | 20 mg | 液体窒素 |
| IVT-M008084 | <i>InVitroCYP</i> H-class Individual Male Human Liver Microsomes | 10 mg | 液体窒素 |
| IVT-M008085 | <i>InVitroCYP</i> H-class Individual Male Human Liver Microsomes | 20 mg | 液体窒素 |

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4. Guidance for Industry: Drug Interaction Studies – Study Design, Data Analysis, and Implications for Dosing and Labeling (Draft, 2006)

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