

Dynabeads® Sheep anti-Rat IgG

For research use only

INDEX

1. PRODUCT DESCRIPTION

- 1.1 Intended Use
- 1.2 Principle of Isolation
- 1.3 Description of Materials

2. PROTOCOLS

- 2.1 Dynabeads Washing Procedure
- 2.2 Sample Preparation
- 2.3 Critical Steps for Cell Isolation
- 2.4 Cell Isolation – Indirect Technique
- 2.5 Cell Isolation – Direct Technique

3. TECHNICAL ADVICE

- 3.1 Indirect versus Direct Technique
- 3.2 Antibody Selection
- 3.3 Labelling Cells with Rat IgG Antibodies
- 3.4 Isolation and Depletion of Target Cells

4. GENERAL INFORMATION

5. REFERENCES

1. PRODUCT DESCRIPTION

1.1 Intended Use

Dynabeads Sheep anti-Rat IgG in combination with any primary rat IgG antibody are ideal for depletion or positive isolation of cells from different species (e.g. mouse) depending on the specificity of the primary antibody. Cells can be directly isolated from any sample such as whole blood, bone marrow, MNC suspensions or tissue digests.

1.2 Principle of Isolation

The primary rat IgG antibody is either added to the cell sample (indirect technique) or pre-coated onto the beads (direct technique) prior to cell isolation (see fig. 1).

Dynabeads are then mixed with the cell sample in a tube. The Dynabeads will bind to the target cells during a short incubation, and then the bead-bound cells are separated by a magnet.

- **Positive isolation** – discard the supernatant and use the bead-bound cells for downstream applications (e.g. molecular analysis or cell culture).
- **Depletion** – discard the bead-bound cells and use the remaining, untouched cells for any application.

1.3 Description of Materials

Dynabeads Sheep anti-Rat IgG are uniform, superparamagnetic polystyrene beads (4.5 µm diameter) coated with polyclonal Sheep anti-Rat IgG antibodies. Cross-reactivity to mouse antibodies is high and cross-reactivity to human antibodies is minimal.

Materials Supplied

- **5 ml Dynabeads Sheep anti-Rat IgG** 4 x 10⁸ beads/ml in phosphate buffered saline (PBS), pH 7.4, containing 0.1% bovine serum albumin (BSA) and 0.02% sodium azide (NaN₃).
- **This product will process up to 2 x 10⁹ cells**

Additional Material Required

- Magnet: (DynaL MPC™) - MPC-S for 20 µl to 2 ml samples, MPC-L for 1-8 ml samples, MPC-15 for 1-15 ml samples and MPC-50 for 15-50 ml samples.
- Mixer allowing both tilting and rotation.
- Buffer 1: PBS (without Ca²⁺ and Mg²⁺) w/0.1% BSA and 2 mM EDTA, pH 7.4.
- Rat IgG antibodies.

Important Notes:

BSA can be replaced by human serum albumin (HSA) or fetal calf serum (FCS).

EDTA can be replaced by sodium citrate.

PBS containing Ca²⁺ or Mg²⁺ is not recommended.

2. PROTOCOLS

2.1 Dynabeads Washing Procedure

Dynabeads should be washed before use. The washing procedure is facilitated by the use of a magnet.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer 1, or at least 1 ml, and mix.
4. Place the tube in a magnet for 1 min and discard the supernatant.
5. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads.

2.2 Sample Preparation

Cells can be directly isolated from any sample such

as whole blood, bone marrow, MNC suspensions or tissue digests.

Please visit www.invitrogen.com/samplepreparation for a list of recommended sample preparation procedures.

2.3 Critical Steps for Cell Isolation

- Use a mixer that provides tilting and rotation of the tubes to ensure Dynabeads do not settle at the bottom of the tube.
- When incubating Dynabeads and cells, the incubation temperature must be 2-8°C to reduce phagocytic activity and other metabolic processes.
- Never use less than 25 µl (1 x 10⁷) Dynabeads per ml cell sample and at least 4 Dynabeads per target cell.

Table 1: Volume of Dynabeads added per ml of cell sample. The volumes can be scaled up as required.

	Positive isolation	Depletion
Sample volume (1 x 10 ⁷ cells/ml*)	1 ml Max 2.5 x 10 ⁶ target cells	1 ml Max 2.5 x 10 ⁶ target cells
Volume of Dynabeads	25 µl	50 µl

* If the concentration of cells is increased or the target cell concentration exceeds 2.5 x 10⁶, the Dynabeads volume must be increased accordingly. Cell concentration can be up to 1 x 10⁸ cells per ml.

2.4 Cell Isolation – Indirect Technique

2.4.1 Labelling Cells with Rat IgG Antibodies

- Use approximately 1 µg antibody per 10⁶ target cells.
 - Recommended cell concentration: 1 x 10⁷ cells/ml.
1. Add rat IgG antibodies to the cell suspension and mix.
 2. Incubate for 20 min at 2-8°C.
 3. Wash the cells by adding 2 ml Buffer 1 per 1 x 10⁷ cells and centrifuge at 300 x g for 8 min. Discard the supernatant.
 4. Resuspend the cells in Buffer 1 at 1 x 10⁷ cells/ml.
 5. Proceed to Isolation or Depletion of Cells (2.4.2)

2.4.2. Isolation or Depletion of Cells

1. Add Dynabeads to the prepared sample according to table 1.
2. Incubate for 20 min (positive isolation) or 30 min (depletion) at 2 - 8°C with gentle tilting and rotation.
3. Double the volume with Buffer 1 to limit trapping of unbound cells (optional).
4. Place the tube in a magnet for 2 min.
5. **Depletion:** Transfer the supernatant containing the unbound cells to a fresh tube for further experiments.
6. **Positive isolation:** Discard the supernatant and gently wash the bead-bound cells 4 times, using the following procedure:
 - i. Add 1 ml Buffer 1 per 1 x 10⁷ Dynabeads.
 - ii. Place the tube in the magnet for 1 min and discard the supernatant.
7. Resuspend the cells in buffer/medium for downstream application.

2.5 Cell Isolation – Direct Technique

2.5.1 Pre-coating Dynabeads

- Use 0.5 – 1.5 µg rat IgG antibodies per 25 µl (1 x 10⁷) Dynabeads. (It is recommended to titrate the amount of antibody.)

1. Transfer washed Dynabeads to a tube.

2. Add rat IgG antibodies.

3. Incubate for ≥ 30 min at 2-8°C with gentle tilting and rotation.
4. Place the tube in a magnet for 1 min and discard the supernatant.
5. Wash the beads twice using 2 ml of Buffer 1.
6. Remove the tube from the magnet and resuspend the Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads.

2.5.2 Isolation or Depletion of Target Cells

1. Add the pre-coated Dynabeads to the cells according to table 1.
2. Incubate for 20 min (positive isolation) or 30 min (depletion) at 2 - 8°C with gentle tilting and rotation.
3. Double the volume with Buffer 1 to limit trapping of unbound cells (optional).
4. Place the tube in a magnet for 2 min.
5. **Depletion:** Transfer the supernatant containing the unbound cells to a fresh tube for further experiments.
6. **Positive isolation:** Discard the supernatant and gently wash the bead-bound cells 4 times, using the following procedure:
 - i. Add 1 ml Buffer 1 per 1 x 10⁷ Dynabeads.
 - ii. Place the tube in the magnet for 1 min and discard the supernatant.
7. Resuspend the cells in buffer/medium for downstream application.

3. TECHNICAL ADVICE

3.1 Indirect versus Direct Technique

Use the indirect technique when:

- A cocktail of rat IgG antibodies is used.
- Very high depletion efficiency is needed.
- The affinities of rat IgG antibodies are low.
- The cells express low number of target antigens.
- The direct technique gives unsatisfactory purity.

The direct technique may be used if:

- The affinity of the primary antibody is high.
- The cells express a high number of target antigens.
- A stock preparation of primary coated Dynabeads is desired.

3.2 Antibody Selection

The choice of primary antibody is the most important factor for successful cell isolation. Note that some antibodies may show reduced antigen-binding efficiency when coated onto beads (direct technique), even though the antibody shows good results in other immunological assays.

3.3 Labelling Cells with Rat IgG Antibodies

- To avoid non-specific binding of cells (e.g. monocytes, B cells), add aggregated IgG to block Fc receptors prior to adding the primary antibodies.
- Excess antibody must be removed by washing before cell isolation.

3.4 Isolation and Depletion of Target Cells

- Remove density gradient media (e.g. Ficoll): Wash cells prior to adding rat IgG antibodies or Dynabeads.
- Remove soluble factors in serum: Serum may contain soluble factors (e.g. antibodies or cell surface antigens), which can interfere with the cell isolation protocol. Washing the cells once may reduce this interference.




Fig. 1: Cell isolation using indirect or direct technique

4. GENERAL INFORMATION

Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

Storage/Stability

This product is stable until the expiry date stated on the label when stored unopened at 2-8°C.

Store opened vials at 2-8°C and avoid bacterial contamination.

Keep Dynabeads in liquid suspension during storage and all handling steps, as drying will result in reduced performance. Resuspend well before use.

Technical Support

Please contact Invitrogen Dynal for further technical information (see contact details).

Warning And Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

Follow appropriate laboratory guidelines.

This product contains 0.02% sodium azide as a preservative, which is cytotoxic. **Avoid pipetting by mouth!** Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. When disposing through plumbing drains, flush with large volumes of water to prevent azide build up.

Certificate of Analysis (CoA) is available upon request.

Material Safety Data Sheet (MSDS) is available at <http://www.invitrogen.com>.

Patents and Trademarks

Dynal®, Dynabeads® and Dynal MPC™ are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal, transportation prepaid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

Limited Use Label License

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,
Invitrogen Corporation,
1600 Faraday Avenue, Carlsbad,
California 92008.
Phone (760) 603-7200.
Fax (760) 602-6500.
Email: outlicensing@invitrogen.com

5. REFERENCES

1. Philippe D *et al.* (2003) Anti-inflammatory properties of the m opioid receptor support its use in the treatment of colon inflammation. *J.Clin.Invest.* 111:1329-1338.
2. Feng B *et al.* (2000) Development of mouse dendritic cells from lineage-negative c-kit^{low} pluripotent hemopoietic stem cells *in vitro*. *Stem Cells.* 18:53-60.
3. Zhou P *et al.* (2004) Expression of dual TCR on DO11.10 T cells allows for ovalbumin-induced oral tolerance to prevent T cell-mediated colitis directed against unrelated enteric antigens. *J. Immunol.* 172:1515-1523.
4. Yeh JRJ *et al.* (2000) The anti-angiogenic agent TNP-470 requires P53 and P21^{CIP/WAF} for endothelial cell growth arrest. *PNAS* 97(23): 12782-12787.
5. Na Nakorn T *et al.* (2002) Characterization of mouse clonogenic megakaryocyte progenitors. *PNAS* 100(1):205-210.
6. Steptoe RJ *et al.* (2003) Transfer of hematopoietic stem cells encoding auto-antigen prevents autoimmune diabetes. *J.Clin.Invest.* 111:1357-1363.
7. Tanchot C *et al.* (2004) Immune regulation by self-reactive T cells is antigen-specific. *J. Immunol.* 172:4285-4291.
8. Weber M *et al.* (2004) Minor histocompatibility antigens on canine hemopoietic progenitor cells. *J. Immunol.* 170:5861-5868.

9. Metwali A *et al.* (2004) Cutting Edge: Hemokinin has Substance P-like function and expression in inflammation. *J. Immunol.* 172: 6528-6532.
10. Tajima F *et al.* (2000) CD34 expression by murine hematopoietic stem cells mobilized by granulocyte colony-stimulating factor. *Blood* 96: 1989-1993.
11. El-Sukkari D *et al.* (2003) The protease inhibitor Cystatin C is differentially expressed among dendritic cell populations, but does not control antigen presentation. *J. Immunol.* 171:5003-5011.
12. Foussat A *et al.* (2003) A comparative study between T Regulatory Type I and CD4⁺CD25⁺ T cells in the control of inflammation. *J. Immunol.* 171: 5018-5026.
13. Joffre O *et al.* (2004) Induction of antigen-specific tolerance to bone marrow allografts with CD4⁺CD25⁺ T lymphocytes. *Blood.* 103: 4216-4221.
14. Hausl C *et al.* (2004) Preventing restimulation of memory B cells in hemophilia A: a potential new strategy for the treatment of antibody-dependent immune disorders. *Blood.* 104: 115-122.

Invitrogen Dynal is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2006 Invitrogen Dynal AS, Oslo, Norway. All rights reserved.

SPEC-05657