

Dynabeads[®] CD3/CD28

(Research version of Dynabeads[®] ClinExVivo[™] CD3/CD28)

For research use only.

INDEX

1. INTENDED USE
2. PRODUCT DESCRIPTION
 - 2.1 Principle of Isolation and Expansion
 - 2.2 Materials Supplied
 - 2.3 Storage and Stability
 - 2.4 Additional Materials Required
 - 2.5 Important Information
3. PROTOCOLS
 - 3.1 Dynabeads Washing Protocol
 - 3.2 Media Preparation
 - 3.3 Starting Material
 - 3.4 Magnetic Separation and Expansion of Cryopreserved CD3⁺ T cells
 - 3.5 Counting and Splitting of Cultures
 - 3.6 Culture Initiation Without Magnetic Concentration of CD3⁺ T Cells
 - 3.7 Separation and Expansion of Fresh (non-cryopreserved) Samples
 - 3.8 Procedures Incorporating Gene Transduction.
4. GENERAL INFORMATION
 - 4.1 Certifications
 - 4.2 Technical Service
 - 4.3 Precautions
5. WARNINGS AND LIMITATIONS
 - 5.1 Patents and Trademarks
 - 5.2 Intellectual Property Disclaimer
 - 5.3 Warranty
 - 5.4 Limited Use Label License
6. REFERENCES

1. INTENDED USE

Dynabeads CD3/CD28 are intended for separation and *in vitro* expansion of human T cells. This product is the research version of Dynabeads ClinExVivo CD3/CD28.

2. PRODUCT DESCRIPTION

Dynabeads CD3/CD28 are uniform 4.5 µm, superparamagnetic polystyrene beads coated with a mixture of monoclonal antibodies against the CD3 and CD28 cell surface molecules of human T cells.

2.1 Principle of Isolation and Expansion

Dynabeads CD3/CD28 offer a simple method for isolation and expansion of human T cells. Firstly, CD3⁺ T cells are separated and concentrated from the apheresis product by magnetic cell separation using Dynabeads CD3/CD28. Following separation, the CD3⁺ T cells are cultured in the presence of the beads. By combining anti-CD3 and anti-CD28 antibodies on Dynabeads, the beads will provide both the primary and co-stimulatory signals that are required for activation and expansion of T cells. Activated T cells produce IL2, GM-CSF, IFN-γ and TNF-α. T cells activated with Dynabeads CD3/CD28 can be expanded 100-1000 fold over a 9-14 day culture period.

2.2 Materials Supplied

• 10 ml of Dynabeads CD3/CD28

Supplied as a suspension containing 1x10⁸ Dynabeads/ml in phosphate buffered saline (PBS), pH 7.4, w/0.1% human serum albumin (HSA).

2.3 Storage and Stability

When stored in unopened vials at 2-8°C, Dynabeads CD3/CD28 are stable until the expiry date stated on the label.

Store opened vials at 2-8°C and avoid bacterial contamination.

Keep Dynabeads in liquid suspension during storage and all handling steps, as drying will result in reduced performance.

Resuspend well before use.

2.4 Additional Materials Required

Materials that are not included, but are recommended to perform the entire protocol:

- Dynabeads[®]M-450 Epoxy –(Cat. No. 140.11)
- Magnet: See www.invitrogen.com/magnets-selection for magnet recommendations
- PBS: Gibco, Cat. No. 20012-027 (Ca²⁺ and Mg²⁺ free)
- Pooled Human Serum (HS) heat inactivated at 56°C for 1 hour
- Human Serum Albumin (HSA)
- OpTmizer[™] T-Cell Expansion SFM (Gibco) serum free 1x formulation designed to support the culture and expansion of human T cells (or equivalent).
- RPMI-1640, Gibco, Cat. No. 12633-012
- L-glutamine, Gibco, Cat. No. 25030-081
- HEPES buffer, Gibco, Cat. No. 15630-080
- Proleukin[®] (IL-2), Chiron, or Biosource IL-2 Cat. No. PHC 0023.
- Appropriate sterile culture vessels.

2.5 Important Information

Follow standard precautions when working with human serum, plasma, or blood products. All human samples must be treated as a potential source of HIV, HBV, and other blood borne pathogens. Gloves and a lab coat must be worn when working with human samples.

Materials contaminated with blood products must be decontaminated by an approved chemical method and disposed of in labeled biohazard containers.

The protocol below is specifically written for the activation and expansion of human T cells from cryopreserved human PBMC obtained from leukapheresis products or Ficoll[®] separated whole blood, cord blood or bone marrow. Cultures may also be initiated from non-cryopreserved fresh samples. As sample sources and T cell or blood collection methods may vary, procedures may require specific modifications to maximize cell recovery, viability, activation and expansion. Such modifications must be determined empirically. Possible protocol variations that may be appropriate under certain experimental conditions are highlighted throughout the protocol, such as gene modification or working with particular disease tissues.

3. PROTOCOLS

Points to consider are listed in brackets ([...]) throughout the protocol. These identify areas where modifications should be considered for specific circumstances.

3.1 Dynabeads Washing Protocol

Dynabeads should be washed before use with the aid of a magnet.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer 1 as the initial volume of Dynabeads, or at least 1 ml, and mix.
4. Place the tube in a magnet for 1 minute until the Dynabeads are separated; discard the supernatant. Remove the tube from the magnet.
5. Resuspend the washed Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads.

3.2 Media Preparation

Buffer 1: PBS, pH 7.4, with 5% heat-inactivated HS or 1% HSA.

Incomplete Medium: OpTmizer[™] T-Cell Expansion SFM (Gibco) serum free 1x formulation designed to support the culture and expansion of human T cells (or equivalent).

Complete Medium: Add 200 IU/ml IL-2 to the incomplete medium.

Store at 2-8°C for a maximum of 7 days.

[Note: Autologous serum may be used, but expect donor to donor variation, IL-2 concentrations can be increased to enhance cell growth, although in static culture there is no significant increase in expansion for normal T cells when IL-2 levels are increased beyond 200-500 IU. Regulatory T cells may benefit from use of IL-2 at ≥500IU/ml. For IL-2 products where activity is not available in IU, typically 6-10ng/ml bioactive IL-2 is sufficient for optimal T cell expansion.]

3.3 Starting Material

The preferred starting material is cryopreserved human PBMC obtained from leukapheresis product or Ficoll[®] separated whole blood. The starting material can also be enriched for specific T cell subsets, such as CD4⁺ T cells or CD25⁺ T cells.

[The cryopreservation and subsequent thawing process facilitates lysis of granulocytes and other cell types that can suppress T cell activation and expansion. This may be particularly important when working with samples from patients in certain disease states with elevated granulocyte counts.]

Alternatively, PBMC from freshly obtained leukapheresis products or Ficoll separated whole blood may be activated and expanded without cryopreservation. If fresh samples are used, deplete monocytes with Dynabeads M-450 Epoxy (see Section 3.7)

For maximum activation and expansion of T cells in PBMC magnetically capture CD3⁺ T cells prior to culture initiation (see Section 3.4). Magnetic concentration is not required if T cells or T cell subsets have been enriched prior to activation and expansion (e.g. purified CD4⁺ T cells).

3.4 Magnetic Separation and Expansion of Cryopreserved CD3⁺ T cells

For fresh samples, see Section 3.7

[Note: Magnetic separation serves two purposes. Firstly, the process tends to synchronize the T cell activation process, accelerating contact between Dynabeads and T cells, as well as promoting early T cell aggregation. Secondly, during the magnetic selection process, T cells are enriched, leaving behind monocytes in the unselected fraction with other cells that might inhibit T cell activation and expansion. Most notably, monocytes can rapidly phagocytose beads at 37°C, thereby reducing the absolute number of beads capable of contacting T cells, in turn reducing the level of T cell activation and expansion.]

1. Determine the percentage of CD3⁺ T cells in the sample by flow cytometry or other suitable methods.
2. For PBMC, resuspend cells in a tube at roughly 2-5 x 10⁷ CD3⁺ cells/ml in Buffer 1, but not exceeding a maximum of 2 x 10⁸ total nucleated cells/ml. For previously enriched/purified T cells (e.g. CD4⁺), resuspend the cells at 5-10 x 10⁷ cells/ml.
3. Add washed Dynabeads CD3/CD28 at a 3:1 bead: CD3⁺ T cell ratio

[Note: Bead:CD3⁺ T cell ratio can be dropped as low as 1:1 while still achieving nearly equivalent T cell capture. Lower capture/activation initiation ratios, such as 1:1, may be preferred when the starting T cell pool is more prone to activation induced cell death or where lower levels of T cell expansion are desired.]

4. Rotate the sample at 1 rpm for 30 min at between 4-25°C. Optimise the mixing temperature between 4-25°C for each application.

[If starting PBMC have relatively low CD3 content, e.g. less than 25%, it may be beneficial to rotate for 1-2 hours at 4-25°C in Incomplete Medium or Buffer 1.]

5. Increase the volume in the tube to sufficient levels for magnetic selection by adding Incomplete Medium or Buffer 1 before placing the tube in a Dynal MPC for 1-2 min.

6. Remove supernatant and stain with anti-CD3 antibody for flow cytometric analysis to calculate depletion efficiency.

7. Gently resuspend bead:cell complexes to an estimated 0.5 x 10⁶ -1 x10⁶ CD3⁺ T cells/ml in Complete Medium (assuming 100% capture of CD3⁺ T cells from PBMC). Typically, 0.5 x 10⁶ CD3⁺ T cells/ml are optimal for culture initiation. Depending upon the particular application (e.g. gene modification), cultures can be initiated with seeding densities up to 1 x 10⁶ CD3⁺ T cells/ml.

8. Plate cells at 2.5-5 x 10⁶ CD3⁺ cells in a total volume of 5 ml/well (0.5-1.0 x 10⁶ cells/ml) of a 6 well culture plate in a humidified incubator at 37°C/ 5% CO₂ for 3 days.

[Note: smaller or larger culture vessels/wells can be used - adjust volumes accordingly. For example, in a 24 well format, total volume should be approximately 1 ml/well; in a 12 well format, total volume should be approximately 2 ml/well. For flasks, culture volumes may vary, however it is important that static cultures do not exceed ~1.2 cm in volume depth, otherwise expansion potential and viability may be compromised due to reduced gas exchange.]

3.5 Counting and Splitting of Cultures

Evaluate cell concentration daily, beginning on day 3 of culture. It may be useful to conduct regular phenotypic analyses by flow cytometry (e.g. days 3, 5, 8, and end of culture) to track T cell purity, phenotype and activation state. Similarly, supernatants may be collected to measure cytokine secretion patterns to further characterize T cell activity.

1. Gently mix contents of wells to dissociate beads and cells. Mixing should be sufficient to disrupt all visible bead:cell complexes.
 2. Remove 50 µl sample and mix in 50 µl trypan blue (do not remove beads before counting).
 3. Count cells manually using a hemocytometer.
- [Caution: Insufficient mixing of bead:cell complexes may result in cell count underestimates as they will not migrate efficiently under hemocytometer coverslips.]*
4. When CD3⁺ T cell density is >1 x 10⁶ cells/ml, dilute cells to approximately 0.5 x 10⁶ CD3⁺ T cells/ml in culture medium.

[Note: The rate of cell growth tends to diminish as cell concentrations exceed 1-2 x 10⁶ CD3⁺/ml]

5. At the end of culture (day 9-14) count cells and remove beads with a magnet.

3.6 Culture Initiation Without Magnetic Concentration of CD3+ T Cells

1. Add 5×10^6 enriched or purified T cells (e.g. purified CD4⁺ T cells) to 1.5×10^7 washed Dynabeads CD3/CD28 in a final volume of 10 ml culture media, placing 5 ml/well in a 6 well tissue culture plate to give a final concentration of 0.5×10^6 CD3⁺ T cells/ml. Mix gently.
2. Culture cells in a humidified incubator at 37°C/5% CO₂ for 3 days
3. Cell concentration should be evaluated daily beginning on day 3 of culture.
4. Counting of cells and splitting of culture; see Section 3.5.

[As described earlier, cells can be cultured in smaller or larger culture vessels by adjusting cell numbers and culture volumes accordingly.]

3.7 Separation and Expansion of Fresh (non-cryopreserved) Samples

[Note: During the thawing process of cryopreserved PBMC, contaminating granulocytes tend to lyse, thereby diminishing any inhibitory effects upon subsequent T cell expansion. If working with fresh, non-cryopreserved samples, take into consideration the level of contaminating granulocytes when designing a culture initiation protocol.

For initiation of cultures with non-cryopreserved samples, use the Magnetic Concentration method described above (see Section 3.4).

Overall T cell activation and expansion can be further improved by specifically depleting phagocytic cells, such as monocytes, using Dynabeads M-450 Epoxy prior to magnetic concentration and culture initiation. This is particularly useful when monocyte levels are >15% of total PBMC.]

1. Dilute PBMC to $5-10 \times 10^6$ cells/ml in Incomplete Medium (pre-warmed to 37°C).
2. Wash Dynabeads M450 Epoxy as described above for Dynabeads CD3/CD28 (see Section 3.1).
3. Add Dynabeads M-450 Epoxy to the PBMC at 1-2 beads per total nucleated cells in a sterile culture vessel. 2 beads per nucleated cell can be used when monocyte levels are >15% of total PBMC. Use any sterile culture vessel, including bags, tissue culture flasks, tissue culture plates or sterile tubes.
4. Place the culture vessel containing beads and cells into a humidified incubator at 37°C/ 5% CO₂ for one hour. If using a culture flask, lay the flask on its side to maximize culture surface area.
5. After 1 hour, remove Dynabeads M-450 Epoxy and any cells that have ingested beads by applying an appropriate DYNAL MPC for 1-2 min.
6. Collect non-selected cells and initiate culture as described in Section 3.4 or 3.6 above, depending upon whether you want to incorporate the magnetic concentration step. For maximal activation and expansion, the magnetic concentration method is recommended. Note: If magnetic concentration is planned after the monocyte/phagocytic cell removal step, keep cells in Incomplete Medium throughout magnetic concentration instead of using Buffer 1.
7. For larger volumes, where PBMC and beads are in bags, beads and bead:cell complexes can be removed on the DYNAL ClinExVivo™ MPC™ (Cat. No. 121-02D). The magnet angle can vary between 0° and 60° and should be determined empirically.
8. Transfer the non-captured cell fraction containing monocyte-depleted T cells to an appropriate culture vessel (plate, flask, or bag).
9. Count cells and measure CD3⁺ T cell content of sample by flow cytometry.
10. Adjust the concentration of cells to $1-5 \times 10^7$ CD3⁺ T cells/ml in Incomplete Medium.

11. Capture CD3⁺ T cells as described above in Section 3.4.
12. Resuspend captured T cells and Dynabeads CD3/CD28 to approximately $0.5-1 \times 10^6$ CD3⁺ T cells/ml in Complete Medium containing IL-2 and transfer to an appropriate culture vessel. Incubate at 37°C/5% CO₂ until Day 3 of culture.
13. Count and split cell culture as described above in Section 3.5.

3.8 Procedures Incorporating Gene Transduction.

Typically, for all culture conditions described earlier, T cells from normal donor tissues begin cycling and start to divide between day 2 and 3 of culture. Days 2, 3 and/or 4 of culture are recommended as optimal days for transduction using moloney-based vectors, whereas day 1, 2 and/or 3 are recommended for lentivirus based vectors. Magnetic removal of beads prior to transduction will diminish overall cell expansion, but should not affect viability. Leaving beads in during the retroviral transduction process should be acceptable for most transduction applications.

[Note: T cells obtained from patients with various diseases and/or undergoing various treatments may be slower to enter cell cycle and cell division may not commence until 1, 2 or even 3 days later than typically observed for tissues from healthy donors. T cells from patients with HIV infection, tissues from patients undergoing chemotherapy, or patients with certain kinds of cancer (e.g. chronic lymphocytic leukemia) may be slower to start cell cycling. Thus it is important to monitor T cell activation markers, such as CD25, as well as cell division to determine optimal splitting schedules and timing for gene modification.]

4. GENERAL INFORMATION

4.1 Certifications

Invitrogen DYNAL AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

4.2 Technical Service

Contact details for further technical information can be found at <http://www.invitrogen.com/contact>.

Certificate of Analysis/Compliance is available upon request.

4.3 Precautions

Material Safety Data Sheet is available upon request.

5. WARNINGS AND LIMITATIONS

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

If mammalian cells are used in this procedure, appropriate laboratory guidelines must be followed. Detailed information on such guidelines can be obtained from: Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body fluids and Tissue: Approved Guideline: M29-A; ISBN 1-56238-339-6; <http://www.nccsl.org>.

Handle all samples as if capable of transmitting disease. All work should be performed wearing gloves and appropriate protection.

5.1 Patents and Trademarks

DYNAL®, Dynabeads®, DETACHaBEAD® and DYNAL MPC™ are either registered trademarks or trademarks of Invitrogen DYNAL AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

ClinExVivo™ is a trademark of Invitrogen DYNAL AS, Oslo, Norway.

Ficoll® is a registered trademark of GE Healthcare Bio-Sciences

Proleukin® is a registered trademark of Chiron Corporation

OpTmizer™ T-cell Expansion SFM is a trademark of Invitrogen Corporation

5.2 Intellectual Property Disclaimer

Invitrogen DYNAL will not be responsible for violations or patent infringements that may occur with the use of our products.

Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen DYNAL's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen DYNAL's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen DYNAL, transportation prepaid, or at Invitrogen DYNAL's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen DYNAL, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen DYNAL's maximum liability is limited in all events to the price of the products sold by Invitrogen DYNAL. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the

purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund.

No. 305: Invitrogen T Cell Expansion Technology – The use of this product for *ex vivo* activation or expansion of human T-cells is covered under various patents owned and licensed from the U.S. Navy, the University of Michigan and the Dana Farber Cancer Institute and is sold for research use only. This product is not to be used for activation or expansion of T-cells modified through gene transfer to specifically modify the T-cells to produce secreted or cell-surface membrane-bound proteins not normally expressed in significant levels by such T-cells, unless the proteins directly enable the selection, or directly modify or preserve the function, of the T cells. Modifications that are acceptable would include, for example, introducing genes that direct expression of cell-surface bound proteins that bind the T-cell to specific target cells.

For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,
Invitrogen Corporation,
1600 Faraday Avenue, Carlsbad,
California 92008.
Phone (760) 603-7200.
Fax (760) 602-6500.
Email: outlicensing@invitrogen.com

6. REFERENCES

1. Kalamasz D *et al* (2004) Optimization of human T-cell expansion *ex vivo* using magnetic beads conjugated with anti-CD3 and Anti-CD28 antibodies. *J Immunother.* 27(5):405-418.
2. Bonyhadi M *et al* (2005) *In vitro* engagement of CD3 and CD28 corrects T cell defects in chronic lymphocytic leukemia. *J Immunol.* 174(4):2366-2375.
3. Earle KE *et al* (2005) *In vitro* expanded human CD4⁺CD25⁺ regulatory T cells suppress effector T cell proliferation. *Clin Immunol.* 115(1):3-9.
4. Bondanza A *et al* (2006) Suicide gene therapy of graft-versus-host disease induced by central memory human T lymphocytes. *Blood* 107(5): 1828-36.
5. Noonan K *et al* (2005) Activated marrow-infiltrating lymphocytes effectively target plasma cells and their clonogenic precursors. *Cancer Res.* 65(5):2026-2034.

Invitrogen DYNAL is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2008 Invitrogen DYNAL AS, Oslo, Norway. All rights reserved.

SPEC-06245