

Dynabeads[®] CD4

For research use only.

1. PRODUCT DESCRIPTION

- 1.1 Intended Use
- 1.2 Principle of Isolation
- 1.3 Description of Materials

2. PROTOCOLS

- 2.1 Dynabeads Washing Procedure
- 2.2 Sample Preparation
- 2.3 Critical Steps for Cell Isolation
- 2.4 Depletion or Positive Isolation of CD4⁺ T Cells
- 2.5 Detachment of Dynabeads from Isolated Cells

3. GENERAL INFORMATION

4. REFERENCES

1. PRODUCT DESCRIPTION

1.1 Intended Use

Isolate or deplete human CD4⁺ T cells directly from whole blood, buffy coat or MNC with Dynabeads CD4. For rapid and consistent results in protein or gene expression analysis, lyse the CD4⁺ T cells while they are still attached to the beads and directly process for further molecular analysis. For downstream cell-based applications (requiring bead-free CD4⁺ cells), Invitrogen

Dynal recommends the Dynabeads[®] FlowComp[™] Human CD4 kit (cat. no. 113.61D).

1.2 Principle of Isolation

Dynabeads are mixed with the cell sample in a tube. The Dynabeads will bind to the target cells during a short incubation, and then the bead-bound cells are separated by a magnet.

- **Positive isolation** – discard the supernatant and use the bead-bound cells for downstream molecular applications.
- **Depletion** – discard the bead-bound cells and use the remaining, untouched cells for any application.

1.3 Description of Materials

Dynabeads CD4 are uniform, superparamagnetic polystyrene beads (4.5 μm diameter) coated with a primary monoclonal antibody specific for the CD4 membrane antigen on human cells.

The primary CD4 antibody is coupled to the Dynabeads via a secondary human IgG4 anti-mouse IgG antibody. The source of the human monoclonal antibody is free of Human Immunodeficiency Virus (HIV), Hepatitis-B Virus (HBV) and Hepatitis-C Virus (HCV).

Materials Supplied

- **5 ml Dynabeads CD4.**
4 x 10⁸ beads/ml in phosphate buffered saline (PBS), pH 7.4, containing 0.1% bovine serum albumin (BSA) and 0.02% sodium azide (NaN₃).
- **This product will process up to 2 x 10⁹ cells.**

Additional Materials Needed

- Magnet (Dynal MPC[™]): See www.invitrogen.com/magnets-selection for magnet recommendations.
- Mixer allowing both tilting and rotation.
- Buffer 1: PBS w/0.1% BSA and 2 mM EDTA, pH 7.4.

Important Notes:

BSA can be replaced by human serum albumin (HSA) or fetal calf serum (FCS).

EDTA can be replaced by sodium citrate. PBS containing Ca²⁺ or Mg²⁺ is not recommended.

2. PROTOCOLS

2.1 Dynabeads Washing Procedure

Dynabeads should be washed before use.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer 1, or at least 1 ml, and mix.
4. Place the tube in a magnet for 1 min and discard the supernatant.
5. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads (step 2).

2.2 Sample Preparation

Cells can be directly isolated from any sample such as whole blood, bone marrow, MNC or tissue digests.

Please visit www.invitrogen.com/cell-isolation and follow our QuickLinks for recommended sample preparation procedures.

Whole Blood and Buffy Coat

Most depletions and positive isolations can use whole blood and buffy coat as a starting sample. Buffy coat is 8-10 times more concentrated than whole blood with regard to number of leucocytes. For this product you have to wash the blood/buffy coat to remove interfering soluble factors.

1. Dilute the whole blood or buffy coat in Buffer 1 (1+2).
2. Centrifuge at 600 x g for 10 min at 2-8°C.
3. Discard the plasma fraction/upper layer. Resuspend blood to the original volume in Buffer 1 and buffy coat 1+1 in Buffer 1 before adding the beads.

2.3 Critical Steps for Cell Isolation

- Use a mixer that provides tilting and rotation of the tubes to ensure Dynabeads do not settle at the bottom of the tube.
- When incubating Dynabeads and cells, the incubation temperature must be 2-8°C to reduce phagocytic activity and other metabolic processes.
- Never use less than 25 μl Dynabeads per ml of cell sample.

Table 1: Volume of Dynabeads added per 10⁷ cells. The volumes can be scaled up as required.

	Positive isolation	Depletion
Sample volume - 1 x 10 ⁷ MNC/ml - washed whole blood - washed buffy coat	1 ml	1 ml
Vol. of Dynabeads per 10 ⁷ cells	25 μl	50 μl
Total no. of cells processed per product	2 x 10 ⁹ cells	1 x 10 ⁹ cells

2.4 Depletion or Positive Isolation of CD4⁺ T Cells

1. Add the appropriate volume of Dynabeads to the prepared sample according to table 1.
2. Incubate for 20 min (positive isolation) or 30 min (depletion) at 2 - 8°C with gentle tilting and rotation.
3. Place the tube in a magnet for 2 min.
4. For depletion, transfer supernatant to a new tube for further use.
5. For positive isolation, discard the supernatant and wash the bead-bound cells 3 times by resuspending in Buffer 1 to the original sample volume, and separate using a magnet for 1 min. Never use less than 1 ml Buffer 1 in each washing step.

For molecular studies, lyse cells while still attached to the beads and transfer supernatant to a new tube for protein or gene expression analysis.

3. GENERAL INFORMATION

Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

Storage/Stability

This product is stable until the expiry date stated on the label when stored unopened at 2-8°C.

Store opened vials at 2-8°C and avoid bacterial contamination.

Keep Dynabeads in liquid suspension during storage and all handling steps, as drying will result in reduced performance. Resuspend well before use.

Technical Support

Please contact Invitrogen Dynal for further technical information (see contact details).

Warning And Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.




Fig. 1 Overview of method

Follow appropriate laboratory guidelines. This product contains 0.02% sodium azide as a preservative, which is cytotoxic. **Avoid pipetting by mouth!** Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. When disposing through plumbing drains, flush with large volumes of water to prevent azide build up.

Certificate of Analysis (CoA) is available upon request.

Material Safety Data Sheet (MSDS) is available at <http://www.invitrogen.com>.

Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

Patents and Trademarks

Dynal[®], Dynabeads[®], Dynal MPC[™] and FlowComp[™] are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its compo-

nents, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: outlicensing@invitrogen.com

Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal, transportation prepaid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

4. REFERENCES

1. Roscic-Mrkic B et al (2003) RANTES (CCL5) uses the proteoglycan CD44 as an auxiliary receptor to mediate cellular activation signals and HIV-1 enhancement. *Blood* 102: 1169-1177.
2. Peng Koh K et al (2004) T cell-mediated vascular dysfunction of human allografts results from IFN-g dysregulation of NO synthase. *J. Clin. Investigation*. 114:846-856.
3. Dai KZ et al (2004) Transcriptional activation of the SH2D2A gene is dependent on a cyclic adenosine 5'-monophosphate-responsive element in the proximal SH2D2A promoter. *J. Immunol*. 172: 6144-6151.

Invitrogen Dynal is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2008 Invitrogen Dynal AS, Oslo, Norway. All rights reserved.

SPEC-06045