

Dynabeads® CD3/CD28 T Cell Expander

For research use only

Index

1. PRODUCT DESCRIPTION

- 1.1 Intended Use
- 1.2 Principle of Activation and Expansion
- 1.3 Description of Materials

2. PROTOCOLS

- 2.1 Technical Advice
- 2.2 Dynabeads Washing Procedure
- 2.3 Isolation of Human T Cells
- 2.4 Short Term Stimulation of Human T Cells (1-3 days)
- 2.5 Polyclonal Expansion of T Cells
- 2.6 Expansion of T Cell Clones

3. GENERAL INFORMATION

4. REFERENCES

1 PRODUCT DESCRIPTION

1.1 Intended Use

This product is intended for activation and expansion of human T cells and T cell clones. T cell

clones (CD28⁺) retain specificity and functional properties after repeated re-stimulation with Dynabeads® CD3/CD28 T Cell Expander.

1.2 Principle of Activation and Expansion

Dynabeads CD3/CD28 T Cell Expander offers a simple method for activation and expansion of T cells that does not require antigen-presenting cells or antigen. Just add the Dynabeads for activation or Dynabeads plus recombinant IL-2 (rIL-2) for expansion of T cells. Cell cultures showing signs of exhaustion can be re-stimulated by adding fresh beads and rIL-2.

Dynabeads CD3/CD28 T Cell Expander offer the first artificial antigen-presenting cells to provide simultaneous signals to TCR/CD3 and CD28 for full activation and expansion of human T cells.

1.3 Description of Materials

Dynabeads CD3/CD28 T Cell Expander are uniform 4.5 µm, superparamagnetic polystyrene beads coated with an optimized mixture of monoclonal antibodies against the CD3 and CD28 cell surface molecules of human T cells.

The CD3 antibody coated on the CD3/CD28 T Cell Expander is specific for the epsilon chain of human CD3, a subunit of the TCR complex. The CD28 antibody is specific for the human CD28 co-stimulatory molecule, which is the receptor for CD80 (B7-1) and CD86 (B7-2). Both antibodies are coupled to the same Dynabead, mimicking *in vivo* stimulation by antigen presenting cells.

The plasma used for human serum albumin (HSA) in this product is tested and found free of Human Immunodeficiency Virus (HIV-1 and HIV-2), Hepatitis-B Virus (HBV) and Hepatitis-C Virus (HCV).

Materials Supplied

2 ml or 5 x 2 ml Dynabeads CD3/CD28 T Cell Expander:

4 x 10⁷ beads/ml in phosphate buffered saline (PBS), pH 7.4, w/0.1% human serum albumin (HSA).

Additional Materials Required

Materials that are not included, but are needed to perform the entire protocol:

- Magnet (DynaL MPC™): See www.invitrogen.com/magnets-selection for magnet recommendations.
- Buffer 1: PBS w/ 0.1% BSA, pH 7.4.
- Culture Medium: OpTmizer™ T-Cell Expansion SFM (Gibco)serum free 1x formulation designed to support the culture and expansion of human T cells (or equivalent).
- rIL-2.

2. PROTOCOLS

2.1 Technical Advice

- Do not use excess volume of Dynabeads as extra beads per well will inhibit cell growth.

2.2 Dynabeads Washing Procedure

Dynabeads should be washed before use.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer 1, or at least 1 ml, and mix.
4. Place the tube in a magnet for 1 minute and discard the supernatant.

5. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads (step 2).

2.3 Isolation of Human T Cells

For isolation of human T cells it is recommended to use the Dynal T Cell Negative Isolation Kit (Cat. no. 113.11D), Dynal CD4 Negative Isolation Kit (Cat. no. 113.17D), Dynal CD8 Negative Isolation Kit (Cat. no. 113.19D), Dynal CD4 Positive Isolation Kit (Cat. no. 113.31D) or Dynal CD8 Positive Isolation Kit (Cat. no. 113.33D).

2.4 Short Term Stimulation of Human T Cells (1-3 days)

For short term culture of T cells for 1 - 3 days Invitrogen Dynal recommends using 1 bead per T cell. Short-term culture can be performed without adding exogenous rIL-2.

MNC or T cells can be used for short term stimulation with Dynabeads CD3/CD28 T Cell Expander.

Day 0:

1. Add MNC or T cells to the culture medium at 1 x 10⁶ cells/ml.
2. Add an equal number of Dynabeads CD3/CD28 T Cell Expander to the cells, 25 µl beads to 1 x 10⁶ T cells.
3. Leave for 1-3 days in an incubator.

2.5 Polyclonal Expansion of T Cells

For expansion of T cells for >3 days Dynal recommends using 3 beads per T cell.

MNC or T cells can be used for expansion with Dynabeads CD3/CD28 T Cell Expander.

Day 0:

1. Add MNC or T cells to the culture medium at 1 x 10⁶ cells/ml.
2. Add 75 µl Dynabeads CD3/CD28 T Cell Expander to 1 x 10⁶ cells.
3. Leave in an incubator.

Day 2-3:

- Add 100-100 U rIL-2/ml.

Day 4-5:

1. Resuspend the cells and beads with a pipette and count the cells in a microscope.
2. Dilute the cells to 0.5 x 10⁶ cells/ml in culture medium containing 10-100 U rIL-2/ml.

Day 7-14:

Split the culture when needed. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures. Count the cells at least twice weekly after thorough resuspension. When the cell density exceeds 2 x 10⁶ cells/ml or when the medium becomes yellow, split cultures to a density of 0.5 x 10⁶ cells/ml in culture medium containing 10-100 U rIL-2/ml.

Day 14: Cells can usually be expanded for up to 2 weeks before restimulation is needed.

Re-stimulation:

1. Resuspend cells and beads before removing beads with a magnet.
2. Count the cells and dilute cells to 0.5 x 10⁶ cells/ml in culture medium containing 10-100 U rIL-2/ml.
3. Re-stimulate by adding new Dynabeads CD3/CD28 T Cell Expander, 3 beads per T cell.
4. Culture cells for 2 - 3 days, count and dilute cells as described above.
5. Re-stimulation with Dynabeads CD3/CD28 T Cell Expander can be performed several times, as long as the T cells express CD28.

2.6 Expansion of T Cell Clones

Day 0:

1. Transfer CD4⁺CD28⁺ or CD8⁺CD28⁺ T cell clones from a full Terasaki well to a 96 well plate.
2. Add 0.125 µl Dynabeads CD3/CD28 T Cell Expander and 10-100 U rIL-2/ml per well. (1 bead per 5 - 10 T cells, higher number of beads will reduce expansion of T cell clones).

3. Leave in an incubator.

Day 5-9: Grow the clones until the well is half-full (~500,000 cells). If wells turn yellow during culture, add 50-100 µl medium with 10-100 U rIL-2/ml. If higher cell numbers are needed, continue with restimulation as described below.

Re-stimulation:

Day 0:

1. Resuspend cells and beads before removing beads with a magnet.
2. Count the cells and dilute cells to 0.25 - 0.5 x 10⁶/ml in culture medium containing 10-100 U rIL-2/ml. Transfer the cloned cells from a 96 well plate to a 24 well plate. If estimated cell number is below 0.25 x 10⁶ cells transfer to a 48 well plate.
3. Re-stimulate by adding 1.25 µl Dynabeads CD3/CD28 T Cell Expander per 24 well (1 bead per 5-10 T cells) or 0.63 µl Dynabeads to a 48 well plate.
4. Leave in an incubator.

Day 3-14:

Split the culture when needed. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures. Count the cells at least twice weekly after thorough resuspension. When the cell density exceeds 2 x 10⁶ cells/ml or when the medium becomes yellow, split cultures to a density of 0.5 x 10⁶ cells/ml in culture medium containing 10-100 U rIL-2/ml.

Further re-stimulation can be performed with fresh Dynabeads CD3/CD28 T Cell Expander as described above.

3. GENERAL INFORMATION

Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

Storage & Precautions

This kit is stable until the expiry date stated on the label, when stored unopened at 2-8°C

Store opened vials at 2-8°C and avoid bacterial contamination

Keep Dynabeads in liquid suspension during storage and all handling steps, since drying will result in reduced performance. Resuspend well before use.

Technical Service

Please contact Invitrogen Dynal for further technical information at www.invitrogen.com.

Certificate of Analysis (CoA) is available upon request.


Warnings and Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

Material Safety Data Sheet (MSDS) is available at <http://www.invitrogen.com>.

Patents and Trademarks

Dynal®, Dynabeads® and Dynal MPC™ are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries. OpTmizer™ T-cell Expansion SFM is a trademark of Invitrogen Corporation.



Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal, transportation prepaid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,
Invitrogen Corporation,
1600 Faraday Avenue, Carlsbad,
California 92008.
Phone (760) 603-7200.
Fax (760) 602-6500.

Email: outlicensing@invitrogen.com

4. REFERENCES

1. Chung HK *et al* (2003) Activation of interleukin-13 expression in T cells from HTLV-1-infected individuals and in chronically infected cell lines. *Blood*. 102: 4130-4136.
2. Danese S *et al* (2004) Cutting Edge: T cells trigger CD40-dependent platelet activation and granular RANTES release: A novel pathway for immune response amplification. *J. Immunol*. 172:2011-2015.
3. Garlie NK *et al* (1999) T cells co-activated with immobilized anti-CD3 and anti-CD28 as potential immunotherapy for cancer. *J. Immunother*. 22: 336-345.
4. Godfrey WR *et al* (2004) *In vitro*-expanded human CD4⁺CD25⁺ T regulatory cells can markedly inhibit allogeneic dendritic cell-stimulated MLR cultures. *Blood*. 104:453-461.
5. Godfrey WR *et al* (2005) Cord blood CD4⁺CD25⁺-derived T regulatory cell lines express FoxP3 protein and manifest potent suppressor function. *Blood*: 105:750-8.
6. Hoffmann P *et al* (2004) Large-scale *in vitro* expansion of polyclonal human CD4⁺CD25^{high} regulatory T cells. *Blood*. 104:895-903.
7. Ishiguro K and Xavier R(2004) Homer-3 regulates activation of serum response element in T cells via its EVH1 domain. *Blood* 103:2248-2256.
8. Jacobsen M *et al* (2002) Oligoclonal expansion of memory CD8⁺ T cells in cerebrospinal fluid from multiple sclerosis patients. *Brain* 125: 538-550.
9. Trickett A *et al* (2003) T cell stimulation and expansion using anti-CD3/CD28 beads. *J. Imm. Methods* 275: 251-255.
10. Zhang R *et al* (2004) CD40 ligand dysregulation in HIV infection: HIV glycoprotein 120 inhibits signaling cascades upstream of CD40 ligand transcription. *J. Immunol*. 172: 2678-2686.
11. von Essen M *et al* (2004) Constitutive and ligand-induced TCR degradation. *J. Immunol*. 173:384-393.

Invitrogen Dynal is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2008 Invitrogen Dynal AS, Oslo, Norway.
All rights reserved.

SPEC-05780