

Keep the buffers cold!

Note:

BSA can be replaced by human serum albumin (HSA) or 2% FBS/FCS.

 EDTA can be replaced by sodium citrate. PBS containing Ca²⁺ or Mg²⁺ is not recommended.

2. PROTOCOL
2.1 Dynabeads Washing Procedure

Dynabeads should be washed before use.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer, or at least 1 ml, and mix.
4. Place the tube in a magnet for 1 min and discard the supernatant.
5. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer as the initial volume of Dynabeads (step 2).

2.2 Sample Preparation

 Prepare MNC from buffy coat. Please visit www.invitrogen.com/cellisolation and follow our QuickLinks for recommended sample preparation procedures.

2.3 Critical Steps For Cell Isolation

- Use a mixer that provides tilting and rotation of the tubes to ensure Dynabeads do not settle at the bottom of the tube.
- Never use less than 100 µl Dynabeads per 1 x 10⁷ MNC sample.
- Dynal MPC[™]-1 (Cat. no. 120.01D) should not be used with this product.

2.4 Isolation of Human B Cells from MNC

 This protocol is based on 1 x 10⁷ MNC. It is scalable from 1 x 10⁷-5 x 10⁸ cells, (see table 1).

1. Transfer 100 µl (1 x 10⁷) MNC in Buffer to a tube.
2. Add 20 µl of Antibody Mix.
3. Mix well and incubate for 20 min. at 2-8°C.
4. Wash the cells by adding 2 ml Buffer. Mix well by tilting the tube several times and centrifuge at 300 x g for 8 min at 2-8°C. Discard the supernatant.
5. Resuspend the cells in 800 µl Buffer.
6. Add 200 µl pre-washed Depletion Dynabeads.

7. Incubate for 15 min at 18-25°C with gentle tilting and rotation.
8. Resuspend the bead-bound cells by **gently** pipetting 5 times using a pipette with a narrow tip opening, (e.g. a 1000 µl pipette tip or a 5 ml serological pipette).
9. Add 1 ml Buffer and mix.
10. Place the tube in the magnet for 2 min.
11. Transfer the supernatant containing the negatively isolated B cells to a new tube.
12. Repeat steps 9-11 and pool the two supernatants.

 Table 1. Volumes for human B cell isolation (1 x 10⁷ to 5 x 10⁸ MNC)

2.5 Downstream Applications

	Working volume per 1 x 10 ⁷ MNC
Cell volume (step 1)	100 µl
Antibody Mix (step 2)	20 µl
Washing (step 4)	2 ml (Max. 40 ml when scaled up)
Resuspension (step 5)	800 µl (Max. 5 ml when scaled up)
Depletion Dynabeads (step 6)	200 µl
Increase volume (step 9)	1 ml (Max. 35 ml when scaled up)

Isolated B cells can be used in applications such as:

- Flow cytometry (for recommended products and protocols visit www.invitrogen.com/immunology.)
- Antigen presentation by B cells and interaction with other cells of the immune system.
- Analysis of B cell immunoglobulin class switching and somatic hypermutation.
- Analysis of B cell activation, proliferation and differentiation.
- B cell signalling pathway studies.

3. GENERAL INFORMATION

Manufactured by Invitrogen Dynal AS. Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

Description of Materials

Depletion Dynabeads are uniform, superparamagnetic polymer beads (4.5 µm diameter) coated with a monoclonal human anti-mouse IgG antibody. The antibody coated onto Depletion Dynabeads recognises all mouse IgG subclasses and is Fc-specific. The antibody on the Dynabeads is a human IgG4 anti-mouse IgG. The source of the human monoclonal antibody is free of

Dynal[®] B Cell Negative Isolation Kit

For research use only.
1. PRODUCT DESCRIPTION

- 1.1 Intended Use
- 1.2 Principle of Isolation
- 1.3 Materials Supplied
- 1.4 Additional Materials Required

2. PROTOCOL

- 2.1 Dynabeads Washing Procedure
- 2.2 Sample Preparation
- 2.3 Critical Steps For Cell Isolation
- 2.4 Isolation of Human B Cells from MNC
- 2.5 Downstream Applications

3. GENERAL INFORMATION
4. REFERENCES
1. PRODUCT DESCRIPTION
1.1 Intended Use

Isolate untouched human naïve B cells by depleting activated B cells, T cells, NK cells, monocytes, macrophages, granulocytes, platelets, plasma cells and erythrocytes from peripheral blood. Isolated B cells are bead- and antibody-free and are suitable for any downstream application.

This kit contains anti-CD43 antibodies. CD43 is expressed on B cells from various lymphomas (e.g. CLL, ALL) as well as subsets of activated B cells.

1.2 Principle of Isolation

Add a mixture of monoclonal antibodies against the non-B cells to the starting sample. Add Dynabeads to bind to the non-B cells during a short incubation. Separate the bead-bound cells by a magnet. Discard the bead-bound cells and use the remaining, untouched human B cells for any application (fig.1).

1.3 Materials Supplied

 Dynal[®] B Cell Negative Isolation Kit for untouched human cells contains Depletion Dynabeads and Antibody Mix.

- **2 x 5 ml Depletion Dynabeads**
Supplied in phosphate buffered saline (PBS), pH 7.4, containing 0.1% bovine serum albumin (BSA) and 0.02% sodium azide.
- **1 ml Antibody Mix**
The antibody mix contains mouse IgG antibodies against CD2, CD14, CD16 (specific for CD16a and CD16b), CD36, CD43 and CD235a (Glycophorin A).

 Supplied in PBS with 0.1% BSA and 0.02% sodium azide (NaN₃).

- **The kit will process up to 5 x 10⁸ cells.**

1.4 Additional Materials Required

- Magnet: See www.invitrogen.com/magnets-selection for magnet recommendations.
- Mixer allowing both tilting and rotation.
- Isolation buffer: Ca²⁺ and Mg²⁺ free phosphate buffered saline (PBS) (e.g. Gibco cat. no. 14190-094) supplemented with 0.1% BSA and 2mM EDTA.
- Lymphoprep[™] for MNC preparation




Fig. 1 Simple method for isolating untouched human B cells.

Human Immunodeficiency Virus (HIV), Hepatitis-B Virus (HBV) and Hepatitis-C Virus (HCV).

Storage/Stability

This kit is stable until the expiry date stated on the label when stored unopened at 2-8°C.

Store opened vials at 2-8°C and avoid bacterial contamination.

Keep Dynabeads in liquid suspension during storage and all handling steps, as drying will result in reduced performance. Resuspend well before use.

Technical Support

Please contact Invitrogen DYNAL for further technical information (see contact details).

Warning And Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use.

Follow appropriate laboratory guidelines.

This product contains 0.02% sodium azide as a preservative, which is cytotoxic. **Avoid pipetting by mouth!** Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. When disposing through plumbing drains, flush with large volumes of water to prevent azide build up.

Certificate of Analysis/Compliance is available upon request.

Material Safety Data Sheet (MSDS) is available at <http://www.invitrogen.com>.

Patents and Trademarks

DYNAL®, DYNABEADS® and DYNAL MPC™ are either registered trademarks or trademarks of Invitrogen DYNAL AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

Lymphoprep® is a registered trademark of Axis-Shield Poc AS, Norway.

Intellectual Property Disclaimer

Invitrogen DYNAL will not be responsible for violations or patent infringements that may occur with the use of our products.

Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or

otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: outlicensing@invitrogen.com

Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen DYNAL's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen DYNAL's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen DYNAL, transportation prepaid, or at

Invitrogen DYNAL's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen DYNAL, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen DYNAL's maximum liability is limited in all events to the price of the products sold by Invitrogen DYNAL. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

4. REFERENCES

1. Kai S *et al* (2004) Indoleamine 2,3-dioxygenase is necessary for cytolytic activity of natural killer cells. *Scand. J. Immunol.* 59:177-182.
2. Pillero LM *et al* (2004) T cell homeostasis with thymic hypoplasia due to chromosome 22q11.2 deletion syndrome. *Blood*: 103:1020-1025.
3. Mulder A *et al* (2003) Identification, isolation and culture of HLA-A2-specific B lymphocytes using MHC Class I tetramers. *J. Immunol.* 171:6599-6603.
4. Cao T *et al* (2001) Hepatitis B Virus Core Antigen Binds and Activates Naive Human B Cells in Vivo: Studies with a Human PBL-NOD/SCID Mouse Model. *J. Virol.* 75 (14): 6359-6366.
5. Cuni S *et al* (2004) A sustained activation of PI3K/NF-κB pathway is critical for the survival of chronic lymphocytic leukemia B cells. *Leukemia* 18(8): 1391-1400.

Invitrogen DYNAL is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2007 Invitrogen DYNAL AS, Oslo, Norway.

All rights reserved.

Revised: 09.2007

Printed: 09.2007