

# Dynal® Mouse B Cell Negative Isolation Kit

## For research use only.

### 1. PRODUCT DESCRIPTION

- 1.1 Intended Use
- 1.2 Principle of Isolation
- 1.3 Description of Materials

### 2. PROTOCOLS

- 2.1 Dynabeads Washing Procedure
- 2.2 Sample Preparation
- 2.3 Critical Steps for Cell Isolation
- 2.4 Isolation of Mouse B Cells from Spleen or Lymph Node Leucocytes
- 2.5 Downstream Applications

### 3. GENERAL INFORMATION

#### 1. PRODUCT DESCRIPTION

##### 1.1 Intended Use

Isolate untouched mouse B cells by depleting non-B cells (T cells, monocytes/macrophages, NK cells, dendritic cells, platelets, plasma cells, erythrocytes and granulocytes) from mouse spleen




Fig. 1 Simple method for isolating untouched mouse B cells.

or lymph node cells. Other sources of mouse B cells can also be used as starting material after optimisation for the particular application. Isolated mouse B cells are bead- and antibody-free and are suitable for downstream applications.

##### 1.2 Principle of Isolation

Add a mixture of monoclonal antibodies against the non-B cells to the starting sample. Add Mouse Depletion Dynabeads to bind to the non-B cells during a short incubation. Separate the bead-bound cells by a magnet. Discard the bead-bound cells and use the remaining, untouched mouse B cells for any application (fig. 1).

##### 1.3 Description of Materials

Dynal® Mouse B Cell Negative Isolation Kit contains uniform, superparamagnetic polystyrene Dynabeads® (4.5 µm diameter) coated with a polyclonal sheep anti-rat IgG antibody.

#### Materials Supplied

Dynal Mouse B Cell Negative Isolation Kit contains Mouse Depletion Dynabeads and Antibody Mix.

##### • 2 x 10 ml Mouse Depletion Dynabeads.

Supplied at  $4 \times 10^8$  beads per ml in phosphate buffered saline (PBS), pH 7.4, containing 0.1% bovine serum albumin (BSA) and 0.02% sodium azide.

##### • 2 ml Antibody Mix.

The Antibody Mix contains a mixture of rat monoclonal antibodies against mouse CD43 (Ly-48), CD4 and Ter-119. CD43 is expressed on IL7-responsive pro-B cells, plasma cells and CD5<sup>+</sup> B cells (B-1 cells). These cells of the B lineage are thus removed during the isolation procedure, leaving resting, conventional B cells.

##### • The kit will process up to $1 \times 10^9$ leucocytes.

#### Additional Materials Required

- Magnet (Dynal MPC™): See [www.invitrogen.com/magnets-selection](http://www.invitrogen.com/magnets-selection) for magnet recommendations.
- Heat inactivated Fetal Calf Serum (FCS).
- Mixer allowing both tilting and rotation.

- Buffer 1: PBS w/0.1% BSA and 2 mM EDTA, pH 7.4.

- Buffer 2: RPMI-1640 w/10% FCS.

Keep the buffers cold!

#### Important Notes:

BSA can be replaced by FCS.

EDTA can be replaced by sodium citrate. PBS containing Ca<sup>2+</sup> or Mg<sup>2+</sup> is not recommended.

### 2. PROTOCOLS

#### 2.1 Dynabeads Washing Procedure

Dynabeads should be washed before use.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer 1, or at least 1 ml, and mix.
4. Place the tube in a magnet for 1 min and discard the supernatant.
5. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads (step 2).

#### 2.2 Sample Preparation

Please visit [www.invitrogen.com/cell-isolation](http://www.invitrogen.com/cell-isolation) and follow our QuickLinks for recommended sample preparation procedures.

#### 2.3 Critical Steps for Cell Isolation

- Use a mixer that provides tilting and rotation of the tubes to ensure Dynabeads do not settle at the bottom of the tube.
- Never use less than 200 µl Dynabeads per  $1 \times 10^7$  leucocytes.
- It is critical to follow the magnet recommendations to ensure a successful isolation.

#### 2.4 Isolation of Mouse B Cells from Spleen or Lymph Node Leucocytes

This protocol is based on  $1 \times 10^7$  leucocytes. It is scalable from  $1 \times 10^7$ - $1 \times 10^9$  cells, (see table 1).

1. Transfer 100 µl ( $1 \times 10^7$ ) leucocytes in Buffer 1 to a tube.
2. Add 20 µl heat inactivated FCS.
3. Add 20 µl of Antibody Mix.
4. Mix well and incubate for 20 min at 2-8°C.
5. Wash the cells by adding 2 ml Buffer 1. Mix well by tilting the tube several times and centrifuge at 300 x g for 8 min at 2-8°C. Discard the supernatant.
6. Resuspend the cells in 800 µl Buffer 1.

7. Add 200 µl pre-washed Mouse Depletion Dynabeads.

8. Incubate for 15 min at 18-25°C with gentle tilting and rotation.

9. Resuspend the bead-bound cells by **gently** pipetting 5 times using a pipette with a narrow tip opening, (e.g. a 1000 µl pipette tip or a 5 ml serological pipette).

10. Add 1 ml Buffer 1.

11. Place the tube in the magnet for 2 min.

12. Transfer the supernatant to a new tube.

The supernatant contains the negatively isolated mouse B cells.

Table 1. Volume requirements for mouse

	Working volume per $1 \times 10^7$ leucocytes
Cell volume (step 1)	100 µl
FCS (step 2)	20 µl
Antibody Mix (step 3)	20 µl
Washing (step 5)	2 ml
Resuspension (step 6)	800 µl
Mouse Depletion Dynabeads (step 7)	200 µl
Volume added before magnet (step 10)	1 ml
Dynal MPC recommended	MPC-L/MPC-15/MPC-50

B cell isolation per  $1 \times 10^7$  starting leucocytes.

When working with higher cell numbers, scale up all reagents and volumes accordingly. Up to  $5 \times 10^7$  leucocytes can be processed in a single 15 ml tube. Up to  $2 \times 10^8$  leucocytes can be processed in a single 50 ml tube.

#### 2.5 Downstream Applications

Isolated mouse B cells can be used in applications such as flow cytometry, functional assays and studies on B cell activation, proliferation and differentiation.

### 3. GENERAL INFORMATION

Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

#### Storage/Stability

This kit is stable until the expiry date stated on the label when stored unopened at 2-8°C.

Store opened vials at 2-8°C and avoid bacterial contamination.

Keep Dynabeads in liquid suspension during storage and all handling steps, as drying will result in reduced performance. Resuspend well before use.

## Technical Support

Please contact Invitrogen Dynal for further technical information (see contact details).

## Warning And Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

Follow appropriate laboratory guidelines.

This product contains 0.02% sodium azide as a preservative, which is cytotoxic. **Avoid pipetting by mouth!** Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. When disposing through plumbing drains, flush with large volumes of water to prevent azide build up.

Certificate of Analysis (CoA) is available upon request.

Material Safety Data Sheet (MSDS) is available at <http://www.invitrogen.com>.

## Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

## Patents and Trademarks

Dynal®, Dynabeads® and Dynal MPC™ are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

MyPure™ is a trademark of Invitrogen Dynal AS, Oslo, Norway.

## Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal, transportation prepaid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRES-

SED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

## Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of

this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,  
Invitrogen Corporation,  
1600 Faraday Avenue, Carlsbad,  
California 92008.  
Phone (760) 603-7200.  
Fax (760) 602-6500.  
Email: [outlicensing@invitrogen.com](mailto:outlicensing@invitrogen.com)

---

---

**Invitrogen Dynal is a part of the Invitrogen Group.**

**Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>**

© Copyright 2007 Invitrogen Dynal AS, Oslo, Norway.

All rights reserved.

Revised: 02.2007

Printed: 02.2007