

Dynabeads® Mouse CD3/CD28 T Cell Expander

For research use only

1. KIT CONTENTS
2. PRODUCT DESCRIPTION
3. PROTOCOL
4. GENERAL INFORMATION
5. REFERENCES

1. KIT CONTENTS

Dynabeads Mouse CD3/CD28 T Cell Expander	
Catalog no.	114.52D/114.53D
Dynabeads Mouse CD3/CD28 T Cell Expander*	1x2 ml (114.52D)
	5x2 ml (114.53D)

* 4 x 10⁷ beads/ml in phosphate buffered saline (PBS), pH 7.4, containing 0.1% bovine serum albumin (BSA).

2. PRODUCT DESCRIPTION

Intended Use

This product is intended for the activation and expansion of mouse T cells, e.g. polyclonal T cells, antigen specific T cells, CD4⁺ T cells, CD8⁺ T cells and regulatory T cells (Treg). The expanded cells can be further utilized either for subsequent experimental *in vitro* manipulation, or for adoptive *in vivo* transfer in mouse models of human disease. Applications may include those for the study of infectious disease, transplantation, auto-immunity and cancer.

Principle of T cell Activation and Expansion

Dynabeads Mouse CD3/CD28 T Cell Expander offers a simple method for activation and expansion of T cells that does not require antigen-presenting cells (APCs) or antigen. Simply add Dynabeads Mouse CD3/CD28 T Cell Expander for activation or Dynabeads Mouse CD3/CD28 T Cell Expander plus recombinant IL-2 (rIL-2) for expansion of T cells. Cell cultures showing signs of exhaustion can be re-stimulated several times by adding fresh Dynabeads Mouse CD3/CD28 T Cell Expander and rIL-2. CD8⁺ T cells remain cytotoxic after repeated re-stimulations. Polyclonal T cells and CD4⁺ T cells have been successfully expanded for at least 4 weeks; CD8⁺ T cell cultures typically resist further expansion after 3 weeks. Treg cells retain FoxP3 expression after 2 weeks expansion.

Description of Materials

Dynabeads Mouse CD3/CD28 T Cell Expander are uniform 4.5 µm, superparamagnetic polymer beads coated with an optimized mixture of monoclonal antibodies against the CD3 and CD28 cell surface molecules of mouse cells. The CD3 antibody is specific for the epsilon chain of mouse CD3, which is considered to be a subunit of the TCR complex. The CD28 antibody is specific for the mouse CD28 co-stimulatory molecule, which

is the receptor for CD80 (B7-1) and CD86 (B7-2). Both antibodies are hamster anti-mouse IgGs coupled to the same bead, mimicking *in vivo* stimulation by APCs. Both the bead size and the covalent antibody coupling technology are critical parameters to allow the simultaneous presentation of optimal stimulatory signals to the T cells in culture, thus allowing their full activation and expansion.

Additional Materials Required

- Dynabeads for isolation of desired T cell subset: see www.invitrogen.com/cellisolation for recommended products
- Buffer: Phosphate buffered saline (e.g. Gibco cat.no 10010-023) w/0.1% bovine serum albumin, pH 7.4 (PBS w/0.1% BSA)
- Magnet (DynaMag™ or Dynal MPC™): See www.invitrogen.com/magnets-selection for magnet recommendations
- Culture medium: RPMI-1640 (e.g. Gibco cat.no 12633-012) with 2mM L-Glutamin, 10% FCS/FBS and 100 U/ml penicillin/streptomycin or equivalent
- Heat inactivated Foetal Calf Serum (FCS)
- rIL-2
- Flat bottom tissue culture plates or tissue culture flasks of suitable size
- Humidified incubator

▲ Critical notes

- Follow the recommended volumes
- Avoid air bubbles during pipetting
- Do not use the Dynal MPC-1 magnet for any cell isolations

3. PROTOCOL

The protocol describes the activation and expansion of polyclonal T cells, CD4⁺ or CD8⁺ T cells isolated from spleen or lymph nodes. In general, it is recommended to use a bead to cell ratio of 1:1, and the cell concentrations given in Table 1. Other bead to cell ratios and cell concentrations can be used following optimization for your particular application. We recommend less intense stimulation for expansion of antigen specific T cells (for example bead to cell ratios 1:3, 1:5 and 1:10 could be tried, choosing an optimal ratio for your application). For Treg expansion, see separate protocol below.

Dynabeads Washing Procedure

Dynabeads should be washed before use.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer, or at least 1 ml, and mix.
4. Place the tube in a magnet for 1 min and discard the supernatant.

5. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer, as the initial volume of Dynabeads (step 2).

Isolation procedure

See www.invitrogen.com/cellisolation for recommended cell isolation products (positive or negative isolation), and refer to the procedure described in the respective package insert. It is critical for isolation of Treg cells (flow sorting or magnetic bead isolation) to use an anti-CD25 antibody that do not block the binding of IL-2.

Note: Before flow cytometric analysis or use in downstream applications, Dynabeads and bead-bound cells should be removed. During the first days of expansion, some cells will bind strongly to the beads. Resuspend the bead/cell suspension thoroughly by pipetting to increase cell recovery, separate on a magnet (either in a culture plate or following transfer to suitable tubes) and collect supernatant containing the T cells. The bead-bound cell fraction can be cultured overnight and the above process repeated to further increase T cell recovery.

Short term stimulation/activation of mouse T cells

1. Start with 1-1.5 x 10⁶ purified T cells/ml culture medium in a suitable tissue culture plate or tissue culture flask.
2. Add Dynabeads Mouse CD3/CD28 T Cell Expander at a bead-to-cell ratio of 1:1.
3. Incubate in a humidified CO₂ incubator at 37°C.
4. Remove the beads at any time before analysis and further use of the cells.

Activation and expansion of polyclonal Mouse T cells

1. Start with 1-1.5 x 10⁶ purified T cells/ml culture medium in a suitable tissue culture plate or tissue culture flask.
2. Add Dynabeads Mouse CD3/CD28 T Cell Expander at a bead-to-cell ratio of 1:1.
3. Add 30 U/ml rIL-2

4. Incubate in a humidified CO₂ incubator at 37°C.
5. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures.
6. Count the cells at least twice weekly after thorough re-suspension.
7. When the cell density exceeds 2.5 x 10⁶ cells/ml or when the medium becomes yellow, split cultures back to a density of 0.5-1 x 10⁶ cells/ml in culture medium containing 30 U/ml rIL-2.

Re-stimulation

Re-stimulation is typically necessary when cell shrinking and reduced rate of proliferation are observed. Guidelines for re-stimulation are provided in Table 2, although we recommend optimization for your particular application.

Do not use an excess volume of Dynabeads Mouse CD3/CD28 T Cell Expander, as excess Dynabeads per cell may inhibit expansion.

Before re-stimulation, remove the used Dynabeads by transferring the cells to a suitable tube. Place the tube in a magnet for 1-2 minutes until the Dynabeads are separated. Transfer the supernatant containing the cells to a new tube. Continue as described in step 1 below.

1. Count the cells and resuspend to a density of 1 x 10⁶ cells/ml in culture medium with 30 U/ml rIL-2 in a suitable culture plate or tissue culture flask.
2. Add Dynabeads Mouse CD3/CD28 T Cell Expander at a bead-to-cell ratio of 1:1.
3. Add 30 U/ml rIL-2.
4. Incubate in a humidified CO₂ incubator at 37°C.
5. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures.
6. Count the cells at least twice weekly after thorough re-suspension.
7. When the cell density exceeds 2.5 x 10⁶ cells/ml or when the medium becomes yellow, split cultures back to a density of 0.5-1 x 10⁶ cells/ml in culture medium with 30 U/ml rIL-2.

Table 1: Volume recommendations for bead-to-cell ratio = 1:1

	1 x 10 ⁵ T cells	1 x 10 ⁶ T cells	35 x 10 ⁶ T cells*
Type of culture plate/flask	Per well in 96-well plate	Per well in 24-well plate	175 cm ² tissue culture flask
Dynabeads Mouse CD3/CD28 T Cell Expander	2.5 µl	25 µl	875 µl
rIL-2	30 U/ml	30 U/ml	30 U/ml
Seeding volume (medium)	100-200 µl	1-2 ml	35-70ml

* Average number of T cells obtained from one mouse spleen.

Table 2: Re-stimulation guidelines for anti-CD3/CD28-expanded cultures

Cell type	First re-stimulation**	Subsequent re-stimulations**
CD4 ⁺ (polyclonal)	8-10 days	8-11 day intervals
CD8 ⁺ (polyclonal)	7-9 days	7-10 day intervals
T cells	7-9 days	10-12 day intervals

** Establish optimal times for your particular cells. Please note that these are only generic guidelines.

Activation and Expansion of Regulatory T cells

1. Start with 1-1.5 x 10⁶ cells/ml culture medium in a suitable tissue culture plate (96 well plate: 10⁵ total cells/well, 24 well plate: 1-1.5 x 10⁶ total cells/well).
2. Add Dynabeads Mouse CD3/CD28 T Cell Expander at a bead-to-cell ratio of 2:1.
3. Add 2000 U/ml rIL-2.
4. Incubate in a humidified CO₂ incubator at 37°C.
5. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures.
6. Count the cells at least twice weekly after thorough re-suspension.
7. When the cell density exceeds 2.5 x 10⁶ cells/ml or when the medium becomes yellow, split cultures back to a density of 0.5-1 x 10⁶ cells/ml in culture medium containing 2000 U/ml rIL-2.

4. GENERAL INFORMATION

Storage and Stability

This kit is stable until the expiry date stated on the label when stored unopened at 2-8°C.

Store opened vials at 2-8°C and avoid bacterial contamination.

Keep Dynabeads in liquid suspension during storage and all handling steps, as drying will result in reduced performance. Resuspend well before use.

Technical Support

Please contact Invitrogen Dynal for further technical information (see contact details).

Warnings and Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

Certificate of Analysis/Compliance is available upon request.

Patents and Trademarks

Dynal®, Dynabeads®, DynaMag™ and Dynal MPC™ are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal,

transportation prepaid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

Limited Use Label Licence

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund.

No. 305: Invitrogen T Cell Expansion Technology – The use of this product for *ex vivo* activation or expansion of human T-cells is

covered under various patents owned and licensed from the U.S. Navy, the University of Michigan and the Dana Farber Cancer Institute and is sold for research use only. This product is not to be used for activation or expansion of T cells modified through gene transfer to specifically modify the T cells to produce secreted or cell-surface membrane-bound proteins not normally expressed in significant levels by such T cells, unless the proteins directly enable the selection, or directly modify or preserve the function, of the T cells. Modifications that are acceptable would include, for example, introducing genes that direct expression of cell-surface bound proteins that bind the T cell to specific target cells.

For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,
Invitrogen Corporation,
1600 Faraday Avenue, Carlsbad,
California 92008.
Phone (760) 603-7200.
Fax (760) 602-6500.
Email: outlicensing@invitrogen.com

5. REFERENCES

1. Ito F *et al.* (2003) Antitumor reactivity of anti-CD3/anti-CD28bead-activated lymphoid cells: implications for cell therapy in a murine model. *J Immunotherapy* 26(3): 222-233.
2. Li J *et al.* (2007) Transfer of in vitro expanded T lymphocytes after activation with dendritomas prolonged survival of mice challenged with EL4 tumor cells. *Int J Oncol* 31(1): 193-197.
3. Masteller EL *et al.* (2005) Expansion of functional endogenous antigen-specific CD4⁺CD25⁺ regulatory T cells from nonobese diabetic mice. *J Immunol* 175: 3053-3059.

Invitrogen Dynal is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2008 Invitrogen Dynal AS, Oslo, Norway. All rights reserved.

SPEC-06299