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
→ 2013 JASTRO

- A. **Do not cover the Calibration Window** - 1.5 cm wide area where scanning starts;
- B. Place film at the center of the scan window perpendicular to the scan direction;
- C. Choose landscape (preferred), or portrait scanning orientation, but **be consistent**;

LANDSCAPE ORIENTATION: SHORT EDGE OF FILM PARALLEL TO SCAN DIRECTION

Figure 6: Film Placement for Flatbed Scanning

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→ 2013 JASTRO

so that the adjustment icons appear gray;  
F. For optimum results use FilmQPro software and the triple-channel dosimetry protocol.

1  
2  
3  
4  
5  
NO YES  
TURN OFF ALL IMAGE ADJUSTMENT FEATURES  
6

7

Figure 6a: Scan Driver Set-up for Epson Scanners


Figure 6b: Scan Driver Set-up for Epson Scanners

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→ 2013 JASTRO

2103 Niroomand-Rad et al. Radiochromic film dosimetry

marked with at least three fiducial marks and then exposed to a uniform dose of  $D_1$  (normally about 15–20 Gy for MD-55-1 film or about 10 Gy for MD-55-2 film). The mean optical density of the film,  $OD_1(i, j)$ , in the area of interest can be determined after 24 h. Then the film is exposed to an unknown dose of  $D_2$  in the desired irradiation field. The film is scanned, again after 24 h, and the two images are correlated using the three fiducial marks. The optical density from the first exposure,  $OD_1(i, j)$ , and the second exposure,  $OD_2(i, j)$ , can be used to obtain a nonuniformity corrected net optical density,  $OD_{net}(i, j)$ , using the following relations:

$$OD_{net}(i, j) = \frac{OD_2(i, j) - OD_1(i, j)}{f(i, j)}, \quad (11)$$

where

$$f(i, j) = \frac{OD_1(i, j) - OD_b(i, j)}{OD_1(i, j) - OD_b(i, j)}, \quad (12)$$

and  $OD_b(i, j)$ , which represents the background (i.e., fog) readings, is the film optical density prior to the first exposure. The double-exposure technique may require writing a computer program for image handling and data analysis. Note that the double-exposure technique has some limitations. For instance, the accuracy of matching the fiducial marks in the two exposures is limited and should be evaluated by the user. Therefore, improvement of the film uniformity by the manufacturer is necessary in order to explore further the full application of radiochromic film.

**M. Film calibration and sensitivity**

Radiochromic film should be calibrated using a large well-characterized uniform radiation field. Film should be placed on the central portion of a large photon beam (such as a 40 × 40 cm) at depth of interest (preferably >5 cm). The characteristics of the calibration beam should be determined by some other dosimeter (such as an ionization chamber). This would allow direct film calibration in terms of absolute dose within the dose range of interest.

The relationship between absorbed dose and film response should be determined. This relationship can be plotted as a curve, often known as a calibration curve. The slope of the calibration curve decreases as dose increases [Fig. 12(A)]. The calibration curve can provide information for conversion of film response to dose and vice versa.

The relationship between dose and film response can also be tabulated. The change in film response per unit absorbed dose can be represented by a single number for a net optical density up to 1.0. This number, defined here as film average sensitivity, is the average change in response (i.e., net OD) per unit absorbed dose calculated over the lower, more linear portion of the calibration curve. This number depends on one or more of the following: (1) the wavelength used for readout, (2) the particular dosimeter used for readout, (3) film batch, (4) the delay between irradiation and readout, (5) beam quality of the calibration source, and (6) other factors (such as temperature and humidity) previously discussed. Table II is an example of a film sensitivity calculation for

MD-55-2 (lot No. 941206), exposed to uniform  $^{60}\text{Co}$  doses and scanned after a week delay with a 633 nm LKB dosimeter at NIST. The average sensitivity (net OD/Gy) for this combination, for the dose range 0–20 Gy, from the data above the line is  $0.023 \pm 0.0007$  (3.2%). Figure 12 is a plot of (A) net OD and (B) net OD per unit absorbed dose measured two days postirradiation. Figure 12(A) has two fairly linear portions, one from 0 to 30 Gy, and the other from 30 to 100 Gy. The average sensitivity for the dose range 30–100 Gy drops about 15% as compared to the one for the dose range 0–30 Gy. This average sensitivity can be used for conversion of film response to dose for these conditions and is clearly inappropriate for higher doses or other conditions. Each user should determine film response and sensitivity for their specific conditions.

**N. Dose calculation using a double-exposure technique**

Klassen *et al.*<sup>50</sup> used a double-exposure technique to evaluate MD-55-2 as a high precision dosimeter using a Gafar

Unknown doses, in the range 3–7 Gy, were given to a set of films and the increase in OD was measured. Then a known dose of 6 Gy was given to each film and the increase in OD was measured. In order to measure the change in sensitivity of the films between the first and second exposures, control

2103

Fig. 12. Comparison of relative optical density for MD-55-2 film exposed to uniform radiation fields with  $^{60}\text{Co}$  gamma ray beam in (a) longitudinal and (b) transverse directions, respectively (Ref. 70).

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
### Response Equivalence Scanner-to-Scanner

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


### Measuring - Film Flatness


Callier effect - Light change from collimated to diffused  
**Avoidable !**

Ashland, Inc.  
A. Micke, D. Lewis, Oregon, November 2012  
www.FilmQAPro.com

### Measuring - Film Flatness

**Rescan solves the problem**  
Multi-channel dosimetry recognizes such disturbances!

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
→ 2013 JASTRO

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### Newton's Rings Pattern

- Interference pattern formed by an air gap of varying thickness between two closely spaced surfaces
  - Constructive/destructive interference between reflected/refracted rays
  - Monochromatic light yields light and dark bands
  - White light yields colored bands
    - Constructive interference for  $\lambda_1$ , destructive interference for  $\lambda_2$

### Avoid Contact Between Film and Glass

- Make a simple frame
  - Approximately 1-2 mm thick
  - Cardboard, plastic, rubber

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