# **GAFCHROMIC 事始め**

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
FILM-SCAN-NEW Ver 3.0 for GAFCHROMIC and X-Ray Film and GAFCHROMIC Film

Max/Min :28702.00 / 2239.00

Display gamma: 0.30

Scan info Scanner Select

**Step 1: Select Scanner**

Kind of Format  
 OLD(Normal)  New(2Ch)  
 CH3(R,G,B)

R G B Select  
 R  G  B

With Median Filter 3  
 With Self Calibration  
 With Base Subtraction Method  
 Base data clip

Number of Scan 1

**Step 2 : Scan Start**

Save data

Kind of Display Format  
 Histogram Display  
 Profile Display

X Profile  
Y Profile

Read Test


Close

C:\Users#R\_TECH\FILM\_SCAN\_NEW#Temp#


The screenshot shows the main interface of the FILM-SCAN-NEW software. At the top, it displays the title and version information. Below this, there are two main panels. The upper panel shows a scan of a film strip with a dashed box indicating the scanning area. The lower panel is a histogram showing the distribution of scan data. The histogram has a red line representing the scan data and a blue line representing the profile data. The x-axis ranges from 57.30 to 544.00, and the y-axis ranges from 3654.93 to 36549.30. On the right side, there is a control panel with various settings and buttons. The control panel includes a 'Step 1: Select Scanner' section with radio buttons for 'Kind of Format' and 'R G B Select', and checkboxes for 'With Median Filter', 'With Self Calibration', 'With Base Subtraction Method', and 'Base data clip'. Below this is a 'Step 2 : Scan Start' section with a 'Save data' button. At the bottom of the control panel, there are 'X Profile' and 'Y Profile' sliders, a 'Read Test' button, and a 'Close' button. The status bar at the bottom shows the file path: C:\Users#R\_TECH\FILM\_SCAN\_NEW#Temp#.

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