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## DEPLETION OF SSEA-4<sup>+</sup> HUMAN EMBRYONAL STEM CELLS USING DYNABEADS

Pauline Lieu,<sup>1</sup> Ying Liu,<sup>1</sup> Steve Spusta,<sup>2</sup> Anette Kullmann,<sup>3</sup> Xianmin Zeng,<sup>2</sup> Berit Reed,<sup>3</sup> John Chesnut,<sup>1</sup> Karoline W. Schjetne,<sup>3</sup>  
<sup>1</sup> Invitrogen Corporate, Life Technologies, Carlsbad, California, US  
<sup>2</sup> Buck Institute of Age Research, Novato, California, US  
<sup>3</sup> Invitrogen Dynal, Life Technologies, Oslo, Norway


AIM: Develop a protocol for gentle depletion of undifferentiated human embryonal stem cells (hESCs)

### Background

Human embryonal stem cells (hESC) research has become one of the most exciting and fastest growing areas in cell biology today and has an enormous potential for clinical applications and in drug discovery. Since embryonic stem cells have unlimited proliferative capacity and can differentiate into multiple cell types, they provide an inexhaustible source of precursors and differentiated cells for many therapeutic applications. Today, scientists have successfully used various methods to differentiate human embryonic stem cells into many different cell types, such as neurons and cardiomyocytes. One of the challenges is to obtain a pure and homogenous population of target cells for drug screening or clinical applications.

### Materials and Methods

Human embryonal stem cells (hESCs) were differentiated for 16 days into neuronal precursor cells using NSC medium. Cells were detached from the culture flask wall and analyzed for SSEA-4 expression. After differentiation, some cells still expressed SSEA-4 indicating that undifferentiated hESC are present (Fig 1a). Dynabeads® SSEA-4 were added to the cell suspension and placed on a magnet after a short incubation. The bead-free cells were transferred to a new tube and analyzed by flow cytometry showing >99% depletion of undifferentiated SSEA-4<sup>+</sup> cells (Fig 1b). Also, stainings of Oct4 (Fig 2) and nestin (Fig 3) were performed on the differentiated NSC demonstrating that the depletion protocol efficiently removed Oct4<sup>+</sup> cells (ESC) leaving a pure nestin positive (NSC) cell population. To show that the differentiated cells were unaffected by the depletion procedure, the cells were further differentiated into motor neurons and glial precursors. The morphology and specific expression of beta tubulin and A2B5 markers of the differentiated cells are visualized in fig 4.

### Dynabeads® SSEA-4 remove hESC from the differentiated NSC populations


Fig 1: Human ESC stem cells (BC01) differentiated into NSC were stained with anti-SSEA-4 antibody pre- and post depletion using Dynabeads® SSEA-4 prior to analysis by flow cytometer.

### Oct4 positive cells are removed after depletion of SSEA-4 depletion


Fig 2: Oct4 positive cells (green) were detected in the differentiated NSC population before depletion (left), demonstrating that undifferentiated ESC were still present. After depletion of SSEA-4<sup>+</sup> cells using only background staining is observed (right).

### Neural stem cells express nestin after SSEA-4 depletion


Fig 3: Differentiated NSC were stained with anti-nestin antibodies (pink) after depletion of SSEA-4 positive cells.

### Neural stem cells depleted of SSEA-4<sup>+</sup> cells differentiate into motor neurons and glial precursor cells


Fig 4: NSC depleted of SSEA-4<sup>+</sup> cells where further differentiated into motor neurons and glial precursors prior to staining with beta actin to identify neural precursor and anti-A2B5 mab to identify oligodendrocytes.

Ordering information  
 Dynabeads® SSEA-4  
 DynaMag™-2  
 DynaMag™-15  
 Anti-SSEA-4 antibody  
 Anti-A2B5 antibody  
 StemPro® hESC SFM

Cat.no  
 111.60D  
 123.21D  
 123.01D  
 41-4000  
 41-1000  
 A1000701  
 Application  
 Depletion of undifferentiated SSEA-4<sup>+</sup> hESCs  
 Holds 2 ml tubes  
 Holds 5 ml and 15 ml tubes  
 Detection of SSEA-4<sup>+</sup> cells  
 Detection of A2B5<sup>+</sup> cells  
 Culture of hESCs

### CONCLUSIONS

- Dynabeads® SSEA-4 efficiently deplete undifferentiated SSEA-4<sup>+</sup> hESCs
- SSEA-4 depleted neural stem cells further differentiate into neural precursor cells such as motor neurons and glial precursor cells
- Isolation is tube-based, fast and gentle to your cells

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