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##### Prepare IgG- Dynabeads complex

1. Take 100  $\mu$ l Dynabeads Protein A (or G, depending on species IgG isotypes) to 10  $\mu$ g IgG.
2. Place the tube on the magnet and discard the supernatant. Wash Dynabeads with PBS-Tween-20 (0.01% Tween-20 can be added to the buffer to prevent aggregation of the Dynabeads Protein A - Ig complex.)
3. Resuspend Dynabeads Protein A with 100  $\mu$ l PBS-Tween20
4. Add antibody (10  $\mu$ g IgG) to Dynabeads Protein A
5. Rotate at 4C more than 1hr
6. Place the tube on the magnet and Remove unbound solution and Wash with PBS-Tween20 twice.

##### Crosslinking by DMP

1. Place the tube on the magnet and remove sup.
2. Add 1 ml 0.2M Na-borate (Borax)/boric acid pH9.0 to the Dynabeads - IgG complex. Wash twice with 0.2M Na-borate (Borax)/boric acid pH9.0.

##### Prepare A and B

A: 0.2M boric acid solution (12.4g in 1000mL)

B: 0.2M sodium borate solution (19.05 g in 1000mL) = 0.05M sodium tetraborate (Borax)

To make 0.2M sodium borate (pH9.0) buffer, Mix 50mL of A and 59mL of B

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3. Resuspend the Dynabeads - IgG complex in 1 ml of 0.2M sodium borate (pH9.0) containing 20 mM DMP (dimethyl pimelimidate dihydro-chloride, Sigma) .  
5.4mg DMP powder is dissolved in 1 ml of 0.2M sodium borate (pH9.0) buffer.  
This solution must be prepared immediately before adding to the Dynabeads - IgG complex.

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4. Incubate with rotator for 30 minutes at 20° C.
5. Place the tube on the magnet and discard the supernatant.
6. Remove the tube from the magnet and stop the reaction by resuspending the Dynabeads - Ig complex in 1 ml of 50 mM Tris pH 7.5/PBS, and incubate for 15 minutes with rotator. This is for quenching of unreacted DMP.

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7. Wash the beads twice with 1 ml of 100mM Glycine (pH 2.5) + 150mM NaCl to remove uncoupled IgG on the beads. Soon after washing, neutralize the beads by adding 1 ml of 0.1M Tris-Cl(pH8.0).
8. Place the tube on the magnet and discard the supernatant.
9. Wash the crosslinked Dynabeads - IgG complex 2 times with 1 ml PBS-Tween20.
10. Resuspend the Dynabeads - IgG complex in 100 μl PBS-Tween20. add 0.1% sodium azide. for store at 4C.

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