


## PRODUCT DESCRIPTION

Gentle Cell Dissociation Reagent is an enzyme-free reagent suitable for the dissociation of human embryonic stem cells or human induced pluripotent stem cells into cell aggregates for routine passaging or into a single cell suspension.

## COMPONENTS

Gentle Cell Dissociation Reagent does not contain enzymes or other proteins. The product is chemically defined and is sterility tested.



## STABILITY AND STORAGE

Product is stable at 15 - 25°C until expiry date as indicated on label.

## DIRECTIONS FOR USE

For complete instructions on how to use Gentle Cell Dissociation Reagent with TeSR™-E8™, mTeSR™1, TeSR™2, the Vitronectin XF™ Kit, the StemAdhere™ Defined Matrix for hPSC Kit, STEMdiff™ Neural Induction Medium, or the STEMdiff™ Definitive Endoderm Kit, please refer to the corresponding documents listed in the table below. Documents are available on our website at [www.stemcell.com](http://www.stemcell.com) or to request a copy, contact STEMCELL Technologies' Technical Support at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

| PRODUCT                                 | CATALOG #               | DOCUMENT # |
|---|-------------------------|------------|
| TeSR™-E8™                               | 05940                   | 29267      |
| mTeSR™1                                 | 05850/05857/05870/05875 | 29106      |
| TeSR™2                                  | 05860/05880             | 29106      |
| Vitronectin XF™ Kit                     | 07190                   | 29285      |
| StemAdhere™ Defined Matrix for hPSC Kit | 07170                   | 28795      |
| STEMdiff™ Neural Induction Medium       | 05831                   | 28782      |
| STEMdiff™ Definitive Endoderm Kit       | 05110                   | 29550      |

Copyright © 2013 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and STEMdiff are trademarks of STEMCELL Technologies Inc. E8, mTeSR and TeSR are trademarks of WARF. Vitronectin XF and StemAdhere are trademarks of Primorigen Biosciences, Inc. All other trademarks are the property of their respective holders.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.  
FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.



TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM • FOR GLOBAL CONTACT DETAILS VISIT [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

VERSION 2.2.0  
DOCUMENT #28798