

Panomics/Affymetrix 社

転写因子関連製品

ディスカウントキャンペーン！！

2010年5月24日 ～ 2010年7月末日まで



期間中 Panomics/Affymetrix 社
の転写因子関連製品は

25%OFF!!

主なキャンペーン対象製品

- EMSA
- TF ELISA
- Luciferase Reporter Vector
- Protein/DNA Array
- Stable Cell Line
- Procarta Transcription Factor Assay



Panomics/Affymetrix 社で取り扱っている転写因子関連試薬は他にも多数ございますので、詳しくはカタログをご請求ください。

また、ご質問等ございましたら

電話 **(03-3593-3385)** 又はメール (**techservice@veritastk.co.jp**)
までお問合せください。

EMSA Kit

製品概要

EMSA (Electrophoretic Mobility Shift Assay) は DNA/TF 複合体を電気泳動上で Shift させてサンプル内に含まれる TF の有無を検証する転写因子研究の基本的なアッセイです。

特長

- ・ 425 EMSA Probe から選択可
- ・ ラジオアイソトープを使用していないため、安全かつ特別な処分が不要
- ・ カスタムデザイン可能
- ・ 1 ウェルあたり 4～8ug のタンパクで検出
- ・ ビオチン化されたプローブ

製品梱包

- ・ 各種反応用溶液
- ・ 各種 Buffer
- ・ コントロールプローブ
- ・ 検出用発光試薬



TF ELISA Kit

製品概要

特定の転写因子のサブユニットを検出できる ELISA キットです。


特長

- ELISA フォーマットによる簡便操作
- プローブと抗体により高い特異性を実現
- 0.5ug の核タンパク質量で測定可能
- 4 時間の実験操作で終了

製品梱包

- Target probe(cold probe も含む)
- 1st, 2nd 抗体
- 各種 Buffer
- 各種反応用溶液
- プレート及びプレートシール

=測定原理=



転写因子リスト

###	転写因子名	180_HNF-3beta	131_PPAR alpha	109_Sp1	104_AP1
510	NFkB Rel-b	170_EGR1	121_PPAR gamma	108_MEF2	103_ER alpha
410	NFkB c-Rel	160_Smad2	112_NFkB p65	107_Stat1	102_HIF
310	NFkB p52	150_Smad4	111_NFkB p50	106_GR	
200	E2F-1	140_FKHR	110_GATA1	105_p53	

TF Reporter Stable Cell Line

製品概要

レポーター Gene が組み込まれた安定された細胞株です。


特長

- Reporter Vector 発現細胞株
- トランスフェクション効率を気にせず実験
- スクリーニングツールに最適

製品梱包

凍結細胞

CHO/GFP-GRα Stable Cell Line



control

dexamethasone

品名	梱包単位
CHO/ GFP-GR α Stable Cell Line	1vial
CHO/ GFP-NF κ Bp65 Stable Cell Line	1vial
TAD Kinase-luciferase HeLa Reporter Stable Cell Line	1vial
HeLa GR-Luciferase Stable Cell Line	1vial
293 GR-Luciferase Stable Cell Line	1vial
HIF-luciferase NIH3T3 Reporter Stable Cell Line	1vial
NFkB-luciferase C2C12 Reporter Stable Cell Line	1vial
NFkB-luciferase NIH3T3 Reporter Stable Cell Line	1vial
NFkB-luciferase 293 Reporter Stable Cell Line	1vial
NFkB-luciferase HeLa Reporter Stable Cell line	1vial
SRF-luciferase HeLa Reporter Stable Cell line	1vial
NF-AT-luciferase K562 Reporter Stable Cell line	1vial
CREB-luciferase CHO Reporter Stable Cell line	1vial
CREB-luciferase 293 Reporter Stable Cell line	1vial
AP1 luciferase 293 Reporter Stable Cell Line	1vial
NF-AT-Luciferase HeLa Reporter Stable Cell line	1vial
Stat1-Luciferase HeLa Reporter Stable Cell Line	1vial
Stat3-Luciferase HeLa Reporter Stable Cell Line	1vial
NFkB-Luciferase A549 Reporter Stable Cell Line	1vial

Protein/DNA Array

製品概要

プロテインと DNA との結合をメンブレンを使って評価する試薬です。

*メンブレンとそのメンブレンにドットされている項目は、TF-TF Interaction Array 試薬に使われるメンブレンと共通しています。

特長


- 複数の転写因子活性を同時評価
- 345 転写因子をラインナップ
- 転写因子スクリーニングに活用

Function Specific Protein/DNA

- 一度に最大 345 多項目の転写因子プロファイルが可能
- ラジオアイソトープを使用していないので安全
- ケミルミを使用した高感度測定

製品梱包

- 各種反応用溶液
- 各種 Buffer
- コントロールプローブ
- 検出用発光試薬



Luciferase Vector


製品概要

特定の遺伝子の発現を Luciferase の活性で評価するためのベクターです。

特長

- In vivo における転写因子の結合活性評価
- 50 以上の転写因子のラインナップ
- Stable cell line の構築にも利用可能

転写因子リスト



###_転写因子名	087_VDR	057_p53	033_GR	019_ELK1
128_HIF1	079_Stat4	052_NFkB(2)	032_GATA4	017_E2F(2)
127_Stat1(2)	077_STAT3(1)	051_NFkB(1)	031_GATA3	016_E2F(1)
117_SP1(2)	075_Stat1 p84/p91	050_NF-AT	030_GATA2	007_AR
116_Smad	073_SRF	045_MEF3	029_GATA1 / GATA2	005_AP3
115_Smad3 / Smad4	072_SRE	044_MEF2	028_GATA1	003_AP1(2)
114_RXR(2)	071_Sp1	043_MEF1	027_GATA	002_AP1(1)
094_E2F1(3)	070_RXR	040_ISRE	026_GAS / ISRE	000_TA-luc (Empty Control)
093_CRE(1)	068_RAR	039_IRF1	025_GAS	
090_YY1	067_PR	038_HSE	020_ER	

Procarta Transcription Factor Assay Kit


製品概要

Luminex/Bio-Plex を用いて、複数の転写因子の活性化を同時に測定することが可能です。転写因子が結合したビオチン標識されたプローブのみを separation plate で精製します。Luminex ビーズに結合した anti-sense プローブと hybridization させ、Luminex/Bio-Plex により検出します。微量の核抽出液から、44 種類の転写因子の活性化を測定することが可能です。

特長

- 最大 44 種類同時測定
- 83 種類のラインナップ
- ユニークな方法によりバックグラウンドの減少 (特許取得済)
- 核タンパク質の抽出は専用 Buffer により容易
- パネル以外のカスタムキットも対応可能
- 内部コントロールによる内部補正


= 実験の流れ =



1) 細胞・核抽出液とビオチン標識された 2 本鎖 cis-elements を反応。

2) 転写因子と結合したプローブの精製 (separation plate)

3) 転写因子 / プローブ複合体を変性させ、プローブの回収



4) 2 本鎖 oligos を変性させ、beads 標識 anti-sense sequence と hybridization

5) Luminex により測定

転写因子リスト


《 Panel I 》

RUNX/AML	CREB	GR/PR	NFAT	PAX3
AP1	E2F1	HIF-1	NF-E1(Y1)	PAX5
AP2	ELK1	HNF1	NF-E2	PPAR
AR	ER	IRF	NFkB1	SMAD
ATF2	ETS/PEA	ISRE	NKX-2.5	STAT1
BRN3	FAST1	MEF-2	NF-Y	STAT3
CEBP	FKHR-1	MYOD	OCT	STAT4
C-MYB	GATA-1	NF-1	P53	STAT5
TF II D				

《 Panel II 》

EFI-1	XBP-1	NPAS2	PUR-1	PDX-1
HAS+HBS	SIE	CCAAT	USF-1	SRE
NFkB	GAG	CEF-1	PIT-1	LF-A1
TR	ALF-1/TAL-1	HSF	AP-4	RB
H4TF	c-MYC	IKAROS	ANTI-RE	XRE
LvF	MRE	TREF-1/2	TR (DR-4)	SRY
TFE-3	HiNF	Neuro-D1	HBS/XBP	CDP
GFI-1	MTF	ELF-1	EGR	COUP-TF
E47	AP-1	KPF-1	TF II D	

データ例



PMA 処理を 4 時間行った Hela 細胞の核抽出液を用いて、転写因子の活性化を測定
AP-1 の活性化を検出 (positive control)

<TAD kinase Assay キット>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-LK0001	AP2alpha TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0002	ATF2 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0003	c-Jun TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0004	c-fos TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0005	CREB1 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0006	CREB2 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0007	CHOP TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0008	E2F1 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0009	ELK1 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0010	MEF2A TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0011	MEF2C TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0012	Myc TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0013	OCT-1 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0014	PPARalpha TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0015	RelA TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0016	SAP-1 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0017	STAT1 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0018	STAT2 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0019	STAT3 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0020	STAT4 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0021	STAT5A TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0022	STAT5B TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800
PN-LK0023	STAT6 TAD in vivo Kinase Assay Kit	10 ug	¥121,000	¥90,800

<Vector 関連>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-LR0000	TA-luc (Empty Control) Luciferase Reporter Vector 20ug	20 ug	¥102,000	¥76,500
PN-LR0002	AP1(1) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0003	AP1(2) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0005	AP3 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0007	AR Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0016	E2F(1) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0017	E2F(2) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0019	ELK1 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0020	ER Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0025	GAS Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0026	GAS/ISRE Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0027	GATA Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0028	GATA1 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0029	GATA1/GATA2 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0030	GATA2 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0031	GATA3 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0032	GATA4 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0033	GR Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0038	HSE Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0039	IRF1 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0040	ISRE Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0043	MEF1 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0044	MEF2 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0045	MEF3 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0050	NF-AT Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0051	NFkB(1) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0052	NFkB(2) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0057	p53 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0067	PR Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0068	RAR Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0070	RXR Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0071	Sp1 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0072	SRE Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0073	SRF Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0075	Stat1 p84/p91 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0077	STAT3(1) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0079	Stat4 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0087	VDR Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0090	YY1 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0093	CRE(1) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0094	E2F1(3) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0114	RXR(2) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0115	Smad3/Smad4 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0116	Smad Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0117	SP1(2) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0127	Stat1(2) Luciferase Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR0128	HIF1 Luciferase Reporter Vector	10 ug	¥102,000	¥76,500

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-LR1001	ABL1 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1002	BRCA1 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1003	CGPR (CALCR) Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1004	GAGE1 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1005	IRF7 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1006	POMC Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1007	ATF2 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1008	STAT1 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1009	14-3-3-sigma (SFN) Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1010	CIITA Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1011	IFNgamma Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1012	IGRP Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1013	IL1alpha Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1014	IL2 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1015	TNFalpha Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1016	SURVIVIN Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1017	ERalpha Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1018	RB1 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1019	APC Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1020	CFTR Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1021	G6PD Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1022	hMLH1 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1023	IL4 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1024	JUNB Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1025	MyoD Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1026	VHL Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1027	v-myc Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1028	MDR1 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1029	COX2 (PTGS2) Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1030	HOXA2 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1031	HOXA5 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1032	FHIT Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1033	THBS2 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-LR1034	TIMP3 Gene Promoter Reporter Vector	10 ug	¥102,000	¥76,500
PN-ME0001	PPARalpha CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0002	PPARbeta (PPARdelta) CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0003	PPARGamma-1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0004	PPARGamma-2 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0005	NFkB p50 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0006	NFkB p65 (RelA) CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0007	GR (NR3C1) CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0010	E2F2 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0014	IRF1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0018	MAFK CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0019	MAX CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0020	CITED1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0022	NF-YC CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0024	NR1I2 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0025	NR1I3 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0029	PREB CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0032	SRY CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0033	STAT1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0035	STAT5B CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0040	USF1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0043	ERalpha CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0044	FOS CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0045	JUN CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0046	NF-ATC1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0047	p53 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0049	GATA1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0054	PBX1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0055	ATF3 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0059	ATF1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0060	CDX2 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0061	NROB1 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0062	HNF4gamma CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0064	ATF2 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0065	BTG2 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0066	GF11 CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0071	SMAD1 (MADH1) CMV Expression Vector	15 ug	¥281,000	¥210,800
PN-ME0052	MYC CMV Expression Vector	個	¥281,000	¥210,800
PN-ME0052	MYC CMV Expression Vector	個	¥281,000	¥210,800

<TF ELISA キット>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-EK1020	HIF ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1021	HIF ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1030	ER alpha ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1031	ER alpha ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1040	AP1 ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1041	AP1 ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1050	p53 ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1051	p53 ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1060	GR ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1061	GR ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1070	Stat1 ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1071	Stat1 ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1080	MEF2 ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1081	MEF2 ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1090	Sp1 ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1091	Sp1 ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1100	GATA1 ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1101	GATA1 ELISA Kit with nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1110	NFkB p50 ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK1111	NFkB p50 ELISA Kit with nuclear extraction kit	1 kit	¥137,000	¥102,800
PN-EK1120	NFkB p65 ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK1121	NFkB p65 ELISA Kit with nuclear extraction kit	1 kit	¥137,000	¥102,800
PN-EK1210	PPAR gamma ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1211	PPAR gamma ELISA Kit with Nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1310	PPAR alpha ELISA Kit 1kit	1 kit	¥126,000	¥94,500
PN-EK1311	PPAR alpha ELISA Kit with Nuclear extraction kit 1kit	1 kit	¥137,000	¥102,800
PN-EK1400	FKHR ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK1401	FKHR ELISA Kit with Nuclear Extraction Kit	1 kit	¥137,000	¥102,800
PN-EK1500	Smad4 ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK1501	Smad4 ELISA Kit with Nuclear Extraction Kit	1 kit	¥137,000	¥102,800
PN-EK1600	Smad2 ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK1601	Smad2 ELISA Kit with Nuclear Extraction Kit	1 kit	¥137,000	¥102,800
PN-EK1700	EGR1 ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK1701	EGR1 ELISA Kit with Nuclear Extraction Kit	1 kit	¥137,000	¥102,800
PN-EK1800	HNF-3beta ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK1801	HNF-3beta ELISA Kit with Nuclear Extraction Kit	1 kit	¥137,000	¥102,800
PN-EK2000	E2F-1 ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK2001	E2F-1 ELISA Kit with Nuclear Extraction Kit	1 kit	¥137,000	¥102,800
PN-EK2010	Phosphotyrosine/PI3K ELISA Kit	1 kit	¥127,000	¥95,300
PN-EK2011	EGFR/PI3K ELISA Kit	1 kit	¥127,000	¥95,300
PN-EK3100	NFkB p52 ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK3101	NFkB p52 ELISA Kit with Nuclear Extraction kit	1 kit	¥137,000	¥102,800
PN-EK4100	NFkB c-Rel ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK4101	NFkB c-Rel ELISA Kit with Nuclear Extraction kit	1 kit	¥137,000	¥102,800
PN-EK5100	NFkB Rel-b ELISA Kit	1 kit	¥126,000	¥94,500
PN-EK5101	NFkB Rel-b ELISA Kit with Nuclear Extraction kit	1 kit	¥137,000	¥102,800
PN-EK0500	HIF-1 Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0501	HIF-1 alpha Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0502	ER Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0503	ER alpha Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0504	AP1 Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0505	AP1 Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0506	p53 Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0507	p53 Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0508	GR Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0509	GR Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0510	STAT1 Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0511	STAT1 Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0512	MEF2 Cold Probe for TF ELISA	本	¥18,000	¥13,500

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-EK0513	MEF2 Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0514	Sp1 Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0515	Sp1 Control Recombinant Protein for TF ELISA	本	¥18,000	¥13,500
PN-EK0516	PPAR gamma Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0517	PPAR gamma Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0518	PPAR alpha Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0519	PPAR alpha Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0520	GATA1 Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0521	GATA1 Control Recombinant Protein for TF ELISA	本	¥18,000	¥13,500
PN-EK0522	NFkB Cold Probe for TF ELISA	本	¥18,000	¥13,500
PN-EK0523	NFkB Control Nuclear Extract for TF ELISA	本	¥18,000	¥13,500
PN-EK0524	Nuclear Extract Dilution Buffer for TF ELISA	本	¥13,000	¥9,800
PN-EK2100	MITF ELISA Kit	個	¥126,000	¥94,500
PN-EK2101	MITF ELISA Kit with Nuclear Extraction Kit	個	¥137,000	¥102,800

<リン酸化 Array 試薬>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-MA3010	SH3 Domain Array I	1 kit	¥137,000	¥102,800
PN-MA3011	SH3 Domain Array II	1 kit	¥137,000	¥102,800
PN-MA3012	SH3 Domain Array III	1 kit	¥137,000	¥102,800
PN-MA3014	SH3 Domain Array IV	1 kit	¥137,000	¥102,800
PN-MA3020	PDZ Domain Array I	1 kit	¥137,000	¥102,800
PN-MA3022	PDZ Domain Array II	1 kit	¥137,000	¥102,800
PN-MA3023	PDZ Domain Array III	1 kit	¥137,000	¥102,800
PN-MA3024	PDZ Domain Array IV	1 kit	¥137,000	¥102,800
PN-MA3030	WW Domain Array I	1 kit	¥137,000	¥102,800
PN-MA3032	WW Domain Array II	1 kit	¥137,000	¥102,800
PN-MA3040	SH2 Domain Array	1 kit	¥137,000	¥102,800
PN-MA3041	Phosphotyrosine profiling array	1 kit	¥137,000	¥102,800
PN-MA3042	EGFR Tyrosine Phosphorylation Array	1 kit	¥137,000	¥102,800
PN-MA3043	PDGFR Tyrosine Phosphorylation Array	1 kit	¥137,000	¥102,800
PN-MA3044	IGF-IR Tyrosine Phosphorylation Array	1 kit	¥137,000	¥102,800

<Protein/DNA Array 関連>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-MA1210	Protein/DNA Array I	1 kit	¥188,000	¥141,000
PN-MA1211	Protein/DNA Array II	1 kit	¥198,000	¥148,500
PN-MA1212	Protein/DNA Array III	1 kit	¥198,000	¥148,500
PN-MA1213	Protein/DNA Array IV	1 kit	¥188,000	¥141,000
PN-MA1214	Protein/DNA Array V	1 kit	¥188,000	¥141,000
PN-MA1215	Combo Protein/DNA Array	1 kit	¥332,000	¥249,000
PN-MA1310	Protein/DNA I Refill	1 kit	¥106,000	¥79,500
PN-MA1311	Protein/DNA II Refill	1 kit	¥116,000	¥87,000
PN-MA1312	Protein/DNA III Refill	1 kit	¥116,000	¥87,000
PN-MA1313	Protein/DNA IV Refill Kit	1 kit	¥106,000	¥79,500
PN-MA1314	Protein/DNA V Refill Kit	1 kit	¥106,000	¥79,500
PN-MA1315	Combo Protein/DNA Refill Kit	1 kit	¥134,000	¥100,500
PN-MA1500	cAMP/Calcium Protein/DNA Array	1 kit	¥119,000	¥89,300
PN-MA1505	Cell Growth Protein/DNA Array	1 kit	¥119,000	¥89,300
PN-MA1510	Nuclear Receptor Protein/DNA Array	1 kit	¥119,000	¥89,300
PN-MA1600	cAMP/Calcium Protein/DNA Array Refill	1 kit	¥66,000	¥49,500
PN-MA1605	Cell Growth Protein/DNA Array Refill	1 kit	¥66,000	¥49,500
PN-MA1610	Nuclear Receptor Protein/DNA Array Refill	1 kit	¥66,000	¥49,500
PN-CS8533	Protein/DNA Array Spin Column Accessory Kit	個	¥38,000	¥28,500

<癌細胞分離キット>

商品コード	品名	梱包	通常定価	キャンペーン価格
PN-CI0002	Cancer cell isolation kit (2) 2assays	2 assay	¥75,000	¥56,300
PN-CI0004	Cancer cell isolation kit (4) 4assays	4 assay	¥131,000	¥98,300
PN-CI0010	Cancer cell isolation kit (10) 10assays	10 assay	¥283,000	¥212,300

<TF Array 試薬>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-MA4010	Transcription Reporter Array I	1 kit	¥135,000	¥101,300
PN-MA4100	Transcription Reporter Array I Refill	1 kit	¥79,000	¥59,300
PN-MA5010	TF/TF Interaction Array I	1 kit	¥185,000	¥138,800
PN-MA5011	TF/TF Interaction Array II	1 kit	¥196,000	¥147,000
PN-MA5012	TF/TF Interaction Array III	1 kit	¥196,000	¥147,000
PN-MA5041	PPARalpha-TF Interaction Array I	1 kit	¥185,000	¥138,800
PN-MA5042	PPARalpha-TF Interaction Array II	1 kit	¥196,000	¥147,000
PN-MA5043	PPARalpha-TF Interaction Array III	1 kit	¥196,000	¥147,000
PN-MA5044	PPARgamma-TF Interaction Array I	1 kit	¥185,000	¥138,800
PN-MA5045	PPARgamma-TF Interaction Array II	1 kit	¥196,000	¥147,000
PN-MA5046	PPARgamma-TF Interaction Array III	1 kit	¥196,000	¥147,000
PN-MA5051	NFkB p50-TF Interaction Array I	1 kit	¥185,000	¥138,800
PN-MA5052	NFkB p50-TF Interaction Array II	1 kit	¥196,000	¥147,000
PN-MA5053	NFkB p50-TF Interaction Array III	1 kit	¥196,000	¥147,000
PN-MA5054	NFkB p65-TF Interaction Array I	1 kit	¥185,000	¥138,800
PN-MA5055	NFkB p65-TF Interaction Array II	1 kit	¥196,000	¥147,000
PN-MA5056	NFkB p65-TF Interaction Array III	1 kit	¥196,000	¥147,000
PN-MA5110	TF/TF Interaction I Refill Kit	1 kit	¥103,000	¥77,300
PN-MA5111	TF/TF Interaction II Refill Kit	1 kit	¥114,000	¥85,500
PN-MA5112	TF/TF Interaction III Refill Kit	1 kit	¥114,000	¥85,500

<Stabe Cell Line>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-RC0002	NFkB-Luciferase A549 Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0003	Stat3-Luciferase HeLa Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0004	Stat1-Luciferase HeLa Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0005	NF-AT-Luciferase HeLa Reporter Stable Cell line	1VI	¥557,000	¥417,800
PN-RC0006	AP1 Luciferase 293 Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0007	CREB-Luciferase 293 Reporter Stable Cell line	1VI	¥557,000	¥417,800
PN-RC0008	CREB-Luciferase CHO Reporter Stable Cell line	1VI	¥557,000	¥417,800
PN-RC0009	NF-AT-Luciferase K562 Reporter Stable Cell line	1VI	¥557,000	¥417,800
PN-RC0012	SRF-Luciferase HeLa Reporter Stable Cell line	1VI	¥557,000	¥417,800
PN-RC0013	NFkB-Luciferase HeLa Reporter Stable Cell line	1VI	¥557,000	¥417,800
PN-RC0014	NFkB-Luciferase 293 Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0015	NFkB-Luciferase NIH3T3 Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0016	NFkB-Luciferase C2C12 Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0017	HIF-Luciferase NIH3T3 Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0018	293 GR-Luciferase Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC0019	HeLa GR-Luciferase Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC1001	TAD Kinase-luciferase HeLa Reporter Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC2001	CHO/GFP-NF κ Bp65 Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC2002	CHO/GFP-GR α Stable Cell Line	1VI	¥557,000	¥417,800
PN-RC5001	NFkB p50 Knockdown HeLa Stable Cell Line	1VI	¥557,000	¥417,800

<Cytokine Array 試薬>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-MA6110	Human Cytokine Antibody Array 1.0 (2 membranes)	1 kit	¥65,000	¥48,800
PN-MA6120	Human Cytokine Antibody Array 1.0 (4 membranes)	1 kit	¥75,000	¥56,300
PN-MA6130	Human Cytokine Antibody Array 1.0 (6 membranes)	1 kit	¥123,000	¥92,300
PN-MA6140	Human Cytokine Antibody Array 1.0 (8 membranes)	1 kit	¥136,000	¥102,000
PN-MA6150	Human Cytokine Antibody Array 3.0 (4 membranes)	1 kit	¥118,000	¥88,500
PN-MA6160	Human Cytokine Antibody Array 3.0 (8 membranes)	1 kit	¥218,000	¥163,500
PN-MA6310	Human Angiogenesis Antibody Array (4 membranes)	1 kit	¥97,000	¥72,800
PN-MA6320	Mouse Angiogenesis Antibody Array (4 membranes)	1 kit	¥97,000	¥72,800
PN-MA6410	Mouse Cytokine Antibody Array 1.0 (4 membranes)	1 kit	¥102,000	¥76,500
PN-MA6412	Mouse Cytokine Antibody Array 1.0 (8 membranes)	1 kit	¥160,000	¥120,000

<Procarta TF 試薬>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-PC5000	Procarta TF Assay Kit, 40-plex	1/4 plate	¥304,000	¥228,000
PN-PC5001	Procarta Transcription Factor Assay Kit, 40-plex	1 plate	¥924,000	¥693,000
PN-PC5010	Procarta TF Assay Kit, 43-plex	1/4 plate	¥300,000	¥225,000
PN-PC5011	Procarta TF Assay Kit, 43-plex	1 plate	¥929,000	¥696,800
PN-PC5012	Procarta Transcription Factor Assay Kit 10 Plex Kit	Kit	¥226,000	¥169,500
PN-PC5102	Procarta TF Whole Cell Lysis Kit	kit	¥63,000	¥47,300
PN-PC5500	Procarta TF Control Nuclear Extracts	kit	¥26,000	¥19,500
PN-PC5501	Procarta TF Control Whole Cell Lysates	kit	¥26,000	¥19,500
PN-PC5503	Procarta TF Binding Buffer	25 mL	¥10,000	¥7,500
PN-PC5504	Procarta TF Separation Buffer	200 mL	¥17,000	¥12,800
PN-PC5505	Procarta TF Wash Buffer	120 mL	¥17,000	¥12,800
PN-PC5506	Procarta TF Assay Buffer	5 mL	¥10,000	¥7,500
PN-PC5507	Procarta TF Reading Buffer	15 mL	¥10,000	¥7,500
PN-PC5508	Procarta TF Streptavidin-PE (SAPE)	11 mL	¥14,000	¥10,500
PN-PC5509	Procarta TF Elution Buffer	7 mL	¥10,000	¥7,500
PN-PC5510	Procarta TF Sample Collection Plate	枚	¥4,000	¥3,000
PN-PC5511	Procarta TF Utility Plate	枚	¥4,000	¥3,000
PN-PC5512	Procarta TF Separation Plate	枚	¥22,000	¥16,500
PN-PC5513	Procarta TF Plate Seals	3枚	¥6,000	¥4,500
PN-PC5515	Procarta TF PCR Plate Seal	3枚	¥6,000	¥4,500
PN-PC5516	Procarta TF PCR Plate	枚	¥4,000	¥3,000
PN-PC5900	Procarta TF Assay Kit, By Request, 3plex	1 plate	¥152,000	¥114,000
PN-PC5901	Procarta TF Assay Kit, By Request, 4plex	1 plate	¥196,000	¥147,000
PN-PC5902	Procarta TF Assay Kit, By Request, 5plex	1 plate	¥239,000	¥179,300
PN-PC5903	Procarta TF Assay Kit, By Request, 6plex	1 plate	¥280,000	¥210,000
PN-PC5904	Procarta TF Assay Kit, By Request, 7plex	1 plate	¥319,000	¥239,300
PN-PC5905	Procarta TF Assay Kit, By Request, 8plex	1 plate	¥357,000	¥267,800
PN-PC5906	Procarta TF Assay Kit, By Request, 9plex	1 plate	¥393,000	¥294,800
PN-PC5907	Procarta TF Assay Kit, By Request, 10plex	1 plate	¥427,000	¥320,300
PN-PC5908	Procarta TF Assay Kit, By Request, 11plex	1 plate	¥460,000	¥345,000
PN-PC5909	Procarta TF Assay Kit, By Request, 12plex	1 plate	¥491,000	¥368,300
PN-PC5910	Procarta TF Assay Kit, By Request, 13plex	1 plate	¥521,000	¥390,800
PN-PC5911	Procarta TF Assay Kit, By Request, 14plex	1 plate	¥549,000	¥411,800
PN-PC5912	Procarta TF Assay Kit, By Request, 15plex	1 plate	¥576,000	¥432,000
PN-PC5913	Procarta TF Assay Kit, By Request, 16plex	1 plate	¥603,000	¥452,300
PN-PC5914	Procarta TF Assay Kit, By Request, 17plex	1 plate	¥627,000	¥470,300
PN-PC5915	Procarta TF Assay Kit, By Request, 18plex	1 plate	¥650,000	¥487,500
PN-PC5916	Procarta TF Assay Kit, By Request, 19plex	1 plate	¥673,000	¥504,800
PN-PC5917	Procarta TF Assay Kit, By Request, 20plex	1 plate	¥694,000	¥520,500
PN-PC5918	Procarta TF Assay Kit, By Request, 21plex	1 plate	¥714,000	¥535,500
PN-PC5919	Procarta TF Assay Kit, By Request, 22plex	1 plate	¥732,000	¥549,000
PN-PC5920	Procarta TF Assay Kit, By Request, 23plex	1 plate	¥750,000	¥562,500
PN-PC5921	Procarta TF Assay Kit, By Request, 24plex	1 plate	¥767,000	¥575,300
PN-PC5922	Procarta TF Assay Kit, By Request, 25plex	1 plate	¥821,000	¥615,800
PN-PC5923	Procarta TF Assay Kit, By Request, 26plex	1 plate	¥836,000	¥627,000
PN-PC5924	Procarta TF Assay Kit, By Request, 27plex	1 plate	¥812,000	¥609,000
PN-PC5925	Procarta TF Assay Kit, By Request, 28plex	1 plate	¥825,000	¥618,800
PN-PC5926	Procarta TF Assay Kit, By Request, 29plex	1 plate	¥837,000	¥627,800
PN-PC5927	Procarta TF Assay Kit, By Request, 30plex	1 plate	¥849,000	¥636,800
PN-PC5928	Procarta TF Assay Kit, By Request, 31plex	1 plate	¥859,000	¥644,300
PN-PC5929	Procarta TF Assay Kit, By Request, 32plex	1 plate	¥869,000	¥651,800
PN-PC5930	Procarta TF Assay Kit, By Request, 33plex	1 plate	¥878,000	¥658,500
PN-PC5931	Procarta TF Assay Kit, By Request, 34plex	1 plate	¥887,000	¥665,300
PN-PC5932	Procarta TF Assay Kit, By Request, 35plex	1 plate	¥895,000	¥671,300
PN-PC5933	Procarta TF Assay Kit, By Request, 36plex	1 plate	¥902,000	¥676,500
PN-PC5934	Procarta TF Assay Kit, By Request, 37plex	1 plate	¥908,000	¥681,000
PN-PC5935	Procarta TF Assay Kit, By Request, 38plex	1 plate	¥914,000	¥685,500
PN-PC5936	Procarta TF Assay Kit, By Request, 39plex	1 plate	¥919,000	¥689,300
PN-PC5937	Procarta TF Assay Kit, By Request, 40plex	1 plate	¥924,000	¥693,000
PN-PC5938	Procarta TF Assay Kit, By Request, 41plex	1 plate	¥928,000	¥696,000
PN-PC5939	Procarta TF Assay Kit, By Request, 42plex	1 plate	¥931,000	¥698,300

〈プロモーターメチル化検出試薬〉

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-MP1100	Promoter Methylation PCR Kit	個	¥104,000	¥78,000
PN-MP1101	14-3-3 Forward PCR Primer (180bp position)	個	¥14,000	¥10,500
PN-MP1102	14-3-3 Reverse PCR Primer (180bp position)	個	¥14,000	¥10,500
PN-MP1103	14-3-3 Forward PCR Primer (600bp position)	個	¥14,000	¥10,500
PN-MP1104	14-3-3 Reverse PCR Primer (600bp position)	個	¥14,000	¥10,500
PN-MP1105	ABL1 Forward PCR Primer (390bp position)	個	¥14,000	¥10,500
PN-MP1106	ABL1 Reverse PCR Primer (390bp position)	個	¥14,000	¥10,500
PN-MP1107	BAGE Forward PCR Primer (230bp position)	個	¥14,000	¥10,500
PN-MP1108	BAGE Reverse PCR Primer (230bp position)	個	¥14,000	¥10,500
PN-MP1109	BAGE Forward PCR Primer (490bp position)	個	¥14,000	¥10,500
PN-MP1110	BAGE Reverse PCR Primer (490bp position)	個	¥14,000	¥10,500
PN-MP1111	BAGE Forward PCR Primer (230bp position)	個	¥14,000	¥10,500
PN-MP1112	BAGE Reverse PCR Primer (230bp position)	個	¥14,000	¥10,500
PN-MP1113	BCRA1 Forward PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1114	BCRA1 Reverse PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1115	BCRA1 Forward PCR Primer (130bp position)	個	¥14,000	¥10,500
PN-MP1116	BCRA1 Reverse PCR Primer (130bp position)	個	¥14,000	¥10,500
PN-MP1117	Calcitonin CGPR Forward PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1118	Calcitonin CGPR Reverse PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1119	Calcitonin CGPR Forward PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1120	Calcitonin CGPR Reverse PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1121	CD14 Forward PCR Primer (400bp position)	個	¥14,000	¥10,500
PN-MP1122	CD14 Reverse PCR Primer (400bp position)	個	¥14,000	¥10,500
PN-MP1123	CDKN2A Forward PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1124	CDKN2A Reverse PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1125	CDKN2A Forward PCR Primer (240bp position)	個	¥14,000	¥10,500
PN-MP1126	CDKN2A Reverse PCR Primer (240bp position)	個	¥14,000	¥10,500
PN-MP1127	COX2 Forward PCR Primer (320bp position)	個	¥14,000	¥10,500
PN-MP1128	COX2 Reverse PCR Primer (320bp position)	個	¥14,000	¥10,500
PN-MP1129	Cyclin D2 Forward PCR Primer (170bp position)	個	¥14,000	¥10,500
PN-MP1130	Cyclin D2 Reverse PCR Primer (170bp position)	個	¥14,000	¥10,500
PN-MP1131	Cyclin D2 Forward PCR Primer (130bp position)	個	¥14,000	¥10,500
PN-MP1132	Cyclin D2 Reverse PCR Primer (130bp position)	個	¥14,000	¥10,500
PN-MP1133	DAPK Forward PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1134	DAPK Reverse PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1135	DAPK Forward PCR Primer (500bp position)	個	¥14,000	¥10,500
PN-MP1136	DAPK Reverse PCR Primer (500bp position)	個	¥14,000	¥10,500
PN-MP1137	glut4 Forward PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1138	glut4 Reverse PCR Primer (200bp position)	個	¥14,000	¥10,500
PN-MP1139	H-ras Forward PCR Primer (400bp position)	個	¥14,000	¥10,500
PN-MP1140	H-ras Reverse PCR Primer (400bp position)	個	¥14,000	¥10,500
PN-MP1141	H-ras Forward PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1142	H-ras Reverse PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1143	INFgamma Forward PCR Primer (100bp position)	個	¥14,000	¥10,500
PN-MP1144	INFgamma Reverse PCR Primer (100bp position)	個	¥14,000	¥10,500
PN-MP1145	GAGE1 Forward PCR Primer (250bp position)	個	¥14,000	¥10,500
PN-MP1146	GAGE1 Reverse PCR Primer (250bp position)	個	¥14,000	¥10,500
PN-MP1147	GAGE1 Forward PCR Primer (400bp position)	個	¥14,000	¥10,500
PN-MP1148	GAGE1 Reverse PCR Primer (400bp position)	個	¥14,000	¥10,500
PN-MP1149	HOXA2 Forward PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1150	HOXA2 Reverse PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1151	IL4 Forward PCR Primer (230bp position)	個	¥14,000	¥10,500
PN-MP1152	IL4 Reverse PCR Primer (230bp position)	個	¥14,000	¥10,500
PN-MP1153	IRF7 Forward PCR Primer (450bp position)	個	¥14,000	¥10,500
PN-MP1154	IRF7 Reverse PCR Primer (450bp position)	個	¥14,000	¥10,500
PN-MP1155	NIS Forward PCR Primer (350bp position)	個	¥14,000	¥10,500
PN-MP1156	NIS Reverse PCR Primer (350bp position)	個	¥14,000	¥10,500
PN-MP1157	NIS Forward PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1158	NIS Reverse PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1159	NME2 Forward PCR Primer (170bp position)	個	¥14,000	¥10,500
PN-MP1160	NME2 Reverse PCR Primer (170bp position)	個	¥14,000	¥10,500
PN-MP1161	NME2 Forward PCR Primer (500bp position)	個	¥14,000	¥10,500
PN-MP1162	NME2 Reverse PCR Primer (500bp position)	個	¥14,000	¥10,500
PN-MP1163	PDGFB Forward PCR Primer (500bp position)	個	¥14,000	¥10,500
PN-MP1164	PDGFB Reverse PCR Primer (500bp position)	個	¥14,000	¥10,500
PN-MP1165	POMC Forward PCR Primer (490bp position)	個	¥14,000	¥10,500
PN-MP1166	POMC Reverse PCR Primer (490bp position)	個	¥14,000	¥10,500
PN-MP1167	NES Forward PCR Primer (290bp position)	個	¥14,000	¥10,500
PN-MP1168	NES Reverse PCR Primer (290bp position)	個	¥14,000	¥10,500
PN-MP1169	NES Forward PCR Primer (350bp position)	個	¥14,000	¥10,500
PN-MP1170	NES Reverse PCR Primer (350bp position)	個	¥14,000	¥10,500
PN-MP1171	PRA2 Forward PCR Primer (250bp position)	個	¥14,000	¥10,500
PN-MP1172	PRA2 Reverse PCR Primer (250bp position)	個	¥14,000	¥10,500
PN-MP1173	TFF1 Forward PCR Primer (160bp position)	個	¥14,000	¥10,500

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-MP1174	TFF1 Reverse PCR Primer (160bp position)	個	¥14,000	¥10,500
PN-MP1175	VHL Forward PCR Primer (290bp position)	個	¥14,000	¥10,500
PN-MP1176	VHL Reverse PCR Primer (290bp position)	個	¥14,000	¥10,500
PN-MP1177	Survivin Forward PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1178	Survivin Reverse PCR Primer (300bp position)	個	¥14,000	¥10,500
PN-MP1179	Caspase 8 Forward PCR Primer (360bp position)	個	¥14,000	¥10,500
PN-MP1180	Caspase 8 Reverse PCR Primer (360bp position)	個	¥14,000	¥10,500

〈タンパク質 ELISA キット〉

商品コード	品名	梱包	通常定価	キャンペーン価格
PN-BC1009	AFP ELISA Kit	1 kit	¥66,000	¥49,500
PN-BC1013	CA 125 ELISA Kit	1 kit	¥77,000	¥57,800
PN-BC1015	CA 15-3 ELISA Kit	1 kit	¥77,000	¥57,800
PN-BC1017	CA 19-9 ELISA Kit	1 kit	¥77,000	¥57,800
PN-K1001-1	Human Adiponectin ELISA Kit		¥141,000	¥105,800
PN-K1002-1	Mouse/Rat Adiponectin ELISA Kit		¥141,000	¥105,800

<EMSA Probe 及びキット>

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1000	EMSA kit without probe	25 rxn	¥58,000	¥43,500
PN-AY1001	AP1(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1001P	AP1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1002	AP2(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1002P	AP2(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1003	ARE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1003P	ARE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1004	Brn3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1004P	Brn3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1005	C/EBP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1005P	C/EBP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1006	CBF(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1006P	CBF(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1007	CDP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1007P	CDP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1008	c-Myb EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1008P	c-Myb EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1009	CREB EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1009P	CREB EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1010	E2F1(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1010P	E2F1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1011	EGR EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1011P	EGR EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1012	ERE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1012P	ERE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1013	Ets EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1013P	Ets EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1014	Ets-1/PEA3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1014P	Ets-1/PEA3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1015	FAST1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1015P	FAST1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1016	GAS/ISRE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1016P	GAS/ISRE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1017	GATA EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1017P	GATA EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1018	GRE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1018P	GRE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1019	HNF4(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1019P	HNF4(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1020	HSE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1020P	HSE EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1021	IRF1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1021P	IRF1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1022	MEF1(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1022P	MEF1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1023	MEF2(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1023P	MEF2(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1024	MRE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1024P	MRE EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1025	Myc-Max EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1025P	Myc-Max EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1026	NF1(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1026P	NF1(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1027	NF-ATc EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1027P	NF-ATc EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1028	NF-E1/Y1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1028P	NF-E1/Y1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1029	NF-E2(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1029P	NF-E2(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1030	NFkB(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1030P	NFkB(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1031	Oct-1(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1031P	Oct-1(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1032	p53(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1032P	p53(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1033	Pax5(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1033P	Pax5(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1034	Pbx1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1034P	Pbx1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1035	Pit1(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1035P	Pit1(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1036	PPAR EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1036P	PPAR EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1037	PRE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1037P	PRE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1038	RAR (DR5) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1038P	RAR (DR5) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1039	RXR (DR1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1039P	RXR (DR1) EMSA Probe Set	25 rxn	¥25,000	¥18,800

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1040	SIE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1040P	SIE EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1041	Smad SBE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1041P	Smad SBE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1042	Smad3/Smad4 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1042P	Smad3/Smad4 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1043	Sp1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1043P	Sp1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1044	SRE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1044P	SRE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1045	Stat1-p84/p91 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1045P	Stat1-p84/p91 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1046	Stat3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1046P	Stat3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1047	Stat4 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1047P	Stat4 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1048	Stat5 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1048P	Stat5 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1049	Stat5/Stat6 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1049P	Stat5/Stat6 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1050	TFIID EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1050P	TFIID EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1051	TR EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1051P	TR EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1052	TR(DR4) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1052P	TR(DR4) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1053	USF1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1053P	USF1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1054	VDR(D3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1054P	VDR(D3) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1055	AP1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1055P	AP1(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1056	AP2(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1056P	AP2(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1059	ADR1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1059P	ADR1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1060	AhR/Arnt EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1060P	AhR/Arnt EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1061	Angiotensinogen ANG-IRE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1061P	Angiotensinogen ANG-IRE EMSA Probe Set	25 rxn	¥23,000	¥17,300
PN-AY1062	Antioxidant RE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1062P	Antioxidant RE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1063	AP3(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1063P	AP3(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1064	AP4(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1064P	AP4(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1065	Beta-response element EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1065P	Beta-response element EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1066	CEF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1066P	CEF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1067	CEF2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1067P	CEF2 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1068	CTEF-CRE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1068P	CTEF-CRE EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1069	COUP-TF EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1069P	COUP-TF EMSA Probe Set	25 rxn	¥23,000	¥17,300
PN-AY1070	CREB-BP1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1070P	CREB-BP1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1071	c-Rel EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1071P	c-Rel EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1073	DR0 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1073P	DR0 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1074	DR1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1074P	DR1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1075	DR2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1075P	DR2 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1076	DR3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1076P	DR3 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1077	DR4 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1077P	DR4 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1078	DR5 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1078P	DR5 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1079	E47 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1079P	E47 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1080	E4BP4(1) EMSA Kit	25 rxn	¥78,000	¥58,500

商品コード	品名	梱包	通常定価	キャンペーン 定価
PN-AY1080P	E4BP(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1081	E4F/ATF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1081P	E4F/ATF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1082	Elk1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1082P	Elk1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1083	EV1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1083P	EV1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1084	Frec2(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1084P	Frec2(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1085	FKHR(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1085P	FKHR(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1086	FKHR(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1086P	FKHR(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1087	Afxh (Foxo4) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1087P	Afxh (Foxo4) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1088	Frec2(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1088P	Frec2(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1089	Frec4 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1089P	Frec4 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1090	Frec7 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1090P	Frec7 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1091	GAG EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1091P	GAG EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1092	GAS EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1092P	GAS EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1093	GATA1(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1093P	GATA1(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1094	GATA1/GATA2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1094P	GATA1/GATA2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1095	GATA2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1095P	GATA2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1096	GATA3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1096P	GATA3 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1097	GATA4 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1097P	GATA4 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1098	GATA6 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1098P	GATA6 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1099	Gfi1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1099P	Gfi1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1100	HBS (hypoxia-inducible) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1100P	HBS (hypoxia-inducible) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1101	HBS/HAS EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1101P	HBS/HAS EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1102	HFH1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1102P	HFH1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1103	HFH2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1103P	HFH2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1104	HFH3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1104P	HFH3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1105	HFH8 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1105P	HFH8 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1107	HIF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1107P	HIF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1108	HAS EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1108P	HAS EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1109	HNF3 alpha / HNF3 beta / HNF3gamma EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1109P	HNF3 alpha / HNF3 beta / HNF3gamma EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1110	HNF3 beta EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1110P	HNF3 beta EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1111	HNF4(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1111P	HNF4(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1112	HNF4/COUP-TF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1112P	HNF4/COUP-TF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1113	MT-Box EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1113P	MT-Box EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1114	Ikaros EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1114P	Ikaros EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1115	IRF1/IRF2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1115P	IRF1/IRF2 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1116	ISRE(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1116P	ISRE(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1117	ISRE(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1117P	ISRE(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1118	ISRE(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1118P	ISRE(3) EMSA Probe Set	25 rxn	¥25,000	¥18,800

商品コード	品名	梱包	通常定価	キャンペーン 定価
PN-AY1119	KKLF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1119P	KKLF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1120	L-III BP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1120P	L-III BP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1121	MEF2(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1121P	MEF2(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1122	MEF2(3) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1122P	MEF2(3) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1123	MEF3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1123P	MEF3 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1124	MSP1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1124P	MSP1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1125	MUSF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1125P	MUSF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1126	MZF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1126P	MZF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1127	NF-Y/Y-box binding factor EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1127P	NF-Y/Y-box binding factor EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1128	Nkx-2.5 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1128P	Nkx-2.5 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1129	NPAS2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1129P	NPAS2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1131	NZF3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1131P	NZF3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1132	ORE (osmotic response element) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1132P	ORE (osmotic response element) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1133	p300 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1133P	p300 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1134	PARP (poly(ADP-ribose) synthetase/polymerase) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1134P	PARP (poly(ADP-ribose) synthetase/polymerase) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1135	Pax2(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1135P	Pax2(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1136	Pax3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1136P	Pax3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1137	Pax4 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1137P	Pax4 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1138	Pax6(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1138P	Pax6(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1139	Pax8 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1139P	Pax8 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1140	PPARalpha EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1140P	PPARalpha EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1141	PPARgamma EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1141P	PPARgamma EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1142	RREB(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1142P	RREB(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1143	RREB(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1143P	RREB(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1144	RSRFC4 (related to serum factor, C4) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1144P	RSRFC4 (related to serum factor, C4) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1145	SAA EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1145P	SAA EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1146	Skn EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1146P	Skn EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1147	SRY(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1147P	SRY(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1148	Tax/CREB EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1148P	Tax/CREB EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1149	TCF/LEF EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1149P	TCF/LEF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1150	XBP1(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1150P	XBP1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1151	XRE(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1151P	XRE(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1152	ZIA EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1152P	ZIA EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1153	ZIB EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1153P	ZIB EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1154	ZIC EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1154P	ZIC EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1155	ZID EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1155P	ZID EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1156	Cdx4/NKX2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1156P	Cdx4/NKX2 EMSA Probe Set	25 rxn	¥25,000	¥18,800

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1157	HBS (HTF binding sequence) (rat, as human xbp-1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1157P	HBS (HTF binding sequence) (rat, as human xbp-1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1158	ADD1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1158P	ADD1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1159	AF1/ARP1/NF-BA EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1159P	AF1/ARP1/NF-BA EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1160	AFP1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1160P	AFP1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1161	AIC/CBF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1161P	AIC/CBF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1162	AML1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1162P	AML1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1163	AP2/YY1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1163P	AP2/YY1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1164	AREB6 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1164P	AREB6 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1165	ARP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1165P	ARP EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1166	ATF/CRE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1166P	ATF/CRE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1167	ATF2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1167P	ATF2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1168	CACC EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1168P	CACC EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1169	CCAAT EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1169P	CCAAT EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1170	CCAC EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1170P	CCAC EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1171	CD28RC EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1171P	CD28RC EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1172	CD28RC/NF-IL2B EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1172P	CD28RC/NF-IL2B EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1174	CEA EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1174P	CEA EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1175	CP1/CTF/CBTF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1175P	CP1/CTF/CBTF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1176	CPE EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1176P	CPE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1177	CREB2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1177P	CREB2 EMSA Probe Set	25 rxn	¥319,000	¥239,300
PN-AY1178	CSBP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1178P	CSBP EMSA Probe Set	25 rxn	¥23,000	¥17,300
PN-AY1179	CTCF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1179P	CTCF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1180	E12(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1180P	E12(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1181	EBP40/EBP45 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1181P	EBP40/EBP45 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1182	EKLF(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1182P	EKLF(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1183	EKLF(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1183P	EKLF(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1184	Elf1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1184P	Elf1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1185	H4TF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1185P	H4TF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1186	HEN1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1186P	HEN1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1187	HFH8/HNF3/LUN EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1187P	HFH8/HNF3/LUN EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1188	HiNF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1188P	HiNF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1189	HLF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1189P	HLF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1190	HMG EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1190P	HMG EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1192	HNF1alpha EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1192P	HNF1alpha EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1193	HOX4C EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1193P	HOX4C EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1194	HOXD8/HOXD9/HOXD10 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1194P	HOXD8/HOXD9/HOXD10 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1195	HOXD8(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1195P	HOXD8(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1196	HOXD8(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1196P	HOXD8(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1197	HOXD9/HOXD10 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1197P	HOXD9/HOXD10 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1198	ICSBP EMSA Kit	25 rxn	¥80,000	¥60,000

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1198P	ICSBP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1199	Isl1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1199P	Isl1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1200	KPF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1200P	KPF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1201	KTP1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1201P	KTP1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1202	LCR-F1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1202P	LCR-F1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1203	LEF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1203P	LEF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1204	LFA1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1204P	LFA1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1205	LR1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1205P	LR1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1206	LyF1(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1206P	LyF1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1207	LyF1(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1207P	LyF1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1208	LyF1(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1208P	LyF1(3) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1209	MAZ(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1209P	MAZ(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1210	MAZ(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1210P	MAZ(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1211	MTB-Zf EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1211P	MTB-Zf EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1212	MTF(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1212P	MTF(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1213	MyoD EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1213P	MyoD EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1214	MyoG EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1214P	MyoG EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1215	MyTi EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1215P	MyTi EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1216	MZF1(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1216P	MZF1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1217	MZF1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1217P	MZF1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1218	MZF1(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1218P	MZF1(3) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1219	NF-4FA EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1219P	NF-4FA EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1220	NF-Gma(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1220P	NF-Gma(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1221	NF-Gma(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1221P	NF-Gma(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1222	NF-IL2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1222P	NF-IL2 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1223	NF-Y EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1223P	NF-Y EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1224	NRF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1224P	NRF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1225	NRF2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1225P	NRF2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1226	p53(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1226P	p53(3) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1227	p55 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1227P	p55 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1228	PEA3/HIP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1228P	PEA3/HIP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1229	PEBP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1229P	PEBP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1230	PEBP2(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1230P	PEBP2(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1231	PPUR(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1231P	PPUR(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1232	PPUR(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1232P	PPUR(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1233	PuF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1233P	PuF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1234	Pur1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1234P	Pur1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1235	PYR EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1235P	PYR EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1236	RB EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1236P	RB EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1237	RFX1/RFX2/RFX3(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1237P	RFX1/RFX2/RFX3(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1238	RFX1/RFX2/RFX3(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1238P	RFX1/RFX2/RFX3(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1239	RIPE3a1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1239P	RIPE3a1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1240	SIF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1240P	SIF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1241	SIF2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1241P	SIF2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1242	SIF3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1242P	SIF3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1243	SP1/ASP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1243P	SP1/ASP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1244	SRF/SAP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1244P	SRF/SAP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1245	SRY(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1245P	SRY(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1246	TCE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1246P	TCE EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1247	TFE3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1247P	TFE3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1248	Tf-LF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1248P	Tf-LF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1249	TGT3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1249P	TGT3 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1250	TIF1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1250P	TIF1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1251	TREF1/TREF2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1251P	TREF1/TREF2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1252	v-Maf EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1252P	v-Maf EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1253	WT1(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1253P	WT1(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1254	WT1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1254P	WT1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1255	WT1(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1255P	WT1(3) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1256	X2BP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1256P	X2BP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1257	XBP1/X2BP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1257P	XBP1/X2BP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1258	ZNF174 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1258P	ZNF174 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1259	AAF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1259P	AAF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1260	ABF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1260P	ABF1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1261	ACF EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1261P	ACF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1262	ACBPB EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1262P	ACBPB EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1263	AIC2/AIC3/AIC4 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1263P	AIC2/AIC3/AIC4 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1264	ALF1B EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1264P	ALF1B EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1265	alpha-PAL EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1265P	alpha-PAL EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1266	AP2(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1266P	AP2(3) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1267	AP3(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1267P	AP3(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1268	AP4(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1268P	AP4(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1269	apoA1 promoter EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1269P	apoA1 promoter EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1270	ATFalpha/ATFdelta EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1270P	ATFalpha/ATFdelta EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1271	ATFalpha EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1271P	ATFalpha EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1272	betaM-globin factor B1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1272P	betaM-globin factor B1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1273	BZP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1273P	BZP EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1274	C/EBPalph/C/EBPgamma EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1274P	C/EBPalph/C/EBPgamma EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1275	C/EBPalph(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1275P	C/EBPalph(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1276	C/EBPalph(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1276P	C/EBPalph(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1277	C/EBPgamma EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1277P	C/EBPgamma EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1278	CBFB EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1278P	CBFB EMSA Probe Set	25 rxn	¥29,000	¥21,800

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1279	CBF(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1279P	CBF(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1280	c-Ets-1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1280P	c-Ets-1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1281	c-fos enhancer-binding protein EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1281P	c-fos enhancer-binding protein EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1282	c-Myb binding protein EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1282P	c-Myb binding protein EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1283	c-Myc EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1283P	c-Myc EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1284	COUP-TF/EAR3/NR2F1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1284P	COUP-TF/EAR3/NR2F1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1285	CP1B EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1285P	CP1B EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1286	CP1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1286P	CP1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1287	CP2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1287P	CP2 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1288	CREB(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1288P	CREB(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1289	CRE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1289P	CRE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1290	CYP1A1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1290P	CYP1A1 EMSA Probe Set	25 rxn	¥23,000	¥17,300
PN-AY1291	DE I EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1291P	DE I EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1292	dioxin receptor EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1292P	dioxin receptor EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1293	E12/E47/MF3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1293P	E12/E47/MF3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1294	E12/E47 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1294P	E12/E47 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1295	E12(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1295P	E12(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1296	E2F1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1296P	E2F1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1297	E2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1297P	E2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1298	E4BP4(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1298P	E4BP4(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1299	EBP-80 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1299P	EBP-80 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1300	EBP-C2/NF-muE3/YEB3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1300P	EBP-C2/NF-muE3/YEB3 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1301	EGF binding protein EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1301P	EGF binding protein EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1302	EGR1/EGR3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1302P	EGR1/EGR3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1303	EGR1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1303P	EGR1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1304	EGR2 BP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1304P	EGR2 BP EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1305	EGR2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1305P	EGR2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1306	EIL1/EIL2/EIL3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1306P	EIL1/EIL2/EIL3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1307	ETF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1307P	ETF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1308	ETV4/E1A-F EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1308P	ETV4/E1A-F EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1309	Fra1/JUN EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1309P	Fra1/JUN EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1310	gastrin EGF-response element EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1310P	gastrin EGF-response element EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1311	GATA1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1311P	GATA1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1312	GBF1/GBF2/GBF3/HY5 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1312P	GBF1/GBF2/GBF3/HY5 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1313	GCF EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1313P	GCF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1314	GKLF EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1314P	GKLF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1315	H1TF2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1315P	H1TF2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1316	H4TF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1316P	H4TF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1317	HFH11beta/HFH11alpha EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1317P	HFH11beta/HFH11alpha EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1318	HFH EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1318P	HFH EMSA Probe Set	25 rxn	¥25,000	¥18,800

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1319	HiNF-A EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1319P	HiNF-A EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1320	HiNF-B/H1TF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1320P	HiNF-B/H1TF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1321	HiNF-C EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1321P	HiNF-C EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1322	HiNF-D3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1322P	HiNF-D3 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1323	HNF1a/HNF1b/HNF1c EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1323P	HNF1a/HNF1b/HNF1c EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1324	HNF1(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1324P	HNF1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1325	HNF4alpha1/HNF4alpha2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1325P	HNF4alpha1/HNF4alpha2 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1326	HOXA4 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1326P	HOXA4 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1327	HSF(2) (distal) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1327P	HSF(2) (distal) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1328	HSF(3) (proximal) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1328P	HSF(3) (proximal) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1329	HSF/HTF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1329P	HSF/HTF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1330	IFI-56KBP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1330P	IFI-56KBP EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1331	IL6 RE-BP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1331P	IL6 RE-BP EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1332	IR4 (EBV) binding protein EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1332P	IR4 (EBV) binding protein EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1333	ISGF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1333P	ISGF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1334	kBF-alpha EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1334P	kBF-alpha EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1335	lactoferrin BP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1335P	lactoferrin BP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1336	LFA1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1336P	LFA1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1337	LFA2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1337P	LFA2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1338	LFB2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1338P	LFB2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1339	LH2/Lim1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1339P	LH2/Lim1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1340	LSF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1340P	LSF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1341	LXRE1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1341P	LXRE1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1342	MASH1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1342P	MASH1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1343	MBP1(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1343P	MBP1(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1344	MBP1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1344P	MBP1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1345	mck promoter E1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1345P	mck promoter E1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1346	mck promoter E2 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1346P	mck promoter E2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1347	MDBP(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1347P	MDBP(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1348	MDBP(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1348P	MDBP(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1349	MEF2alpha EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1349P	MEF2alpha EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1350	Mfh1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1350P	Mfh1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1351	MHC gene promoter W box EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1351P	MHC gene promoter W box EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1352	msx1/msx2/msx3 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1352P	msx1/msx2/msx3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1353	Myb(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1353P	Myb(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1354	Myc-CF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1354P	Myc-CF1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1355	Myc-PRF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1355P	Myc-PRF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1356	NCAM binding protein EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1356P	NCAM binding protein EMSA Probe Set	25 rxn	¥25,000	¥18,800

商品コード	品名	梱包	通常定価	キャンペーン定価
PN-AY1357	NF1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1357P	NF1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1358	NF1/L EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1358P	NF1/L EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1359	NF-A3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1359P	NF-A3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1360	NF-A EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1360P	NF-A EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1361	NF-AT1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1361P	NF-AT1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1362	NF-Atp EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1362P	NF-Atp EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1363	NF-Atx EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1363P	NF-Atx EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1364	NF-E2(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1364P	NF-E2(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1365	NFE-6/CP1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1365P	NFE-6/CP1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1366	NFkB(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1366P	NFkB(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1367	NFkBp65 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1367P	NFkBp65 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1368	N-ras binding protein EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1368P	N-ras binding protein EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1369	OBP EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1369P	OBP EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1370	Oct-1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1370P	Oct-1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1371	Oct-4/Oct-3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1371P	Oct-4/Oct-3 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1372	Oct-B2 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1372P	Oct-B2 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1373	ODC EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1373P	ODC EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1374	p53(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1374P	p53(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1375	Pax1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1375P	Pax1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1376	Pax2(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1376P	Pax2(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1377	Pax5(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1377P	Pax5(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1378	Pax6(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1378P	Pax6(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1379	PBGD binding protein EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1379P	PBGD binding protein EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1380	PCF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1380P	PCF EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1381	PEBP2(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1381P	PEBP2(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1382	PEPCK promoter GR EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1382P	PEPCK promoter GR EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1383	Pit1(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1383P	Pit1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1384	Pit1(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1384P	Pit1(3) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1385	PO-B EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1385P	PO-B EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1386	PRDI-BFc EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1386P	PRDI-BFc EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1387	PRDII-BF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1387P	PRDII-BF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1388	PREB EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1388P	PREB EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1389	PTF1-beta EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1389P	PTF1-beta EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1390	PTF1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1390P	PTF1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1391	Pu.1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1391P	Pu.1 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1392	PUR EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1392P	PUR EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1393	RAR/RXR/CF1/MB67 EMSA Kit	25 rxn	¥80,000	¥60,000

商品コード	品名	梱包	通常定価	キャンペーン 価格
PN-AY1393P	RAR/RXR/CF1/MB67 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1394	RFX1/RFX2/RFX3(3) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1394P	RFX1/RFX2/RFX3(3) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1395	RORE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1395P	RORE EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1396	RVF EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1396P	RVF EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1397	SAP-1b EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1397P	SAP-1b EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1398	Snail EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1398P	Snail EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1399	Sp1(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1399P	Sp1(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1400	SPERM1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1400P	SPERM1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1401	SRF(2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1401P	SRF(2) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1402	SRF(1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1402P	SRF(1) EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1403	SSAP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1403P	SSAP EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1404	Stat1/Stat3 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1404P	Stat1/Stat3 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1405	Stat5b EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1405P	Stat5b EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1406	Surf-2 regulating protein (1) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1406P	Surf-2 regulating protein (1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1407	Surf-2 regulating protein (2) EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1407P	Surf-2 regulating protein (2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1408	T3R EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1408P	T3R EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1409	Tat EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1409P	Tat EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1410	TEF1/AP5 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1410P	TEF1/AP5 EMSA Probe Set	25 rxn	¥25,000	¥18,800
PN-AY1411	TEF1 EMSA Kit	25 rxn	¥80,000	¥60,000
PN-AY1411P	TEF1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1412	TFE3-L EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1412P	TFE3-L EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1413	TFEB EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1413P	TFEB EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1414	Thy1 binding protein EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1414P	Thy1 binding protein EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1415	tPA binding protein EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1415P	tPA binding protein EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1416	transferrin BP EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1416P	transferrin BP EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1417	TTF1(1) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1417P	TTF1(1) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1418	TTF1(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1418P	TTF1(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1419	TxREF/NF-III EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1419P	TxREF/NF-III EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1420	URE EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1420P	URE EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1421	USF(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1421P	USF(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1422	v-rel 50-55K EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1422P	v-rel 50-55K EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1423	WAP binding protein EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1423P	WAP binding protein EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1424	XBP1(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1424P	XBP1(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1425	XRE(2) EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1425P	XRE(2) EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1426	YB1 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1426P	YB1 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1427	Yi EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1427P	Yi EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1428	ZEB EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1428P	ZEB EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY1429	Zn15 EMSA Kit	25 rxn	¥84,000	¥63,000
PN-AY1429P	Zn15 EMSA Probe Set	25 rxn	¥29,000	¥21,800
PN-AY2002	Nuclear Extraction kit	1 kit	¥38,000	¥28,500
PN-AY2005	Genomic DNA Extraction Kit 1kit	1 kit	¥40,000	¥30,000
PN-AY9999	EMSA Combo Kit (any three probes)	25 rxn	¥120,000	¥90,000

商品コード	品名	梱包	通常定価	キャンペーン 定価
PN-CS7010	Control Nuclear Extract for EMSA Kit	本	¥15,000	¥11,300
PN-CS7026	Loading Dye for EMSA Kit	本	¥12,000	¥9,000
PN-CS7149	Control Probe for EMSA Kit	本	¥12,000	¥9,000
PN-CS7150	Cold Control Probe for EMSA Kit	本	¥12,000	¥9,000
PN-CS7526	Hybridization Buffer for EMSA Kit	本	¥15,000	¥11,300
PN-CS7527	20X SSC for Protein/DNA Arrays	本	¥15,000	¥11,300
PN-CS7528	20% SDS for Protein/DNA Arrays	本	¥15,000	¥11,300
PN-CS7535	2X Blocking Buffer for EMSA Kit	本	¥18,000	¥13,500
PN-CS7536	Streptavidin-HRP Conjugate for EMSA Kit	本	¥12,000	¥9,000
PN-CS7537	10X Wash Buffer for EMSA Kit	本	¥15,000	¥11,300
PN-CS7538	10X Detection Buffer for EMSA Kit	本	¥15,000	¥11,300
PN-CS7539	Substrate Solution I for EMSA Kit	本	¥21,000	¥15,800
PN-CS7540	Substrate Solution II for EMSA Kit	本	¥21,000	¥15,800
PN-CS7541	Substrate Solution III for EMSA Kit	本	¥15,000	¥11,300
PN-CS7542	5X Binding Buffer for EMSA Kit	本	¥12,000	¥9,000
PN-CS7543	Poly d(I-C) for EMSA Kit	本	¥15,000	¥11,300
PN-AY1191	HNF1 EMSA Kit	個	¥84,000	¥63,000
PN-AY1191P	NF1 EMSA Probe Set	個	¥29,000	¥21,800



日本総代理店

株式会社

ベリタス

〒105-0001 東京都港区虎ノ門 2-7-14 八洲ビル
TEL.03-3593-3211(代) FAX.03-3593-3216
E-mail: veritas@veritastk.co.jp

<http://www.veritastk.co.jp/>

RPNF-3045