

**P A C K A G E   N O T I C E**

**Catalog ID: PEC1Q**  
**Product Name: C1qScreen (RUO)**  
**LOT# 018**

Based on our internal testing data, please note the following information regarding Lot# 018 of C1qScreen (PEC1Q, RUO):

- Lot# 018 used with C1Q-NC Lot# 011 may exhibit background values averaging 500 MFI.
- Lot# 018 demonstrates an average 37% decrease in signal strength compared to Lot# 017.

**For additional guidance please contact One Lambda Technical Support:**  
**E-mail [1lambda-techsupport@thermofisher.com](mailto:1lambda-techsupport@thermofisher.com); Toll-free: 800.822.8824**