

		typical area-intensity values of product-groups with a PMT voltage at 600 V.					
excitation wavelenght		488 nm	488 nm	488 nm	488 nm	633 nm	633 nm
emission filter (nm)		530(30)	575(26)	695(40)	780(60)	660(20)	780(60)
particle	voltage (V)	Fitc	Pe	PerCP-Cy5.5	Pe-Cy7	APC	APC-Cy7
CB	600	250	30	15	2	100	3
C570	600	570	80	25	4	50	0
M-270 analogs	600	3000	760	1000	400	200	50
M-280 with antibody	600	3100	690	600	190	<b>500</b>	20
M-280 streptavidine	600	3100	690	600	190	<b>200</b>	20
M 450 analogs	600	2300	740	1400	600	400	150
ACA 100 analogs	600	120	30	10	4	5	0