

INVITROGRO™ CP MEDIUM

Product No.	Description	Size
Z99029	INVITROGRO CP Medium	250 mL
Z990003	INVITROGRO CP Medium	500 mL
Z990004	INVITROGRO CP Medium	1 L

Product Description:

INVITROGRO CP Medium is a hepatocyte incubation medium used for culturing primary hepatocytes¹. INVITROGRO CP Medium is formulated without antibiotics to increase its shelf life. BioIVT recommends the use of the TORPEDO™ Antibiotic Mix with INVITROGRO CP Medium.

Stability: 6 months for INVITROGRO CP Medium without TORPEDO Antibiotic Mix
7 days for Complete INVITROGRO CP Medium

Storage: 4°C

Related Products:

Product No.	Size
Z99000	TORPEDO Antibiotic Mix 5.5 ml
Z990007	TORPEDO Antibiotic Mix 11 mL
Z990008	TORPEDO Antibiotic Mix 22 mL

Reference:

1. Roymans, D.; Van Looveren, C.; Leone, A.; Parker, J. B.; McMillan, M.; Johnson, M. D.; Koganti, A.; Gilissen, R.; Silber, P.; Mannens, G.; Meuldermans, W. Determination of cytochrome P450 1A2 and cytochrome P450 3A4 induction in cryopreserved human hepatocytes. *Biochem. Pharmacol.* **2004**, 67(3), 427–437.

Caution: This product is being sold for research and/or manufacturing purposes only. The biological samples supplied by BioIVT, or any material isolated from the samples, are for in-vitro research use only and are not to be used as a source of material for clinical therapies. Human material may be used in vivo in animals. The user assumes all responsibility for its usage and disposal, in accordance with all regulations.

As of March 5th, 2018, "BioIVT" will replace the current brand of BioreclamationIVT and those of each of its affiliates and subsidiaries: IVT Holdings Inc., Asterand Bioscience, Inc., Asterand UK Acquisition limited, Asterand US Acquisition Corporation, Labquote.com LLC, ILSBio LLC, Seratrials LLC, Sera Laboratories International Ltd., Translational Cell Science LLC, Bioreclamation-IVT India Private Limited, In Vitro Inc., In Vitro GmbH, Qualyst Holdings LLC, Qualyst Transporter Solutions LLC.